VersoCit-2


#50300114

2016.11.01

Instructions

Anleitung

Mode d'emploi

 **Struers**

English

Description

VersoCit-2 is a two component acrylic mounting system, consisting of one powder and one liquid. VersoCit-2 is translucent and used for routine examination. VersoCit-2 cures at room temperature.

Safety Data Sheet

Please download a copy of the **Safety Data Sheet** from the Library section on www.struers.com and read carefully before use.

Access the Struers Library from the top menu on the website.

Please check the Library periodically for updates to the Safety Data information.

Mounting of Specimens

Measuring of Components

The ratio for VersoCit-2 Powder : VersoCit-2 Liquid is **15:10** by weight and **20:10** by volume. Weighing is recommended to obtain the best mounts and for reproducibility.

Important!

Do NOT pour excess liquid back into the bottle as this will significantly reduce the shelf-life.

Mixing of Components

- Measure out the amount of VersoCit-2 Liquid into a mixing cup (max. 40 g / 40 ml).
- Measure out the amount of VersoCit-2 Powder into a mixing cup (max. 60 g / 80 ml).
- Pour the powder into the liquid.
- Stir thoroughly for **at least 30 seconds** until the mixture is homogeneous, without whisking air into the mixture.

Selecting the Mounting Cup

All Struers mounting cups can be used. FlexiForm is recommended for easy unmounting of the mount. If possible, choose a mounting cup where the distance between the specimen and cup wall is approx. 5 mm/0.2".

Pre-treatment of specimens

To promote the best possible adhesion of VersoCit-2 to the specimen ensure that the specimen is thoroughly cleaned and dry prior to mounting. Each surface must be free from grease, moisture and other contaminants.

Mounting

- Place the cleaned specimen in the centre of the mounting cup.
- Mix the required amount of mixture.
- Pour the mixture carefully over the specimen.

Curing

VersoCit-2 is a room temperature curing system. For a 30 mm dia. mount, the curing time is about 10 minutes. We recommend keeping the mount in the mounting cup until the temperature is below 40°C/104°F.

After Curing

After removal from the mounting cup, the mount is ready for preparation.

Technical Data

	Components	Ratio by weight	Ratio by volume
VersoCit	VersoCit-2 Powder	15	20
	VersoCit-2 Liquid	10	10

NB For the best mount use weighing method.

Density

VersoCit-2 Powder: 1.16 g/ml

VersoCit-2 Liquid: 1.03 g/ml

Soluble in

VersoCit-2 Liquid is soluble in organic solvents.

Curing

Can be prepared (ground and polished) after:	10 min.
Peak temperature *:	100°C / 212°F
Pot life:	3 min.

*30 mm dia. mounting cup without specimen at 21°C / 70°F.

Properties after Curing

Hardness: 82 Shore D

Density: 1.2 g/ml

Chemical resistance: Most water diluted acids and bases.

Colouring

VersoCit-2 can be coloured with AcryDye.

Tips

- If the mixture is not miscible add a few drops of liquid.
- To avoid a sticky surface of the mount, cover the mounting cup during curing

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