# AbraPol-10

## Instruction Manual

Manual No.: 15097001

Date of Release 15.06.2005

*AbraPol-10*  
*Instruction Manual*

<b>Table of Contents</b>	<b>Page</b>
User's Guide .....	1
Reference Guide.....	31
Quick Reference Guide.....	63

---

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the Serial No. and Voltage on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

**Instruction Manuals:** Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

**Service Manuals:** Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2005.

**Struers A/S**  
Pederstrupvej 84  
DK-2750 Ballerup  
Denmark  
Telephone +45 44 600 800  
Fax +45 44 600 801

---

## **AbraPol-10**

### **Safety Precaution Sheet**

#### **To be read carefully before use**

1. The operator should be fully aware of the use of the machine according to the instruction manual.
2. The machine must be installed in compliance with local safety regulations.
3. The machine must be placed on a safe and stable support. The machine must be leveled by means of the adjustable legs.
4. The machine must be placed in an adequate working position.
5. Make sure that the actual voltage corresponds to the voltage stated on the side of the machine. The machine must be earthed.
6. Make sure that the water connections are without leaks. The main water supply should be turned off if you leave the machine unattended.
7. The machine must be disconnected from the mains prior to any service.
8. Make sure that the specimens in the specimen holder are securely fixed when operating.
9. Alcohol based consumables: follow the current safety rules for handling, mixing, filling, emptying and disposal of the alcohol-based liquids.  
Struers recommend the use of an external exhaust system.
10. Observe the current safety regulations for handling, mixing, filling, emptying and disposal of the additive for cooling fluid.
11. To achieve maximum safety and lifetime of the machine, use only original Struers consumables.

*AbraPol-10*  
*Instruction Manual*

- 12.** Ensure that the splash ring is in place when operating.
- 13.** Keep clear of the preparation disc when operating.
- 14.** If you observe malfunctions or hear unusual noises - stop the machine and call technical service.
- 15.** Make sure that the screws for the specimen holder have the right length and do not protrude.

---

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

---

# User's Guide

<b>Table of Contents</b>	<b>Page</b>
<b>1. Getting Started</b>	
Checking the Contents of the Packing Box .....	3
Recirculation Cooling Unit TRECI (accessory) .....	3
Additional Dosing Unit ABTDO (accessory) .....	3
Placing AbraPol-10 .....	3
Getting Acquainted with AbraPol-10 .....	4
Front of AbraPol-10 .....	4
Supplying Power .....	5
Direction of Rotation .....	5
Supplying Compressed Air .....	6
Connection to an External Exhaust System .....	6
Mounting a Recirculation Cooling Unit (accessory) .....	7
Placing Bottles in the Dosing Unit .....	8
Mounting an Additional Dosing Unit (accessory) .....	8
Mounting the Outlet Kit (accessory) .....	10
Mounting the Water Supply Kit (accessory) .....	12
Mounting the Disc Cooling Kit (accessory) .....	13
Mounting the Stock Removal Sensor (accessory) .....	14
Mounting the Water Level Sensor (accessory) .....	15
<b>2. Operation</b>	
Front Panel .....	16
Front Panel Controls .....	17
Software Settings .....	18
Setting the Language .....	19
Setup Bottle Configuration .....	20
Reading the Display .....	22
Sleep Mode .....	22
Changing/Editing Values .....	23
Numeric Values .....	23
Alphanumeric Values .....	24
Programming a Preparation Step .....	25
Mounting a Preparation Disc (300 or 350 mm) .....	25
Inserting the Specimen Holder .....	25
Adjusting the Specimen Holder Position .....	25
Starting the Preparation Process (Struers Method) .....	25
Stopping the Preparation Process .....	26

*AbraPol-10*  
*Instruction Manual*

Manual Functions..... 26  
    When Grinding with Water ..... 26  
    When Grinding/Polishing with Lubricant..... 26

**3. Maintenance**

Weekly..... 27  
    Coolant ..... 27  
General Cleaning..... 27  
    Cooling Tank..... 27  
    Painted Surfaces..... 27  
    Cleaning of Tubes ..... 28

## **1. Getting Started**

### **Checking the Contents of the Packing Box**

In the packing box you should find the following parts:

- 1 AbraPol-10 complete
- 1 Splash ring
- 1 Drain angle, 90°
- 1 Pipe angle connection ø50
- 1 Pipe angle connection ø32
- 1 Hose for compressed air
- 3 Hose clamps
- 1 Allen key 4 mm
- 2 Bottles, ½ l
- 2 Bottle lids
- 1 Bottle, 1 l
- 1 Set of Instruction Manuals

### *Recirculation Cooling Unit TRECI (accessory)*

- 1 TRECI complete
- 1 Drain angle, 45°
- 1 Funnel
- 2 Disposable plastic inserts

### *Additional Dosing Unit ABTDO (accessory)*


Three additional pumps for dosing of OP-Suspension, lubricants and suspensions.

### **Placing AbraPol-10**

AbraPol-10 should be placed on a plane and horizontal floor. The machine must be placed close to the power supply, compressed air, water mains and water outlet facilities.

**Getting Acquainted with  
AbraPol-10**


Take a moment to familiarise yourself with the location and names of the AbraPol-10 components.


*Front of AbraPol-10*

- |   |                    |
|---|--------------------|
| ① Specimen holder motor                         | ⑥ Turn/push button |
| ② Handle for adjusting specimen holder position | ⑦ Display          |
| ③ Quick coupling                                | ⑧ Front panel keys |
| ④ Double start buttons                          | ⑨ Nozzle block     |
| ⑤ Bottle unit                                   | ⑩ Emergency stop   |

## Supplying Power

- Remove the 4 Allen screws, on the right hand side of the machine, securing the electrical panel, and let the electrical panel rest on the two tabs.
- Lead the cable through the relief angle and connect the 3 phases and earth according to local regulations.
- In case of an external Recirculation Cooling Unit, connect the electrical cable from the pump, according to the diagram inside the panel.

### **IMPORTANT**

Check that the mains voltage corresponds to the voltage stated on the type plate on the side of the machine.

## *Direction of Rotation*

The preparation disc should rotate counter-clockwise.  
If the direction of rotation is clockwise, switch off AbraPol-10 and unplug the machine.

- Change two of the phases.
- Repeat the rotation check.

### **Supplying Compressed Air**

- Connect the compressed air supply with the inlet on the rear, left side of the machine by means of the air hose and the hose connection delivered with the machine.
- Fasten the air hose with a hose clamp.

The pressure supply should be 6-10 bar and should be supplied either from a central compressor, portable compressor with compressed air reservoir or compressed-air bottle. A capacity of 20 l/min at atmospheric pressure is sufficient. Please refer to the section on Technical Data for recommended air quality.


### **Connection to an External Exhaust System**

An exhaust system can be connected when using alcohol based suspensions or lubricants. Slot the 50 mm pipe angle connection into the hole on the left hand side of the splash guard and connect to the exhaust system. Recommended capacity for exhaust system: 180m<sup>3</sup>/h at 0mm water gauge.

## Mounting a Recirculation Cooling Unit (accessory)

Inlet hose      90° drain angle

### **Warning**

Disconnect the machine from the mains power supply before connecting or disconnecting the recirculation pump.

- Place the unit where you find it convenient, either inside the machine or outside, on the left.

### **Note!**

The Recirculation Cooling Unit cannot be positioned outside the machine when an outlet kit (ABTDR) is also used.

- Connect the inlet hose from the dosing arm to the pump of the Recirculation Cooling Unit and tighten the clamp. (If the cooling unit is placed outside the machine, pass the hose under the cabinet).
- Lead the electrical cable into the electrical panel and connect to the terminal points. (If the cooling unit is placed outside, lead through the angle on the panel).
- Slot the 90° drain angle into the pipe in the top of the reservoir lid. (Note: The 45° elbow pipe delivered with the Recirculation Cooling Unit is not used with this machine).
- Cut the drain hose to an appropriate length, and twist it into the hole located on top of the unit. (If the cooling unit is placed outside, lead the drain hose through the hole on the rear left of the machine).
- Check that there is a steady fall on the whole course of the outlet hose.
- Check that when the power is turned on, the pump rotates in the direction indicated by an arrow on the top of the pump. If not, change two of the phases on the pump motor cable.
- Replace the panel/ back plate.
- Place a disposable plastic insert in the tank and fold it over the rim.
- Fill the tank with 29.1 l water and 900 ml Struers Additive for Cooling Fluid. The water level should be 8-10 cm below the upper rim of the tank.
- Fit the tank lid and sieve.
- Replace the front cover.

### **IMPORTANT**

Too high a level of coolant in the tank might damage the pump. To avoid this, place the disposable insert so that the pre-punched hole is in front of the overflow aperture in the tank.

Always maintain the correct concentration of Struers Additive in the cooling water (percentage stated on the Additive container).

Remember to top up with Struers Additive each time you refill with water.

### Placing Bottles in the Dosing Unit

- Place the filled bottles in the dosing unit and connect the tubing.
- Enter the bottle details in the Bottle Configuration menu to make them available for preparation methods. See “Setup Bottle Configuration”.



The tubes can easily go through the bottle unit to reach larger containers placed on the floor e.g. Lubricants.

**Nb!**

OP-X can not be placed on the floor.


### Mounting an Additional Dosing Unit (accessory) ABTDO

- Disconnect main power.
- Open the doors to the bottle unit.

- Remove the rectangular cover plate ① with the 4 mm Allen key supplied.
- Remove the cover plate located on the left side of the machine with the supplied 4 mm Allen key supplied, and loosen the plug.
- Lead the reinforced water hose from the hole on the left, through the hole in the reinforcement beam on the inside of the bottle unit, and through the rectangular hole in the dosing unit.
- Connect the hose to the additional dosing unit, and tighten the clamp.
- Place the ½" - ¾" fitting in the hole (¾" on outside) and tighten the ½" nut.
- Connect the reinforced hose to the ½" thread.

- Connect the plug on the additional dosing unit to the plug located in the bottom of the dosing unit.

**NB**

Remember the earth connections.



- Replace the two cover plates.
- Connect the white water hose to the  $\frac{3}{4}$ " thread.
- Mount the straight end of the pressure hose onto the water inlet tube on the back of AbraPol-10:
  - Insert the filter gasket in the coupling nut with the flat side against the pressure hose.
  - Tighten the coupling nut completely.
- Mount the other end of the pressure hose on the water mains tap for cold water:
  - Mount the reduction ring with gasket on the water mains tap, if necessary.
  - Introduce the gasket and tighten the coupling nut completely.
- Connect the silicone tubes to the inlets/outlets of the dosing pumps.
- Check that the reinforced hose does not interfere with the belts.
- Reconnect main power and allocate the bottles.

**NB**

OP-X suspension can only be used in bottle 4.


**Mounting the Outlet Kit  
(accessory)  
ABTDR**

- Disconnect main power.
- Remove the lower front cover plate.
- Remove the Recirculation Cooling Unit, if there is one.
- Take the entire unit, and place it in position by letting the two Allen bolts slide into the slots located in the left rear side.
- Tighten the Allen bolts.

- Strip the steel spring from approx. 3 cm of the drain hose and cut. Bend the cut end of the steel spring inwards towards the centre of the hose. Mount the drain hose onto pipe **A** and clamp the stripped section using a hose clamp.
- Take another piece of drain hose and check that it is long enough to reach the Recirculation Cooling Unit. Mount onto pipe **B**. (Strip some steel spring from the hose and secure the hose as described earlier).
- Connect to the Recirculation Cooling Unit.
- Connect a piece of drain hose from the outlet at the rear (**C**) and lead to the drain.
- Reposition the Recirculation Cooling Unit.
- Make sure that the hoses do not interfere with the belts.



- Locate the electrical cable for the air solenoid valve (located on the inner left side of the machine) and remove the cap covering the stripped wires.
- Open the cover of the DIN terminal by removing the screw and connect the wires to the terminals.
- Replace the cover for the DIN terminal.

**Note:**

It is not important to which terminal the wires are connected, each wire can be connected to either terminal.



- Find the air supply hose and cut off the plug (stopper). Fix the cut end into the air supply inlet.
- Replace the front cover and reconnect main power.

**Mounting the Water Supply Kit  
(accessory)  
ABWAT**

- Disconnect main power.
- Remove the lower front cover plate.
- Remove the cover plate located on the left side of the machine with the supplied Allen key, and break loose the plug.
- Mount the water solenoid valve with the supplied Allen bolts and key according to the type plate.
- Connect the two blades (cable W23) to the solenoid.
- Connect the reinforced hose from the doser arm to the hose nipple, and tighten the clamp.
- Be careful that the hose does not interfere with the belts.
- Remount the cover plate.
- Connect the white hose to the solenoid valve.
- Mount the straight end of the pressure hose onto the water inlet tube on the back of AbraPol-10:
  - Insert the filter gasket in the coupling nut with the flat side against the pressure hose.
  - Tighten the coupling nut completely.
- Mount the other end of the pressure hose on the water mains tap for cold water:
  - Mount the reduction ring with gasket on the water mains tap, if necessary.
  - Introduce the gasket and tighten the coupling nut completely.
- Replace the front cover and reconnect the main power.

**Mounting the Disc Cooling Kit (accessory)**

- Disconnect main power.
- Remove the lower front cover plate.
- Remove the cover plate located on the left side of the machine with the supplied Allen key, and break loose the plug.
- Mount the water solenoid valve and the throttle valve with the Allen bolts supplied.
- Connect the two blades (cable W24) to the solenoid.
- Cooling is factory installed; simply connect the grey hose (ø4) from the nozzle to the single banjo located on the throttle valve.
- Remount the cover plate.
- Connect the white hose to the solenoid valve.
- Mount the straight end of the pressure hose onto the water inlet tube on the back of AbraPol-10:
  - Insert the filter gasket in the coupling nut with the flat side against the pressure hose.
  - Tighten the coupling nut completely.
- Mount the other end of the pressure hose on the water mains tap for cold water:
  - Mount the reduction ring with gasket on the water mains tap, if necessary.
  - Introduce the gasket and tighten the coupling nut completely.
- Replace the front cover and reconnect main power

**Adjusting Disc Cooling**

Should the rate of disc cooling require adjustment:

- Loosen the contra nut and turn the regulator screw to the required rate.
- Re-tighten the contra nut.


### Mounting the Stock Removal Sensor (accessory)

ABMEU

Plug


- Disconnect main power.
- Remove the lower front cover.
- Connect the plug on the removal sensor with the plug located in the recess to the left of the control panel.


Stock Removal Sensor


- Take the removal sensor and gently slide it up in the hole located in the bottom of the cast iron beam.
- Let the bracket slide over the factory installed support bracket.
- From the left, slide the two screws through the brackets, and tighten with a 10mm spanner.
- Reconnect main power and follow the instructions on the screen, using the tools supplied and the calibrating sticks.

**Mounting the Water Level  
Sensor (accessory)**  
PAMWA


- Switch the machine off.
- Pull out the Recirculation Cooling Unit.
- Remove the pump, drain hose and the top cover.
- Remove the red cover from the hole for the sensor or remove the existing Water Level Indicator.
- Remove the locking nut on the bottom of the sensor and slide off the black float.
- Mount the new Water Level Sensor on the Recirculation Cooling Unit cover and secure using the nut supplied (flat side of nut up).
- Replace the float (with the clear surface up) and secure the locking nut in place.
- On the right hand side of the machine, remove the 4 Allen screws securing the electrical panel, and let the electrical panel hang on the two tabs (as illustrated in Supplying Power).
- Remove the dummy jack plug located on the electrical panel by pressing down the red tab.
- Plug in the jack plug from the sensor.
- Reassemble the Recirculation Cooling Unit and slide it back into position.
- Replace the front cover.
- Switch the machine on.









## 2. Operation

### Front Panel

For more information also see page on Front Panel Controls.


**Front Panel Controls**

Name	Key	Function	Name	Key	Function
<b>1</b> FUNCTION KEY		Controls for various purposes. See the bottom line of the individual screens.	<b>7</b> Push/Turn Knob		Used for entering and changing steps and parameters. Combined cursor and enter key.
<b>2</b> FUNCTION KEY		Controls for various purposes. See the bottom line of the individual screens.	<b>8</b> DISC		Starts the rotation of the disc.
<b>3</b> FUNCTION KEY		Controls for various purposes. See the bottom line of the individual screens.	<b>9</b> Water		Activates the water flow. Remember to adjust the flow on the water tap.
<b>4</b> FUNCTION KEY		Controls for various purposes. See the bottom line of the individual screens.	<b>10</b> LUBRICANT		Manual dosing of lubricant
<b>5</b> Esc		Leaves the present menu or aborts functions/changes.	<b>11</b> ABRASIVE		Manual dosing of abrasive
<b>6</b> STOP		Stops the preparation process.	EMERGENCY STOP		- Push the red button to stop. - Pull the red button to release.
			MAIN SWITCH		The main switch is located on the right side of the machine.


## Software Settings

Switch on the power at the main switch located at the right hand side of the machine. The following display will appear briefly:


The display will then change to the last screen shown before AbraPol-10 was switched off, usually a preparation method. When switching on AbraPol-10 for the first time, the MAIN MENU display will appear. If the heading in the display is different, press Esc, until the MAIN MENU appears.

The MAIN MENU is the highest level in the menu structure. Pre-defined Struers methods, your own methods, manual functions or the configuration menu can then be selected.

Turn knob to select CONFIGURATION.


Push knob to activate the CONFIGURATION MENU.

*Setting the Language*

 Turn knob to select Language.

 Push knob to activate the Language pop-up menu.

 Turn knob to select the language you prefer.

 Push knob to accept the language.


The CONFIGURATION Menu now appears in the language you have chosen.

Press Esc to return to the MAIN MENU.


*Setup Bottle Configuration*

 Turn knob to select Bottle Configuration.

↓

 Push knob to activate EDIT BOTTLE CONFIGURATION.

↓

 Turn knob to select LUB/SUSP.

↓

 Push knob to activate pop-up menu.

↓

 Turn knob to toggle between Lubricant and Suspension.

↓


 Push knob to select Lubricant or Suspension.

Turn knob to move to TYPE.

Push knob to activate pop-up menu.

Turn knob to choose between the different options.

Push knob to select an option.

Repeat this procedure for all the bottles in use.


Press Esc twice to return to the Main Menu.

**NB**

The different suspensions and lubricants used must always be defined here so that AbraPol-10 can locate the correct suspension or lubricant.


**Reading the Display**

The display can show various kinds of information, for example parameters of a preparation method, or the ongoing process. An example of the screen for a preparation method is shown below:

- A Heading.
- B Inverted text: cursor position.
- C Help text.
- D Function key options.

During the process the screen could look as follows:

- A Time remaining in currently running step.
- B Preparation step currently running.

*Sleep Mode*

To increase the lifetime of the display, the backlight is automatically switched off if AbraPol-10 has not been used for 15 min. Press any key to re-activate the backlight.

**Please Note**

The sample screens in this Instruction Manual show a number of possible texts. The actual display screen may differ from the samples in the Manual.

## Changing/Editing Values

Depending on the type of value, there are two different ways of editing.

STRUERS METHODS can not be edited or changed in any way, this is only valid for USER METHODS and CONFIGURATION.

### Numeric Values

 Turn knob to select the value to be changed, e.g. Force:

↓

Step No. : 2		EDIT METHOD STEP	
Surface	: MD-Allegro		, 150 rpm
Grit/Grain	: User 1 Susp.		
Lubricant	: DP-Blue		P. dose: 6s
Step Mode	: Removal		
S. Mover	: 150 rpm		Direc.: ▷ ▷
<hr/>			
Stock Removal	:	MAIN	FINAL
		90 μm	110 μm
Force	:	<b>[100]</b> N	200 N
Abrasive Dosing	:	5	6
Lubricant Dosing	:	7	8
<hr/>			
F1	↗	F3	F4
Prev.	Next		

↓

 Push knob to edit the value.

↓

Two square brackets [ ] appear around the value.

↓

Step No. : 2		EDIT METHOD STEP	
Surface	: MD-Allegro		, 150 rpm
Grit/Grain	: User 1 Susp.		
Lubricant	: DP-Blue		P. dose: 6s
Step Mode	: Removal		
S. Mover	: 150 rpm		Direc.: ▷ ▷
<hr/>			
Stock Removal	:	MAIN	FINAL
		90 μm	110 μm
Force	:	<b>[100]</b> N	200 N
Abrasive Dosing	:	5	6
Lubricant Dosing	:	7	8
<hr/>			
F1	↗	F3	F4
Prev.	Next		

↓


 Turn knob to increase or decrease the numeric value.

↓

 Push knob to accept the new value. (Pressing Esc aborts the changes, preserving the original value.)

Alphanumeric Values

 Turn knob to select the alphanumeric value to be changed, e.g. Surface

↓  
 Push knob to edit the value.

↓ A pop-up menu appears.

↓  
 Turn knob to select the correct surface.

↓  
 Push knob to accept the new preparation surface and to return to the previous screen. (Pressing Esc aborts the changes, preserving the original setting.)

### **Programming a Preparation Step**

Different preparation steps can be programmed individually in regard of time, speed of disc, water, type and dosing of lubricant and suspension. These steps can be combined in a number of different ways so that different materials can be prepared without having to repeat the programming.

### **Mounting a Preparation Disc (300 or 350 mm)**

- Move the nozzle block to the side.
- Place the disc on the turntable and move it until the three pins engage with the holes in the turntable.

### **Inserting the Specimen Holder**

- Place the specimen holder under the quick coupling.
- Press and hold the black flange of the column down with your thumbs while guiding the pressure tap of the specimen holder into the coupling.
- Release the black flange.
- Turn the specimen holder until the three pins engage with the corresponding holes.

### **Adjusting the Specimen Holder Position**

The position of the specimen holder has to be adjusted correctly in relation to the preparation disc to get the best possible preparation results.

- Loosen the two handles placed on each side of the specimen holder motor.
- The specimen holder motor can now be moved manually to the sides. The eccentricity can be read on the scale placed on the specimen holder motor.
- When the correct position is found fasten the handles again.

### **Starting the Preparation Process (Struers Method)**

- Check that the splash ring is in place.
- Press ESC to get to Main Menu.
- Select Struers Methods in Main Menu by turning the knob.
- Push the knob for entering Struers Methods.
- Turn the knob to select a Struers Method.
- View Method by pushing the knob.
- Start the preparation by pressing the two double start buttons simultaneously.
- Keep the buttons pressed until the sample motor starts.

## **Stopping the Preparation Process**

When the time has elapsed the preparation disc will automatically stop and the specimen holder will return to its initial position. The machine is now ready for the next step.

- A preparation step can be stopped at any time during the process by pressing the key STOP ⊖.
- If a process has been stopped using the emergency stop button – you will get a message when pulling the button outwards again to release it. The sample mover will return to its initial position before starting again.

## **Manual Functions**

Manual operations are possible as the disc can run independently of the specimen holder.

The parameters for manual functions disc speed, lubricant and suspension can be altered in the menu Manual Functions, see Advanced Operations.

- Start the disc separately by pressing the key DISC ○. The disc will operate at low speed (150 rpm).

### *When Grinding with Water*

- Start the flow of water independently of a preparation step by pressing the key WATER ☒.

### *When Grinding/Polishing with Lubricant*

- To be able to give manual dosing of lubricant and suspension you should first choose a preparation step where required lubricant and suspension is selected.
- You can now press the key LUBRICANT ☒ and the key SUSPENSION ☒. The dosing will be at maximum for as long as the key is pressed.
- Stop the disc and water by pressing the key STOP ⊖ or by pressing the key DISC ○ and the key WATER ☒.

### **3. Maintenance**

#### **General Cleaning**

The machine and the Recirculation Cooling System should be kept as clean as possible in order to avoid contamination of the specimens.

#### *Cooling Tank*

Every time water is changed, the cooling tank should be cleaned carefully to remove all grinding waste. Cleaning will be easier if Struers plastic inserts are used.

***IMPORTANT***

Always maintain the correct concentration of Struers Additive in the cooling water (percentage stated on the container of the Additive). Remember to add Struers Additive each time you refill with water.

***IMPORTANT***

Do not use benzene or petroleum for cleaning when the additive is used.

#### *Painted Surfaces*

Painted surfaces and keyboard should be cleaned with a moist cloth and common household detergents.

***IMPORTANT***

Never use alcohol, acetone or similar solvents.

#### **Weekly Service**

##### *Checking the Recirculation Cooling Unit*

The cooling unit should be checked for cooling water after 8 hours use or at least every week. The unit must be refilled if the flushing pump cannot reach the cooling water or if the cooling water is too contaminated. Remember to add Struers Additive for Cooling Fluid: One part of Additive for 33 parts of water.

To check the concentration of additive, use a refractometer. Concentration = Brix value. The concentration of additive should always lie between 2.7 and 3.3 %. Add Struers Additive for Cooling Fluid if the concentration is too low.

#### **Monthly service**

##### *Replacing the Cooling Water*

- Replace the cooling water in the Recirculation Cooling Unit at least once a month.
- Place a disposable plastic insert (TREPO) in the tank and fold it over the edge.

***Please Note***

Too high a level of coolant in the tank might damage the pump.

- Fill the tank with 29.1 l of water and 900 ml Struers Additive.


*Cleaning of Tubes*

The cleaning procedure for CLEANING OF TUBES can be used weekly or when necessary.

Clean the tubes every time bottles are changed or replaced, thus avoiding lubricant/suspension left in tubes interfering with the preparation procedures.

Push knob to select MANUAL FUNCTIONS.

Push knob to activate CLEANING OF TUBES.

Cleaning of Selected Tubes

Turn knob to select a bottle.

CLEANING OF TUBES		
BOTTLE NO.	STATUS	SELECT
1, DP-Green	Clean	No
2, DP-Red	Used	Yes
3, 6Pm	Clean	No
4, None	Clean	No
5, None	Clean	No
6, None	Clean	No

Turn knob to select tube for cleaning  
Push knob to toggle parameter

**F1** Start cleaning of selected tubes      **F2** Select all used tubes

Push knob to toggle between YES and NO in the select column. One or more bottles can be selected.

Select YES and press **F1**, the cleaning process starts. After a while a pop-up message will appear.

Follow the instructions in the pop-up message.

```
[ TUBE CLEANING MESSAGE ]
The selected tube(s) are empty,
disconnect the tube(s) from the
bottles, and place the tube(s) in
a bowl filled with water. Press F1
when you are ready.
```

When tube cleaning is finished you get the following pop-up message.

```
[ TUBE CLEANING MESSAGE ]
Tube cleaning is finished, remove
tube(s) from bowl, and connect to
bottles again.
```

Cleaning of All Used Tubes

The tubes in use should be cleaned from time to time.

**F2** When in the CLEANING OF TUBES menu press **F2** to select all used tubes.

↓

CLEANING OF TUBES		
BOTTLE NO.	STATUS	SELECT
1, DP-Green	Clean	No
2, DP-Red	Used	Yes
3, 6PM	Used	Yes
4, None	Clean	No
5, None	Clean	No
6, None	Clean	No

Turn knob to select tube for cleaning  
Push knob to toggle parameter

**F1** Start cleaning of selected tubes      **F2** Select all used tubes

↓

**F1** Press **F1** to start the cleaning process and follow the pop-up instructions.

↓

[ TUBE CLEANING MESSAGE ]  
The selected tube(s) are empty, disconnect the tube(s) from the bottles, and place the tube(s) in a bowl filled with water. Press F1 when you are ready.

↓ When tube cleaning is finished you get the following pop-up message.

[ TUBE CLEANING MESSAGE ]  
Tube cleaning is finished, remove tube(s) from bowl, and connect to bottles again.

# Reference Guide


Table of Contents	Page
<b>1. Advanced Operations</b>	
Configuration Menu .....	33
Setup Bottle Configuration .....	33
Setup User Consumables .....	34
User Lubricant/Suspension.....	35
User Surface .....	36
Disc Cooling (optional) .....	37
Display Contrast.....	37
Language.....	37
<b>F1</b> , Default Value.....	37
Method Options.....	38
Struers Methods .....	38
Storing Capacity .....	38
Copying a Method.....	38
Inserting a Method .....	39
Renaming a Method .....	40
Name Editing Principles .....	41
Alphabetical Sorting of Names .....	42
Saving a Method .....	43
Deleting a Method .....	44
Step Options.....	45
Copying a Step.....	45
Inserting a Step .....	46
Deleting a Step .....	47
Method Parameters .....	48
Surface .....	48
Disc speed .....	48
Grit/grain .....	48
Lubricant.....	48
Pre-dosing .....	49
Step Mode .....	49
Sample Mover .....	49
Direction.....	49
Process Time; Stock Removal; Removal/Time .....	50
Force.....	50
Abrasive dosing.....	50
Lubricant dosing.....	50
300 rpm on Sample Mover and Disc.....	50

Manual Functions.....	51
Cleaning of Tubes .....	51
Manual Preparation .....	52
Disc Speed.....	52
Lubricant.....	52
Suspension .....	52
Starting Manual Preparation .....	52
<b>2. Accessories and Consumables .....</b>	<b>53</b>
<b>3. Struers Metalog Guide™.....</b>	<b>54</b>
<b>4. Trouble shooting</b>	
Error Messages .....	56
Messages .....	56
Errors .....	56
Fatal Errors .....	56
<b>5. Technical Data.....</b>	<b>60</b>
<b>6. Menu Overview.....</b>	<b>62</b>

## 1. Advanced Operations

### Configuration Menu


Press Esc until the Main Menu appears (a long beep can be heard).

Turn knob to select CONFIGURATION.

Push knob to activate the CONFIGURATION menu.

### *Setup Bottle Configuration*

Please look in the User Guide for detailed information on Bottle Configuration.


*Setup User Consumables*

The menus User Lub./Susp. Configuration and User Surface Configuration allow the user to add five suspensions, five lubricants and eight grinding/polishing surfaces.

A name can contain up to 16 characters and a name should only be used once.

When a lubricant or suspension has been named it is important to configure the bottle(s) i.e. a bottle number must be allocated for each lubricant and suspension name. For more information please see *Setup Bottle Configuration* in the User Guide.

The tree pictures below show the three pop-ups shown when in EDIT METHOD STEP to choose a suspension/lubricant and a surface.

The user defined names will always appear at the bottom of the pop-ups and are not sorted in any way.


User defined consumables can be selected in all user methods.

<b>IMPORTANT</b>
If consumable names change this will have influence on all user methods containing the consumable in question.


The following pages show how user consumables are created.

User Lubricant/Suspension

↓  
 Push knob to EDIT NAME OF USER CONSUMABLES.


↓  
 Push knob to select the name you want to change.


**F2** Edit name of lubricant/suspension and press **F2** to accept the new name.

For detailed information please see “Name Editing Principles”.


User Surface

↓  
 Push knob to EDIT SURFACE NAME AND STRUCTURE.

↓  
 Push knob to select the name you want to change.

**F2** Change name and press **F2** to accept new name.

For detailed information please see “Name Editing Principles”.

↓

 Turn knob to go the column  
STRUCTURE.

If the, user defined, surface is abrasive (contains  
abrasive grain) select Abrasive.  
If the surface does not contain abrasive grain. e.g. a  
polishing cloth select Non-abrasive.

 Push knob to chose one of the two  
options.

If ABRASIVE was selected, abrasive dosing will be set to OFF.

If NON-ABRASIVE was selected, the suspension to be used can  
be selected under Grit/Grain, a lubricant can be chosen under  
Lubricant, and the level for both Abrasive and Lubricant dosing  
can be adjusted.

*Disc Cooling (optional)*

Disc cooling can be adjusted in steps from 0 to 5. If cooling is  
turned off the setting will be 0 and the maximum setting is 5.

*Display Contrast*

The contrast settings of the display can be adjusted for maximum  
readability.

*Language*

The language can be set to English, German, French, Spanish or  
Japanese.

*F1, Default Value*

The factory setting of some parameters can be restored by  
pressing the function key **F1** when the appropriate value is  
highlighted.

**Method Options**

Methods can be adjusted to suit the particular requirements of the user.

*Struers Methods*

The AbraPol-10 software includes 20 Metalog Guide Methods, 10 for the 300 mm disc and 10 new methods which have been adapted to suit the new 350 mm disc.

As the Struers methods can not be changed or deleted, the first step will often be to copy a Struers method into the User Methods database and then adjust it to suit individual requirements.

*Storing Capacity*


The storing capacity for new/adapted methods is 250 units – a method name is one unit and a method step is also one unit.

Example: a method with four steps takes up five units  
a method with five steps takes up six units etc.


The typical method has approximately 3 to 5 steps and this allows room for 40 to 50 methods.

*Copying a Method*

In the Main Menu select STRUERS METHODS.

Push knob to activate the STRUERS METHODS menu

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↓ When asked for confirmation,

**F1** Press **F1**.

The method is now resident in the buffer. To save the method in a different location.


↓

*Inserting a Method*


**Esc** Press Esc to return to the Main Menu.

↓ In the Main Menu select **USER METHODS**.

**F2** Press **F2** to insert the copied method.

↓ When asked for confirmation,

**F2** Press **F2**.

*Renaming a Method*

The names of the methods in the USER METHODS database can be edited and changed. After copying a method from the Struers Methods the name can be changed to a name of your choice.

Select the method you want to rename.

Press **F4**: RENAME.

The following screen will appear:

Use **F1** or **F4** to select the character to be changed.

Use **F3** to delete the selected character.

Turn knob to move the cursor in the character set.

Push knob to insert new character.

For detailed information please see

“*Name Editing Principles*”.

**NB!**  
If a method name starts with “Copy of .....” the text “Copy of” can be deleted by pressing F3 once.

Name Editing Principles


Place the auxiliary cursor on the character you want to select, using **F1**: LEFT or **F4**: RIGHT. Turn the knob to move to the character cursor to select a character.

- Write the new name using the following keys:
  - F1** Moves the cursor to the left
  - F2** Accepts the new name
  - F3** Deletes one character in the text
  - F4** Moves the cursor to the right
  -  Moves character cursor when turned. A quick turn of the knob moves cursor to next/previous line.
  -  Enters characters (from character set) when pressed.

- A Auxiliary cursor
- B Character cursor
- C Character set

Write the new name using the name editing procedures.

*Alphabetical Sorting of Names*

Names are sorted alphabetically but EMPTY METHOD is always placed as the last method in the list.

The list can show eight methods at a time – if you want to move down in the list quickly a quick turn of the knob allow you to go to the next/previous page (eight methods at a time).


An arrow at the start and/or end of the list indicates more methods.

*Saving a Method*

While working with a user method, changes made in the database can be saved.

- As soon as a parameter is changed, **F4: SAVE** will be shown on the display and changes can be saved. Changes can be saved at both step and method level.

**F4** Press **F4: SAVE** to save the changes after all necessary parameters have been altered.

**Important**

When saving changes, the original method will be overwritten. To preserve the original method, a copy should be made and renamed. The new method can then be changed as required.  
For details, see *Copying a Method*.

If a new method is created from EMPTY METHOD the name will automatically change name from Empty Method to Unnamed Method 01 when saved.

The method can then be renamed and altered as required.

A new empty method is always created automatically.

*Deleting a Method*

If a method is no longer used it can be deleted.

Select the method you want to delete.

↓  
**F3** Press **F3:DELETE**.

↓ When asked for confirmation,  
**F3** Press **F3**.

The storing capacity is 250 units (see Storing Capacity).

If a lot of different user methods are stored it can be necessary to delete some user methods from time to time.

**NB!**  
Struers Methods can not be deleted.

## Step Options

### Copying a Step

Single steps can also be copied to customise preparation methods.

Select the method and the preparation step to be copied.

Method B 300 mm			VIEW	METHOD
SURFACE	GRIT GRAIN	LUBRICANT		REMOVAL TIME
1 SiC-Paper	#320	Water		1m 00s
2 MD-Largo	9µm	DP-Green		2m 20s
3 PoliSat1	6µm	DP-Green		3m 00s
4 PoliNat1	3µm	DP-Red		2m 00s
5 PoliFloc4	OP-U	None		1m 00s

Turn/Push knob to view step or select  
step for START

F1  F2 F3 F4  
F1 Copy

**F1** Press **F1**: Copy.  
A pop-up will appear.

↓

```
[ COPY ]
Are you sure you want
to copy step no. 4 to
buffer ?
Press F1 or F4
```

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*Inserting a Step*


My own method			VIEW METHOD
SURFACE	GRIT	LUBRICANT	REMOVAL TIME
1 SiC-Paper	#320	Water	2m 00s
2 MD-Largo	9µm	DP-Green	3m 00s
3 MD/DP-Oac	6µm	DP-Green	2m 30s
4 MD/DP-Mol	3µm	DP-Red	2m 00s
5 MD/DP-Chem	OP-U	None	1m 00s

↓  
 Turn/Push knob to edit step or select step for START

F1 Copy    F2 Insert    F3 Delete    F4

↓ Select the method and step where the copied step is to be inserted.

**F2** Press **F2**: Insert.  
 A pop-up will appear.

When asked for confirmation,

↓  
**F2** Press **F2**: Insert.  
 The pasted step is now inserted.

↓  
**F4** Press **F4**: Save.  
 The changed method is saved.

### *Deleting a Step*

When a step is inserted it is added to the method.  
Replacing a step therefore means inserting a step and afterwards deleting a step.

Turn knob to move to step you want to delete.

**F3** Press **F3** Delete.  
A pop-up will appear.

When asked for confirmation,


**F3** Press **F3**  
The step will be deleted.

**F4** Press **F4** Save  
To save the changed method.


**Method Parameters**

For every method step a series of parameters can be changed:

*Surface*

Pushing the knob allows you to select between all the surfaces in the AbraPol-10 database:

*Disc speed*

The speed can be set to 150 or 300 rpm.


*Grit/grain*

Pushing the knob allows you to select between the different grain sizes available:

*Lubricant*

Pushing the knob allows you to select between the different lubricants available:

*Pre-dosing*

Pushing the knob allows you to adjust the pre-dosing between 0 and 30 seconds.

*Step Mode*

Pushing the knob  $\leftarrow \sqcap$  allows you to select between the 3 different modes:

*Time*

The preparation time can be specified in both a main and a final phase. Usually the main phase is used to remove a large amount of material, by using high forces and longer times, whereas the final phase is used to improve the scratch pattern by applying lower forces for a short while at the end of the preparation step.

*Removal*

The preparation removal can be specified in both a main and a final phase. Using only one phase gives better measurement accuracy. The stock removal sensor must be installed to use this parameter.

*Removal/Time*

A combination of the above. First a specified amount of material is removed using a relatively high force. Thereafter the specimens are run for a short while at a low force to improve the scratch pattern. The stock removal sensor must be installed to use this parameter.

*Sample Mover*

Pushing the knob allows you to select between 150-300 rpm.

*Direction*

Turning the knob allows you to toggle between co- or counter rotation. Push the knob to choose

The direction of rotation of the specimen mover can be set to either co- (►►) or counter (►◄) rotation. Co-rotation will give the best preparation result and the most uniform polishing. Counter-rotation can be used when polishing using oxide polishing suspensions to keep the suspension on the polishing cloth.

*Process Time; Stock Removal;  
Removal/Time*

Depending on the mode selected either the process time; the stock removal or a combination of both can be adjusted.  
The maximum preparation time is 20 minutes in each phase.  
The maximum stock removal is 1000 µm in each phase.

*Force*

Pushing the knob allows you to adjust the force between 50 and 700 Newton.

*Abrasive dosing*

Pushing the knob allows you to adjust the abrasive dosing level between 0 and 20.

**IMPORTANT**

If a grinding or polishing disc with a surface containing abrasive is selected, dosing is automatically set to OFF.

*Lubricant dosing*

Pushing the knob allows you to adjust the lubricant dosing level between 0 and 20. It is possible to toggle between ON and OF if water is selected.

*300 rpm on Sample Mover  
and Disc*

When using 300 rpm on the sample mover and 300 rpm on the disc please note the following recommendations:



- Use disc cooling.
- Pre-dose sufficiently when working with new preparation surfaces.
- Use a 200 mm sample holder on a 350 mm surface to distribute the lubricant and abrasive adequately.
- Use DP-Lubricant, Green or DP-Lubricant, Red.  
DP-Lubricant, Blue is alcohol-based and will evaporate too quickly.

## Manual Functions

Several manual functions are included in the software of AbraPol-10. Press Esc until you reach the Main Menu.

 Push knob to select MANUAL FUNCTIONS.


### *Cleaning of Tubes*

Cleaning of tubes is described in depth earlier in the manual in section 3 Maintenance.

*Manual Preparation*

Manual functions are possible as the disc can run independently of the specimen holder.


*Disc Speed*

The speed can be set to 150 or 300 rpm.

*Lubricant*

The different lubricants can be chosen here.


*Dosing Level*

Pushing the knob and then turning it allows you to adjust the dosing level between 0 and 20.

*Suspension*

The different suspensions can be chosen here.

*Dosing Level*

Pushing the knob and then turning it allows you to adjust the dosing level between 0 and 20.

*Starting Manual Preparation*

After defining all parameters the manual preparation can be started by pressing **F1**.

**NB!**

The suspension and lubricant selected in the Manual Preparation menu must correspond to the available bottles – see also Bottle Configuration.

## 2. Accessories and Consumables

<b>Specification</b>	<b>Code</b>
<i>Disc for Magnetic Fixation</i> MD-Disc, 350 dia. Gross weight 4.6 kg.	DEMFI
<i>Additional Dosing Unit</i> Three additional pumps for dosing of OP-Suspension, lubricants and suspensions - 3 pumps	ABTDO
<i>Outlet Kit</i> Shift valve for recirculation cooling unit (TRECI) for drain	ABTDR
<i>Water Supply Kit</i> For control of tap water supply. For AbraPol-10 without recirculation cooling unit	ABWAT
<i>Disc Cooling</i> Disc cooling for AbraPol-10. Only to be used with tap water. Require water supply kit (ABWAT)	ABTCO
<i>Accurate Material Removal Unit</i> Stops process after required amount of material has been removed (Stock Removal Sensor)	ABPEU
<i>Recirculation Cooling Unit</i> Recirculation cooling unit, complete, 3-phased	TRECI
<i>Inserts for Recirculation Cooling Unit</i> 20 disposable plastic inserts	TREPO
<i>Water Level Indicator</i> Sensor indicates the level of water in the recirculation cooling unit (TRECI)	PAMWA
<i>Additive for Cooling Fluid</i> 1 litre 5 litre	ADDUN ADDFI

***Remember...***

Struers offers a comprehensive range of consumables for grinding and polishing. Please ask for separate leaflets.

### **3. Struers Metalog Guide™**

In Struers Metalog Guide™ you will find a detailed description of grinding/polishing methods for automated mechanical specimen preparation.

Struers Metalog Guide™ offers preparation methods for the most common materials, based on a simple analysis of two key properties: hardness and ductility. Finding the right method is easy, including choice of consumables. Always consult Struers Metalog Guide™ for the correct preparation method for the actual specimens.

Struers Metalog Guide™ contains 6 useful chapters:

Metalogram: a quick and safe guide to the right preparation method.

Metalog Methods: a complete catalogue of preparation methods, based on Struers' vast store of materialographic experience, and employing Struers' range of consumables.

Preparation Philosophy: the basics of modern specimen preparation, seen from a professional point of view.

Metalog Process: the materialographic preparation process from start to finish, logically explained.

Metalog Master: a combined trouble-shooting guide and supply of in-depth information on the processes of mechanical preparation, including an expert system for the solving of preparation problems.

Consumables Specification: quick access to the relevant consumables for the chosen preparation methods.

A complete guide to materialographic specimen preparation.  
Contact your local dealer for a free copy of Metalog Guide™.

## 4. Trouble shooting

<b>Error</b>	<b>Explanation</b>	<b>Action</b>
Machine Problems		
AbraPol-10 is switched on but the display is blank	The backlight of the display has been switched off. AbraPol-10 is in Sleep Mode.	Press any key to re-activate the backlight.
Display text not clear	The display is sensitive to temperature changes.	Change the display contrast in the CONFIGURATION Menu.
Display defective	Defect in the circuit of the display or defect in the CPU circuit.	Call a technician.
Disc rotates in the wrong direction.	Wrong connection to the power supply.	Switch two of the phases.
Water not draining away.	Drain hose squeezed.	Straighten the hose.
	Drain hose clogged.	Clean the hose.
	Drain hose does not slope downward.	Adjust the hose to an even slope.

## **Error Messages**

Error messages are divided into three classes:

Messages

Errors

Fatal Errors

### *Messages*

Messages are intended to inform the operator of the machine's progress and advise about minor operational errors.

### *Errors*

Errors must be rectified before process can continue.

### *Fatal Errors*

In case of Fatal Errors, process cannot continue before an authorised technician has rectified the error. Turn off the unit at the main switch immediately. Do not attempt to operate unit before technician has rectified problem.

<b>Message</b>	<b>Explanation</b>	<b>Action</b>
Process in Progress	Start is not possible as other process in progress.	
Process already stopping	Appears if Stop is pressed repeatedly.	
Water level in tank too low	Water level in recirculation tank is too low.	Refill circulation water if needed.
Process is stopping	Stop has been pressed.	
Process is stopping, await flushing please	The process is stopping but flushing removes residues of OP-X.	After waiting a moment you can continue.
Emergency Stop is active	The emergency button was pressed and the preparation was interrupted.	Pull the button outwards to release it and follow the instructions in the pop-up text.
Database is full !	The database storing capacity has been reached.  For details please see Method Options.	Delete one or more of the methods – this will free space to store new methods.  Nb! it is not possible to delete Struers Methods.
10 steps is maximum !	The maximum of 10 steps has been reached.	Two or more methods can be used alternately.
2-hand-start released too early	The two start buttons have not been pressed long enough to start the process.	Keep the start buttons pressed until the sample motor starts.
Process is finished	The process is finished.	You can now start another process.

*AbraPol-10*  
*Instruction Manual*

<b>Message</b>	<b>Explanation</b>	<b>Action</b>
START denied method not selected	A method has not been select.	Select a method and press start.
Manual dosing denied method step not selected	Manual dosing has been denied because step has not been selected.	Choose a preparation step that contains the suspension or lubricant in question.
The following lubricant is not configured: " xx "	The lubricant chosen for a certain preparation has not been configured and the preparation can not proceed.	In the CONFIGURATION MENU it is possible change the set-op for the different lubricants. (See Setup Bottle Configuration)
You can not dose with water	The manual lubricant key has been pressed in a step where water was selected.	Use the  key for manual dosing of water.
Process time or stock removal not specified (zero)	The step can not start because neither time nor removal have been specified.	In EDIT METOD STEP specify time or stock removal.
START denied, selected method or step is Empty	It is not possible to start a preparation process from an empty method or step.	Chose a method or a step which is not empty before pressing start.
Tube(s) not selected for cleaning push knob to select tube(s)"	The tube has not been selected and can therefore not be cleaned.	Push knob to toggle between yes and no in the select column – yes must be selected if you want to clean the tube.
The following abrasive is not configured: " xx "	The abrasive chosen for a certain preparation has not been configured and the preparation can not proceed.	In the CONFIGURATION MENU it is possible change the set-op for the different abrasives. (See Setup Bottle Configuration)
Method name already in use, please rename the method	A method name can only be used once.	Use a different name.

*AbraPol-10*  
*Instruction Manual*

<b>Message</b>	<b>Explanation</b>	<b>Action</b>
Disc motor overloaded	The disc motor was overloaded and therefore overheated.	Wait till the motor has cooled off and reduce the force before restarting the process.
Sample motor overloaded	The sample motor was overloaded and therefore overheated.	Wait till the motor has cooled off and reduce the force before restarting the process.
Pump motor overloaded	Water hose blocked or too much swarf in coolant.	Replace coolant. Inspect pump.
Main supply voltage too low	Supply voltage insufficient.	Await normalisation of power grid.
Removal rate too low, phase time exceeded	It has not been possible to grind/polish the samples to the expected level within 15 minutes.	Select a different surface or/and increase the force.
Sample holder unable to move, upwards after process	It was not possible to move the sample holder upwards after finishing a process.	The pressure in the compressed air system is too low.
		Force regulation error – contact Struers service technician.
Force regulation error	It is not possible to reach the selected force.	The pressure in the compressed air system is too low.
		Force regulation error – contact Struers service technician.
Air pressure too low	Missing or incomplete air supply.	Check that the air hoses are tight and properly clamped.
Stock removal sensor not installed	Stock removal sensor has not been installed or has not been calibrated.	If stock removal sensor (accessory) is not installed this must be done before the stock removal mode can be used.
		If the stock removal sensor is installed it is not calibrated correctly, please look in the manual for the stock removal sensor.

*AbraPol-10*  
*Instruction Manual*

<b>Fatal Errors</b>	<b>Explanation</b>	<b>Action</b>
15 V DC supply in pcb missing	Fault in the internal power supply.	Turn off unit at Main Switch. Contact Struers service technician.
24 V DC supply in pcb missing	Fault in the internal power supply.	
No serial communication	No contact between machine system and operating panel.	
Program version mismatch	Conflict between software in machine control system and operating panel.	


## 5. Technical Data

Subject		Specifications	
		Metric/International	US
Disc	Rotational speed	150/300 rpm	150/300 rpm
	Size	300/350 mm	11.8"/13.7"
	Power Consumption	1.5 kW	2 Hp
Specimens	Speed	166/333 rpm	
	Direction	CW/CCW	
	Force	0-700 N	0-150 lbf
	Motor power consumption	0.37 kW	0.5 Hp
Dosing Unit	Level	0-20	
	No. of pumps	3	
Cooling water Recirculation Cooling Unit	Power supply	3 (3L+PE)	
	Power Consumption	140 W	0.2 Hp
	Width	500 mm	
	Depth	400 mm	
	Height	410 mm	
	Weight	7.2 kg	
Software and electronics	LC Display	240x128 pixels	
	Controls	Touch pads/Push-turn knob	
	Memory	EPROM/RAM/NV-RAM	
Environment	Safety standards	98/37/EØF EN292-1:1991; EN292-2:1991+A1; EN418:1992; EN574:1996; EN953:1997; EN954-1:1997; EN1050:1997; NFPA70:1996; UL508; NFPA79:1997 FCC part15, subpart	
	EMC	EN50081-1:1992; EN50082-1:1997	
	Noise level(idle 300/300)	77 dbA	
	Surrounding temperatur	5-40°C	
	Humidity	Max. 95%RH	

*AbraPol-10*  
*Instruction Manual*

Subject		Specifications	
		Metric/International	US
Supply Voltage	Power consumption	2 kW	
	No. of phases	3 (3L+PE)	
	Output, main motor	1.5 kW	
	Voltage/frequency:		
	3 x 200 V / 50Hz 3 x 200 - 210 V / 60Hz CSA 3 x 220 - 230 V / 50 Hz 3 x 220 - 240 V / 60 Hz 3 x 380 - 415 V / 50 Hz 3 x 380 - 415 V / 60 Hz 3 x 460 - 480 V / 60 Hz CSA		
Water mains	Water supply	max. 10 bar	
Compressed air	Compressed air supply Recommended quality, ISO 8573-1, class 5.6.4 for normal use	6-10 bar	
Fittings	Splash ring		
Accessories (opt.)	Water mains		
	Disc cooling(water mains)		
	Shift valve incl. external drain		
	3 dosing pumps incl. OP-X unit		
	Guard		
	Stock removal sensor		
Dimensions and weight	Width	670 mm	26.5"
	Depth	930 mm	36.6"
	Height	1500 mm	59"
	Weight	250 kg	

## 6. Menu Overview

# Quick Reference Guide

## Description

AbraPol-10 is an automated machine ideal for laboratories or workplaces carrying out quality control testing on a large number of specimens. AbraPol-10 has the unique possibility of preparation on a 350 mm disc. Furthermore, the Speed of specimen holder and speed of disc can be set to 300 rpm, resulting in reduced grinding and polishing time.

## Main Controls

### *Main Switch*

The main switch is located at the right side of the machine.

### *Emergency Stop*

Red push button on the front of the machine. Stops all ongoing operations. Release the emergency stop by pulling the button.

### *Double Start Buttons*

Start the actual preparation by pressing the two buttons simultaneously.

## Front Panel Controls

### *Manual Functions*

Start the Rotation of the disc

Switch the Water on

Manual dosing of Lubricant

Manual dosing of Abrasive

### *Function Keys*

Controls for various purposes.  
See the bottom line of the individual screens.

### *Stop and ESC*


Stops the grinding or polishing process

Leaves the present menu or aborts Edit functions/changes

### *Turn / Push Knob*

The Knob is used for entering and changing steps and parameters.

Combined cursor and enter key.

### **Running a Struers Method**

- Press ESC to get up to Main Menu.
- Select Struers Methods in Main Menu, by turning the knob.
- Push the knob for entering Struers Methods.
- Turn the knob to select the chosen Struers Method.
- View Method by pushing the knob.
- Start the preparation by pressing the two double start buttons simultaneously.

### *Editing process*

- For editing the process, push knob during the process.

### **Copying a Preparation Method**

- Turn the knob to select the preparation method you want to copy.
- Press **F1** for copying a Preparation Method
- Press **F1** to accept copying
- Enter the screen USER METHODS press **F2** to insert the method

### **Editing Names**

- Select the method name you want to edit.
- Edit name by pressing **F4**.
- Use knob to select/insert character.
- Move the cursor left or right by pressing **F1** or **F4**.
- Press **F2** to accept the new name.

### **Modifying a Preparation Method**

- Select User Methods in Main Menu, by turning the knob.
- Push the knob for entering User Methods.
- Select the Preparation Method and the step you want to change.
- Turn the knob to the parameter you want to change.
- Push the knob and select the new parameter.
- Save changes by pressing **F4**.

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*AbraPol-10*  
*Spare Parts and Diagrams*

---

**Always state *Serial No* and *Voltage/frequency*  
if you have technical questions or when ordering spare parts.**

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Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to change without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2005.

**Struers A/S**  
Pederstrupvej 84  
DK-2750 Ballerup  
Denmark  
Telephone +45 44 600 800  
Telefax: +45 44 600 801

---

## **Spare Parts and Diagrams**

### **Table of contents**

### **Drawing**

## **AbraPol-10**

### *Drawings*

AbraPol-10, complete.....	15090001I
Windows for doser module .....	15090004
Turntable with quad rings .....	15090005
Quick-release coupling, complete.....	15090007
Plate for electrical connection, complete .....	15090015C
Stand, complete.....	15090010D
Dosing arm, complete .....	15090055B
Plate with process print and pneumatic.....	15090031D
Front plate, complete.....	15090060B
Air filter/regulator unit, complete.....	15090032A
Contactator box, complete .....	15090033C
1L Bottle for OPS, complete.....	15090092B
1L Bottle for lubricant/suspension, complete.....	15090091B
Cabinet for drain tray.....	15090050D
Dosing module.....	15090051A
Pump, complete.....	15090053
Column guide, complete .....	15090020B
Main motor, complete.....	15090021A
Bearing unit sub assembly.....	15090040A
Sample mover motor, complete.....	15090045A
Outlet kit, complete .....	15090084C
Outlet kit for cooling unit/drain.....	15090082B
Outlet kit, assembled.....	15090094C
Additional dosing unit .....	15090085C
Accurate material removal unit.....	15090083F
Disc cooling .....	15090087B
Base unit assembly.....	15090023A

Some of the drawings may contain position numbers  
not used in connection with this manual.

*AbraPol-10*  
*Spare Parts and Diagrams*

*Diagrams*

Complete Block Diagram.....	15093050A
Contactor plate, layout + wiring (4 pages) .....	15093110B
Control voltage diagram, CPU-board, A1 (3 pages).....	15093120D
Main supply circuit.....	15093100C
Functional diagram, safety circuit.....	15093130D
Control panel.....	15093200D
Machine control board, A2 (8 pages).....	15093210E
Removal Measurement board .....	15093220D
Wiring diagram (5 pages).....	15093450F
Transformer connections (3 pages).....	15093452A

*AbraPol-10*  
*Spare Parts and Diagrams*

<b>Drawing</b>	<b>Spare Part list for AbraPol-10</b>	
	<b>Pos.</b>	<b>Cat no.</b>
15090001	<b>AbraPol-10, complete</b>	
	20	Cabinet for drain tray, complete 15090050
	50	Drain tray, complete R5090057
	70	Column guide, complete 15090020
	110	Rubber bushing 15090690
	140	Kipp-handle M10 w. parts RGH30020
	170	Sliding adjusting disc 14220368
	180	M2 Sample Mover motor, complete 15090045
	183	Safety switch 1NO+1NC (4R) RSS00008
	190	12.00A Fuse, midget 10x38 15092921
	190	8.00A Fuse, midget 10x38 15092921
	195	1.00A Fuse, midget 10x38 15092921
	196	2.00A Fuse, midget 10x38 15092921
	220	Process print + pneumatic mont. 15090031
	225	Contacto box, complete 15090033
	250	Plate for electrical connection, complete 15090015
	280	Front panel, complete 15090060
	390	Turntable w. quad rings (ø350 mm) 15090005
	405	Dosing arm, complete 15090055
	425	Lifting iron bar with screws R5494100
	450	Right/left window for doser, complete 15090004
	475-477	Screw & key for quick-release coupling, complete 15092920
	470	Quick-release coupling, complete 15090007
	540	Splash ring 15090900
	601	Bottle ½ l for DP-suspension 71000119
	602	Lid, complete for ½ l R4600603
	603	Bottle 1L for lubricant/suspension, complete 15090091
	604	Bottle 1L for OPS, complete 15090092
	605	Non-return valve elbow, complete 12600575
	606	Pressure hose PVC ¼", armoured 2NU12403
	607	Hose clamp S12/9 2NS11209
	608	Hose clamp S17/9 2NS11709
	609	Hose nipple 13688045
	611	Allen key 4, hardened DIN 911 2GR00040
	612	Fork wrench NV13 2GR00213
	613	Fork wrench NV24 DIN 894 2GR00224
	614	ABTDR Outlet kit 05096901
	615	ABTDO Additional dosing unit 05096902
	616	ABMEU Accurate material removal unit 05096903
	617	ABTCO disc cooling 05096904
	618	ABWAT water supply kit 05096905
	619	PAMWA water level sensor 03756912

*AbraPol-10*  
*Spare Parts and Diagrams*

<b>Drawing</b>	<b>Pos.</b>	<b>Spare Part list for AbraPol-10</b>	<b>Cat no.</b>
15090001		<b>AbraPol-10, complete</b>	
	621	TRECI Cooling unit 3X415-480V/50-60Hz	02536154
	621	TRECI Cooling unit 3X345-415V/50-60Hz	02536149
	621	TRECI Cooling unit 3X200-240V/50-60Hz	02536135
15090004		<b>Right/left window for doser, complete</b>	
	10	Left window for dosing module	15092901
	10	Right window for dosing module	15092902
15090005		<b>Turntable w. quad rings (ø350 mm)</b>	
	10	Turntable w. quad rings (ø350 mm)	15090005
	20	Shell w. screws for bearing disc	R5098001
	50	Quad-ring 4326-366Y	2IQ04326
	60	Quad-ring 4441-366Y	2IQ04441
	70	Quad-ring 4450-366Y	2IQ04450
15090007		<b>Quick-release coupling, complete</b>	
	5	Quick-release coupling, complete	260MP027
	6	Screw & key for quick-release coupling, complete	15092920
	30	Pressure spring ø4.3x0.7	12600718
	40	Driving pin	12600717
	60	Key /fitted key, type A	12600723
	65	Ball KU 5.556 (ø7/32"+D145) RS	2BA00055
15090015		<b>Plate for electric connection, complete</b>	
	10	Plate for electric connection, complete	R5090015
	215	Terminal strip SAK 6/35 PA	2XK41012
	216	Ground terminal	2XR42617
15090010		<b>Stand, complete</b>	
	20	Machine leg M75-16-070	2GB40075
	150	Air filter/regulator unit, complete	15090032
15090055		<b>Dosing arm, complete</b>	
	10	Dosing arm, complete	15090055
	30	INA-Sealing ring SD 14X20X3	2II21420
	100	Silicone hose 7/ø10	2NU19207
	110	Water pipe with spray nozzle	R5090180
	150	Nozzle block	R5090140
	180	Tygon hose ENFT 21 ø2.06	2NU91221
	220	Silicone hose ø4/ø6	2NU11454
	270	Grifl-hose k. arm PVC ½"-ø12.5	2NU29316
	290	Spring for dosing arm	15090151

*AbraPol-10*  
*Spare Parts and Diagrams*

<b>Drawing</b>	<b>Spare Part list for AbraPol-10</b>		
	<b>Pos.</b>	<b>Cat no.</b>	
15090031		<b>Plate with process print and pneumatic, complete</b>	
	20	PCB A2 process print ABRAPOL-10	15093002
	21	Set of Proms, latest version	15092900
	22	2.5AT Fuse, sub-miniature, TR5	2FU32900
	23	4.0AT Fuse, sub-miniature, TR4	2FU32999
	60	Pressure governor, 5-8.5 bar 1/4"	2YR00001
	70	Gasket, PVC 1/8"	2IF00011
	100	Gasket, PVC 1/4"	2IF00012
	150	3/2 solenoid valve 24V DC NORM.Å.	2YM10120
	190	Gasket, PVC M5	2IF00010
	220	3/2 solenoid valve 24VDC	2YM10124
	240	Angle for quick-coupling, rotatable ø5-M5	2NF10081
	270	Gasket, PVC M5	2IF00010
	280	Clamping ring	14590110
	285	PU hose PU-3 blue ø2.9/ø4.3	2NU14031
	290	Nylon hose, SUPERFLEX ø5/ø3.2	2NU12445
	310	Neoprene nipple ø36/ø47/ø54-2.5	2GK90457
	15090060		<b>Front plate, complete</b>
		10	Front plate with foil and glass
60		Display with plug	15343045
80		PCB A1 control panel ABRAPOL-10	15093001
81		Set of Proms, latest version	15092900
110		Light diode PCB	15022396
150		2 channel opt. encoder	2HR12411
160		Turn/push-button	R5090600
15090032		<b>Air filter/regulator unit, complete</b>	
	40	Air filter, air regulator EAW3000-F02D-6	2YF00005
	90	Gasket, PVC 1/4"	2IF00012
	110	Gasket, PVC 1/8"	2IF00011
	140	PVC-Hose, clear 13/32"-ø10	2NU19313
	150	Nylon tube, SUPERFLEX ø5/ø3.2	2NU12445
180	Gasket, PVC 1/8"	2IF00011	

*AbraPol-10*  
*Spare Parts and Diagrams*

<b>Drawing</b>	<b>Spare Part list for AbraPol-10</b>			
	<b>Pos.</b>	<b>Cat no.</b>		
15090033	<b>Contactor box, complete</b>			
	20	Safety relay two-hand 24V	2KS10020	
	30	Contactor CI 12 37H0032/13	2KM10232	
	40	Auxiliary BL.CB-NO 037H0111	2KH00111	
	50	Mechanical blocking f. C109-25	2KH10100	
	65	12.00A Fuse, midget 10x38 (F1-F3)	2FC10120	
	65	8.00A Fuse, midget 10x38 (F1-F3)	2FC10080	
	70	1.00A Fuse, midget 10x38 (F4-F8)	2FC10010	
	110	Contactor CA4-5-0, 24V-50/60HZ	2KM04501	
	120	Auxiliary contact BL.CBM-40 37H3149	2KH03149	
	125	Auxiliary contact BL.CBM-22 37H3152	2KH03152	
	130	Mechanical blocking CM 4	2KH10021	
	160	TRAFO 200-460V/24V+24V/200VA	2MT72034	
	15090092	<b>Bottle, 1L, for OP-S, complete</b>		
		10	Bottle, 1L, for OP-S, complete	15090092
		20	PVC-Hose, clear 5/32"-ø4	2NU19305
30		Lid for OP-S bottle, complete	15092917	
40		O-ring 42.86-3.53 72 NBR 872	2IO35355	
50		Neoprene nipple ø5/ ø7/ø11-1.5	2GK90108	
15090091	<b>Bottle, 1L, for Lubricant/Suspension, complete</b>			
	10	Bottle, 1L, for Lubricant/Suspension, complete	15090091	
	30	Lid for LUB/SUSP bottle, complete	15092918	
	40	O-Ring 42.86-3.53 72 NBR 872	2IO35355	
	50	Neoprene nipple ø2/ø6.5/ø11-1.5	2GK90102	
15090050	<b>Cabinet for drain tray, complete</b>			
	20	Dosing module, complete	15090051	
	80	Emergency stop, complete	260MP007	
	110	Two-hands release push-button, complete	15092903	
	180	Magnetic catch, black	RGL30035	
15090051	<b>Dosing module, complete</b>			
	20	BOSCH motor w. print + axle extension	14600550	
	90	Tygon hose ENFT 21 ø2.06	2NU91221	
	140	Silicone hose ø2/ø4	2NU11452	

*AbraPol-10*  
*Spare Parts and Diagrams*

<b>Drawing</b>	<b>Spare Part list for AbraPol-10</b>		<b>Cat no.</b>
	<b>Pos.</b>		
15090053		<b>Pump, complete</b>	
	20-70	Pump roller, complete	14600601
	50	Pump, assembled	15090053
	100	Lid for pump	14600602
	110	Locking spring	14600024
15090020		<b>Column guide, complete</b>	
	80	Cylinder EDCQ 2B 50-OD	2YC50101
	81	REED-Contact D-A73L	2KR30177
	120	Gasket, PVC 1/8"	2IF00011
	130	Gasket, PVC 1/4"	2IF00012
	160	Coupling FLEXO FK-M10	2JH50010
	220	INA-Sealing ring G 40x52x5	2II04052
	240	Buffer AA-315 ø30x15 m. M8x10	2GS10315
	400	Bearing unit sub assembly	15090040
	440	M1 Main motor, complete	15090021
	450	Narrow V-belt pulley XPZ/3V 1320mm	2JD01320
	470	50 Hz V-belt SPZ-2 ø315	2JE00315
	470	60 Hz V-belt SPZ-2 ø355	2JE00355
	480	Adapter 2012/ø35	2JE92035
15090021		<b>M1 Main motor, complete</b>	
	10	3x200V/50Hz 1,5/0,75kW-1440/710 rpm	15092904
	10	3x220-230V/50Hz 1,5/0,75kW-1440/710 rpm	15092905
	10	M1 CSA Motor 3x200-240V/60Hz 1,5/0,75kW-1440/710 rpm	15092906
	10	3x400V/50Hz 1,5/0,75kW-1440/710 rpm	15092907
	10	M1 CSA Motor 3x460-480V/60Hz 1,5/0,75kW-1440/710 rpm	15092908
	10	3x380-415V/60Hz 1,5/0,75kW-1440/710 rpm	15092909
	30	Key B 8X5X30	2ZF20830
	40	V-belt SPZ-2 ø63	2JE00063
	50	Adapter 1108/ø28	2JE91128
	15090040		<b>Bearing unit sub assembly</b>
20		Bearing unit sub assembly	R5090040
50		Ball bearing 6208-2RS1 ø40/80	2BK00120
60		Angle contact bearing ø50/ø90x20	2BK30050
70		Nilos-ring 7210AVH	2BK97210

*AbraPol-10*  
*Spare Parts and Diagrams*

<b>Drawing</b>	<b>Spare Part list for AbraPol-10</b>		<b>Cat no.</b>
	<b>Pos.</b>		
15090045		<b>M2 Sample mover motor, complete</b>	
	10	M2 Gear-motor 3x200V/50Hz 0,37/0,25kW-300/150rpm	15092910
	10	M2 Gear-motor 3x220-230V/50Hz 0,37/0,25kW-300/150rpm	15092911
	10	M2 CSA Gear-motor 3x200-210V/60Hz 0,37/0,25kW-300/150rpm	15092912
	10	M2 Gear-motor 3x380-415V/50Hz 0,37/0,25kW-300/150rpm	15092913
	10	M2 CSA Gear-motor 3x460-480V/60Hz 0,37/0,25kW-300/150rpm	15092914
	10	M2 Gear-motor 3x220-240V/60Hz 0,37/0,25kW-300/150rpm	15092915
	10	M2 Gear-motor 3x380-415V/60Hz 0,37/0,25kW-300/150rpm	15099016
		<b>ABTDR Outlet kit</b>	05096901
	5	Outlet kit, complete	15090084
20	Cooling coil standard ø32 PVC	2NU29401	
15090084		<b>Outlet kit, complete</b>	
	20	Outlet, complete	15090094
	50	Outlet for cooling unit/drain	15090082
	80	Nylon tube, SUPERFLEX ø5/ø3.D4512	2NU12445
15090082		<b>Outlet kit for cooling unit/drain</b>	
	05	Outlet kit for cooling unit/drain	R5090082
	10	Solenoid valve 5/2-way 1/8"	2YM10010
	40	Gasket, PVC O-1/8	2IF00011
15090094		<b>Outlet kit, complete</b>	
	10	Actuator CRB1BW30-90S	2MK73090
	80-140	Outlet kit complete	R5090094
	100	Sealing ring SD 50x58x4	2IS30050

*AbraPol-10*  
*Spare Parts and Diagrams*








Drawing	Spare Part list for AbraPol-10	
	Pos.	Cat no.
	<b>ABTDO Additional Dosing Unit</b>	
		05096902
	10	Allen key 2.5 hardened DIN 911
	20	Tygon hose ENFT 21 ø2.06
	30	Silicone hose ø4/ø6
	40	Pressure water hose, white nr.880 2x3/4"
	60	Gasket ø11/ø24 x 1.5
	70	Gasket with filter 3/4"
	80	Bottle, 1L, for OP-S, complete
	90-100	Bottle, ½ l, for DP-suspension
	110	Lid for ½ l bottle, complete
		R4600603
15090085	<b>Additional Dosing Unit</b>	
	20	Solenoid valve 2/2 NC-NF-NC G1/8"
	30	Gasket, PVC O-1/8
	70	Gasket with filter 3/8"
	80	O-ring 11.50-1.50 72 NBR 872
	160	BOSCH motor w. print + axle extension
	230	Tygon hose ENFT 21 ø2.06
	270	Silicone hose ø2/ø4
	300	Silicone hose ø4/ø6
	330	Elbow pipe II
	340	Silicone hose ø4/ø6
	350	Silicone hose ø4/ø6
	430	Clear arm SL. M.½+D519" 90°
		2YM10140
		2IF00011
		2IX20405
		2IO15040
		14600550
		2NU91221
		2NU11452
		2NU11454
		14600037
		2NU11454
		2NU11454
		2NU94033
	<b>ABMEU Accurate Material Removal unit</b>	
	05	Accurate Material Removal unit
	10-30	Test push-buttons (3 sizes)
	60	Allen key 5, hardened DIN 911
	70	Angle pipe wrench 289 N 17
		05096903
		15099083
		15092919
		2GR00050
		2GR00317
15099083	<b>Accurate Material Removal unit</b>	
	10	1K 0.5W potentiometer
		2RP10020
	<b>ABTCO Disc Cooling</b>	
	10	Disc Cooling
	20	Pressure water hose, white, no. 880 2x3/4"
	30	Reducing branch, ½" internal 3/4" external
	40	Gasket ø11/ø24 x 1.5
	50	Gasket with filter 3/4"
	60	Allen key 4, hardened DIN 911
		05096904
		15090087
		2NU93020
		2NG30013
		13590359
		2IX20410
		2GR00040
15090087	<b>Disc Cooling</b>	
	10	Solenoid valve 24VAC yellow 337
	120	Grifl-hose k. arm PVC 3/8-ø10
		2YM12337
		2NU29312

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










- 110 Omnifit 200M
- 120 Omnifit 200M
- 20 Omnifit 50M
- 40 High Vacuum Grease
- 90 Omnifit Seal 50H

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














ABRAPOL 10, CIRCUIT DIAGRAM  
MAIN VOLTAGE, MAIN SUPPLY CIRCUIT.


FILE NAME: 1 5093100.DWG	Size A2	CAGE Code	DWG NO 15093100	Rev I
Friday, November 09, 2001	Scale	SLJ / SLJ	Sheet 1 of 1	c


- COLOR CODE:
- BK = black
  - BN = brown
  - RD = red
  - OG = orange
  - YE = yellow
  - GN = green
  - BU = blue
  - VT = violet
  - GY = grey
  - WH = white
  - PK = pink
  - BE = beige
  - RO = rose

Rev. A: Baan PDM Rev. B: Connections at X11 changed Rev. C: Wire no. Changed (DEM 09-11-01) Rev. D: Connection between L1 and L2 removed at 380-420V and 460-480V		STRUERS A/S VALHOEJES ALLÉ 176 DK-2810 ROEDOVRE DENMARK PHONE: + 45 3670 3500			
FILE NAME.: 5093130.DSN PAGE1.SCH		Size A3	CAGE Code  	DWG NO 15093130	Rev D
Thursday, April 24, 2003		Scale  	SLN / KFC		Sheet 1 of 1
<b>ABRAPOL 10, FUNCTIONAL DIAGRAM, SAFETY CIRCUIT</b>					


SUPPLY IN AND RS485 COMMUNICATION


INCREMENTAL DIAL KNOB


TOUCH PAD AND PUSH BUTTONS INPUT


OUTPUT TO POWER SIGN LEDs (AND DECOUPLING)


CPU AND MEMORY CIRCUIT


RESET CIRCUIT


ADDRESS SELECT & R/W LOGIC


JOYSTICK INPUT


RS232 INTERFACE


DISPLAY INTERFACE


CCFT INVERTER


TRANSFORMER INPUT

MOUNTING CONNECTIONS


GROUND CONNECTIONS


UAC-01; +5V, 1AT  
 +15V, 300mAT  
 -15V, 100mAT  
 +24V, 1,5AT

UAC-1; 24Vac, 4AT (CONTACTOR MAIN, DISC MOTOR, SAMPLE MOVER, RECIRK. AND WATER-VALVE)

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











































Se speciel tegning



M1  
Connection box on sample mover motor

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













Herewith declares that

<i>Product Name:</i>	Abrapol-10
<i>Type No:</i>	509
<i>Machine Type:</i>	Grinding and polishing machine

is in conformity with the provisions of the following directives:

**Safety of Machinery** 98/37/ EEC according to the following standard(s):  
EN ISO 12100-1:2003, EN ISO 12100-2:2003, EN ISO 13850:2006, EN574:1996; EN953:1997;  
EN954-1:1997; EN1050:1997; EN60204-1:1997.

**EMC-Directive** 2004/108/EEC according to the following standard(s):  
EN61000-6-1:2002, EN61000-6-3:2002.

**Low Voltage Directive** 2006/95/EEC according to the following standard(s):  
EN60204-1:1997.

**Supplementary Information** The equipment complies with the American standards:  
FCC part 15, subpart B og UL508 ; NFPA 70:1996; NFPA 79:1997

**The above has been declared according to the global method, module A**

Date: 04.12.2007

Christian Skjold Heyde,  
Vice President, R& D and Production, Struers A/S

**Hersteller**  
Struers A/S  
Pederstrupvej 84  
DK-2750 Ballerup, Danmark  
Telefon +45 44 600 800

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Christian Skjold Heyde,  
Vice- President, R& D et Production, Struers A/S

# AbraPol-10

**A high-tech  
automatic floor-  
based grinding and  
polishing machine  
with a high  
throughput**

*AbraPol-10, extremely sturdy and reliable*

**AbraPol-10 is a fast and tough grinding and polishing machine, ideal for laboratories or workplaces carrying out quality control testing of very large samples or a high volume of specimens**


**Unique Advantages:**

- High-speed preparation shortens your preparation time by up to 50%
- Consumables cost is significantly reduced
- User-friendly operations due to built-in method database
- Sturdy design for high volume preparation year after year
- Fits in any production facility
- Optional: Accurate Material Removal and Disc Cooling

Tough floor-based grinding and polishing machine

Low carbon steel

AbraPol-10 is a powerful, high capacity machine setting new standards for grinding and polishing. Never before has a preparation machine been so fast performing thus saving valuable time for the user.

A user-friendly interface makes it easy to operate the machine and establish new preparation methods. Metallog Guide methods are also included, covering the preparation of most materials. Thus, even operators with limited preparation experience can prepare all kinds of materials immediately.

AbraPol-10 handles quickly and easily just about any sample and fits in any production environment. Year after year it lives up to the demand of the daily high volume metallographic sample preparation.

### Reduced Preparation Time

AbraPol-10 is the fastest preparation machine on the market, setting new standards in metallographic sample preparation.

New movable dosing arm

### Preparation Method for low carbon steel:

	The original set-up	New set-up
Disc size (mm)	300	350
Sample holder size (mm)	160	200
Number of samples	6 x 30 mm dia.	8 x 40 mm dia.
	The original set-up	New set-up
Plane Grinding	MD-Piano 220	MD-Piano 220
Speed - disc/specimen holder (rpm)	300/150	300/300
Time (min)	Until plane	Until plane
Fine Grinding 1	MD-Allegro	MD-Allegro
Speed - disc/specimen holder (rpm)	150/150	300/300
Time	4 min.	2 min. 20 sec.
Polishing	MD-Dac	SP-PoliSat 1
Speed - disc/specimen holder (rpm)	150/150	300/300
Time	4 min.	2 min. 30 sec.
Oxide Polishing	MD-Chem	SP-PoliCel 1
Speed - disc/specimen holder (rpm)	150/150	150/150
Time (min)	2	1
<b>Total preparation time</b>	<b>10 min.</b>	<b>5 min. 50 sec.</b>

It is now possible to reduce the total preparation time by up to 50% compared to current technology available on the market and at the same time increase the number and area of samples by 30%.

Using AbraPol-10, the effect will be outstanding: double the throughput of large samples and save valuable time for the user.

The keys to these remarkable results are to be found through a combination of high rotation speed and the increased area of the preparation disc.

AbraPol-10 is designed for a new and larger 350 mm disc. The increased area of the disc allows an efficient utilisation of 200 mm specimen holders. Combining this with a 300 rpm rotation speed, the total preparation time is reduced by 50%. The total preparation time after plane grinding, e.g. for WC/Co sintered carbides will be reduced from 17 min

to 7 min, and for low carbon steel the reduction will be from 10 min to 5 min. 50 sec.

### Reduced Preparation Cost

AbraPol-10 is not only saving valuable time for the user, but also substantially reduces preparation cost. As a result of the reduced preparation time, the consumption of diamond suspension will decrease with 20%.

The large surface of the 350 mm disc significantly increases in the lifetime of the grinding and polishing surfaces. Tests have shown that the grinding and polishing surface lifetime increases by up to 70%. In short, the preparation time and the consumption of diamond suspension is reduced while the lifetime of the preparation surfaces is increased. That makes AbraPol-10 the most economic high volume preparation machine available on the market.

Large display

### Powerful Action

The power of a grinding and polishing machine has never been so strong, making it a sturdy and tough preparation machine. The main motor is built into a strong, welded frame. The whole machine is designed for vibration free preparation. Double start buttons is standard on the AbraPol-10 securing the user the maximum safety.

### Interface

The AbraPol-10 is an automatic machine controlled by a user-friendly interface. A large display makes it easy to select or create a preparation method. Only few operation keys are necessary to start or create a preparation method.

### Metallog Guide Methods Included

The Metallog guide methods which are covering the preparation of most materials are incorporated in the menu program. A suitable method can easily be selected when new materials have to be prepared. The method can even be copied, changed and saved for further use.

### Integrated bottle unit

The placement of the consumables is integrated in the front of the machine making the access for changing bottles or changing tubes in pumps extremely easy. The standard version comes with 3 pumps for suspension and lubricants. It can be expanded with 3 additional pumps of which one is for OP Suspensions.

The bottle unit is protected by two strong clear plastic doors allowing an easy monitoring of the consumable levels. The tubes can easily go through the bottle unit to reach larger containers placed on the floor e.g. lubricants.

### Peristaltic pumps for dosing

Peristaltic pumps always supplies the correct amount of suspensions, lubricants and oxide polishing suspensions, just where it is needed. A unique feature of the peristaltic pumps is that they are entirely aerosol-free and thus represent an environment-friendly alternative.

### Optional

A material removal unit can be added for precise control of material removal rate. A new and advanced potentiometer is the basic principle in the material removal unit. It is a quick and easy operation that is completely software monitored. The removal can be adjusted with 10 µm intervals. This can be used with advantage in the following instances :

- To remove the specified amount of material when preparing brittle materials (e.g. plasma sprayed coatings)
- In three dimensional analysis where material has to be removed in successive steps.
- Or simply to make sure to remove sufficient material from a damaged surface (e.g. a bad sectioning)

The user saves time and consumables, and it allows an entire preparation cycle to be performed on an adequate ground specimen. Another accessory for AbraPol-10 is the disc cooling. A combination of high-speed preparation and high throughput will generate heat and increase the temperature of the preparation disc. The increase in temperature will boost the evaporation

of lubricant and the cloth might run dry which can result in damage of both cloth and specimens. For users with high throughput it is therefore highly recommended to use the disc cooling system.

### Abra System

AbraPol-10 is a highly adaptable machine that will meet all your needs. In its basic form it comes with a standard tap water cooling system and a dosing unit with three pumps. But it can also be supplied with a recirculation cooling unit, for example if the user wants to reduce the risk of corrosion in connection with the preparation of ferrous materials.

The Abra System also includes the plane grinding machine, Abraplan. A company with high capacity requirements will find the ideal solution is an Abraplan for plane grinding, an AbraPol-10 for fine grinding and a second AbraPol-10 for polishing.

AbraPol-10 includes an integrated dosing unit for Diamond Suspension and Lubricants

**Struers A/S**  
 Pederstrupvej 84  
 DK-2750 Ballerup, Denmark  
 Phone +45 44 600 800  
 Fax +45 44 600 801  
 struers@struers.dk

Technical Data		
Subject		Specifikations Metric/International
<b>Disc</b>	Rotational speed Size	150/300 rpm 300/350 mm (12"/14")
<b>Specimen holder</b>	Speed Direction Force Max size of sample holder	150/300 rpm CW/CCW 50-700 N 200 mm (8")
<b>Dosing Unit</b>	Level No. of pumps	0-20 3 standard + 3 optional
<b>Software and Electronics</b>	LC Display Controls	240x128 characters Touch pads/Push-turn knob
<b>Environment</b>	Noise level (idle 300/300) Surrounding temperature HumidityMax.	70 dB(A) 5-40°C 95%RH
<b>Supply Voltage</b>	Power consumption No. of phases Output, main motor	1,9 kW 3 (3L+PE) 1,5 kW
<b>Water mains</b>	Water supply	Max. 10 bar
<b>Compressed air</b>	Compressed air supply	6-10 bar
<b>AbraPol-10 dimensions and weight</b>	Width Depth Height Weight	670 mm (26.5") 930 mm (36.6") 1500 mm (59") 250 kg

Struers' equipment is in conformity with the provisions of the applicable International Directives and their appurtenant Standards. (Please contact your local supplier for details)

## Specifications

AbraPol-10	Code
Automatic microprocessor controlled machine for grinding and polishing. Recirculation cooling unit (TREC) or water supply kit (ABWAT) are ordered separately	ABTEN

## Accessories

Water Supply kit	Code
For control of tap water supply. For AbraPol-10 without Recirculation cooling unit (TREC)	ABWAT

Recirculation Cooling Unit	Code
For water cooling on various equipment. Capacity 30 l. 3-phase. For AbraPol-10 without water supply kit (ABWAT)	TRECI

Rolling Pallet for TREC	Code
	TREJA

Disc for Magnetic Fixation	Code
MD-Disc, 350 mm dia. Gross weight 4.6 kg	DEMFI
MD-Disc, 300 mm dia. Gross weight 3.9 kg	DEMAL

Water Level Sensor	Code
Indicating the level of water in the recirculation cooling unit (TREC)	PAMWA

Disc Cooling	Code
Disc cooling for AbraPol-10. Only to be used with tap water. In connection with recirculation cooling unit (TREC) the outlet kit (ABTDR) is required.	ABTCO

Additional Dosing Unit	Code
Three additional pumps for dosing of OP-Suspensions, Lubricants and Suspensions. In connection with recirculation cooling unit (TREC) the outlet kit (ABTDR) is required.	ABTDO

Outlet Kit	Code
Shift valve for recirculation cooling unit (TREC) or drain	ABTDR

Accurate Material Removal Unit	Code
For accurate removal of material. Stops process after required amount of material has been removed.	ABMEU

Struers' products are subject to constant product development. Therefore, we reserve the right to introduce changes in our products without notice.

**USA and CANADA**  
**Struers Inc.**  
 24766 Detroit Road  
 Westlake, OH 44145-1598  
 Phone +1 440 871 0071  
 Fax +1 440 871 8188  
 info@struers.com

**SWEDEN**  
**Struers A/S**  
 Smältvägen 1  
 P.O. Box 11085  
 SE-161 11 Bromma  
 Telefon +46 (0)8 447 53 90  
 Telefax +46 (0)8 447 53 99  
 info@struers.dk

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**Struers GmbH**