# Accutom-5

## Instruction Manual

Manual No.: 14947001

Date of Release 15.11.2012

*Accutom-5*  
*Instruction Manual*

<b>Table of Contents</b>	<b>Page</b>
User's Guide .....	1
Reference Guide.....	25
Quick Reference Guide .....	47

---

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the *Serial No.* and *Voltage* on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

**Instruction Manuals:** Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

**Service Manuals:** Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

**Original instructions.** The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2012.

**Struers A/S**  
Pederstrupvej 84  
DK-2750 Ballerup  
Telephone +45 44 600 800  
Fax +45 44 600 801

---

## **Accutom-5 Safety Precaution Sheet**

### **To be read carefully before use**

1. The operator should be fully instructed in the use of the machine and its cut-off wheels according to the Instruction Manual and the instructions for the cut-off wheels.
2. The machine must be placed on a safe and stable support table.
3. Be sure that the actual voltage corresponds to the voltage stated on the back of the machine. The machine must be earthed.
4. Use only intact cut-off wheels. The cut-off wheels must be approved for min. 3000 rpm. If other cut-off wheels or saw blades are used, make sure that the speed setting of Accutom-5 complies with the max. speed for the cut-off wheels or saw blades.
5. Observe the current safety regulations for handling, mixing, filling, emptying and disposal of the additive for cooling fluid.
6. The sample must be securely fixed in the specimen holder.
7. Do not touch the sample, the specimen holder head or the cut-off wheel while positioning the sample with the POSITION controls.
8. Never try to open the cover before the cut-off wheel has stopped completely.

---

The equipment should only be used for its intended purpose and as detailed in the Instruction Manual.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

---

## Disposal

Equipment marked with a WEEE symbol  contain electrical and electronic components and must not be disposed of as general waste.

Please contact your local authorities for information on the correct method of disposal in accordance with national legislation.

# User's Guide

Table of Contents	Page
<b>1. Getting Started 3</b>	
Checking the Contents of Packing .....	3
Placing Accutom-5 .....	3
Getting Acquainted with Accutom-5 .....	3
Supplying Power .....	4
Changing the voltage setting .....	4
Recirculation Unit .....	6
Software Settings .....	7
Configuration Menu.....	7
Setting the Language.....	8
<b>2. Basic Operations 9</b>	
Using the Controls.....	9
Front Panel Controls of Accutom-5.....	9
Groups of Keys .....	9
Acoustic Signals.....	9
Location of Main Switch .....	9
Front Panel Controls .....	10
Display .....	11
Reading the Display .....	12
Changing/Editing Values .....	13
Numeric Values.....	13
Alphanumeric Values .....	14
Positioning the Sample .....	15
Reference Position.....	15
Absolute Position .....	15
Relative Position .....	15
Relative Zero.....	15
Stop Position.....	15
Cutting.....	16
Changing the Cut-off Wheel.....	16
Clamping the Sample and Specimen Holder .....	17
Positioning the Sample .....	17
Setting the Cutting Parameters .....	18
Wheel.....	18
Speed.....	18
Feed.....	18
Force.....	18
Rotation.....	19
Water .....	19

Starting the Cutting .....	19
During Cutting .....	20
Changing the Feed Speed .....	20
Retracting the Sample .....	20
Force Limit Reached.....	20
Stopping the Cutting.....	21
Automatic Stop.....	21
Manual Stop.....	21

### **3. Maintenance**

Daily Service .....	22
Checking the Recirculation Unit.....	22
Weekly Service .....	23
Refilling the Cooling Water Tank.....	23
Emptying and Cleaning the Tank.....	23
Refilling the Tank .....	23
Maintenance of Cut-off Wheels.....	24
Maintenance of Diamond and CBN Cut-off Wheels.....	24
Storing of Abrasive Cut-off Wheels.....	24

## 1. Getting Started

### Checking the Contents of Packing

In the packing box you should find the following parts:


- 1 Accutom-5
- 2 Mains cables
- 1 Specimen Holder with parallel vice
- 1 Flange for cut-off wheels (42 mm dia.)
- 1 Stop pin
- 1 Spanner, 17mm
- 1 Small grate
- 1 Large grate
- 1 Allen key, 2 mm
- 1 Allen key, 2.5 mm
- 1 Allen key, 3 mm
- 1 Allen key, 4 mm
- 1 Allen key 5 mm
- 2 Screws M4x20
- 2 Screws M4x35
- 1 Set of Instruction Manuals

### Placing Accutom-5

Accutom-5 should be placed on a stable and plane (tolerance  $\pm 1$  mm) table. The table must be able to carry a weight of min. 50 kg.

### Getting Acquainted with Accutom-5

Take a moment to familiarise yourself with the location and names of the Accutom-5 components.

- A Front panel/Front panel control(s)
- B Main switch
- C Recirculation Cooling Unit
- D Cut-off wheel
- E Specimen holder head

## Supplying Power

Always remember to switch the power off when installing electrical equipment!

**IMPORTANT**  
Check that the mains voltage corresponds to the voltage stated on the type plate on the back of the machine.

## Changing the voltage setting

The factory setting for Accutom-5 is 240V.

If the factory setting is not the correct setting for your mains supply the setting can be changed from 240V to 220V:

- Pull out the fuse holder from the cable terminal at the back of the machine.
- Pull out the voltage switch and turn it to the correct setting.

Voltage Required	Setting
230 or 240V	240V
200 to 220V	200V


Note: The two additional settings, 110V and 120V are not to be used.

- Reinsert the voltage switch.
- Put the fuse holder back into the cable terminal.

Accutom-5 is shipped with 2 types of Mains cables:

Single-phase Supply

The 2-pin (European Schuko) plug is for use on single-phase connections.

If the plug supplied on this cable is not approved in your country, then the plug must be replaced with an approved plug. The leads must be connected as follows:

- Yellow/green: earth
- Brown: line (live)
- Blue: neutral

3-phase Supply

The 3-pin (North American NEMA) plug is for use on 3-phase power connections.

If the plug supplied on this cable is not approved in your country, then the plug must be replaced with an approved plug. The leads must be connected as follows:

- Green: earth
- Black: line (live)
- White: line (live)

Connection to the Machine

Both cables are on the other end equipped with an IEC 320 cable connector that has to be connected to the Accutom.

**WARNING!**

The output voltage from this cable is 200 – 240V and not 110V. DO NOT use this cable to connect equipment that use a 110V power supply. Failure to adhere to this may result in material damage.

## Recirculation Unit

- Pull out the recirculation tank.
- Fill the tank with 3.88 l water and 120 ml Struers Additive for cooling fluid. The water level should be 5 mm below the edge of the front hole in the tank cover.

**IMPORTANT**

Always ensure that there is sufficient water in the tank as the recirculation pump will be damaged if it is run dry.

- Check that the cover is fitted securely into place in the recirculation tank and push the drawer with the recirculation tank into place again.
- Check that the end of the inlet tube has dropped into position.

**Note:**

Changing of cooling water should be done at least once a month.

**IMPORTANT**

Always maintain the correct concentration of Struers Additive in the cooling water (percentage stated on the container of the Additive). Remember to add Struers Additive for cooling fluid each time you refill with water. Do not use any oil, petrol, or turpentine-based additives, only Struers additive.





## Software Settings Configuration Menu

Switch on the power at the main switch located at the back of the machine. The following display will appear briefly:


Afterwards the display will change to the same screen which was shown before Accutom-5 was switched off, usually a cutting method. When switching Accutom-5 on for the first time, the display to appear should be the MAIN MENU. If the heading in the display is different, press Esc, until the MAIN MENU appears. (A long beep can be heard).

The MAIN MENU is the highest level in the menu structure. From here you can go to configuration, manual functions and cutting methods.


Setting the Language

 Press MENU UP/DOWN ▲▼ to select OPTIONS.

 Press ENTER  to activate the OPTIONS menu.

 Press ENTER  to activate the LANGUAGE Menu.


 Press MENU UP/DOWN ▲▼ to select the Language you prefer.

 Press ENTER  to accept the language. The OPTIONS Menu now appears in the language you have chosen.

## 2. Basic Operations

### Using the Controls Front Panel Controls of Accutom-5

#### Groups of Keys

- A** Programming and monitoring functions
- B** Positioning of the specimen holder
- C** Start/stop of Accutom-5

#### Acoustic Signals

-  **Short Beep:** when a key is pressed, a short beep indicates that the command has been accepted.
-  **Long Beep:** a long beep indicates that the key is inactive at that moment.










#### Location of Main Switch

The main switch is located on the back of the machine.


**Please Note...**

The contents of the program memory are not lost when the main switch is turned off.

## Front Panel Controls

Name	Key	Function	Name	Key	Function
① MAIN SWITCH		The main switch is located on the back of the machine.	⑦ ESC		Leaves the present menu or specimen holder position screen.
② START		Starts the cutting process according to the pre-set method. The recirculation water, if selected, is turned on.	⑧ FAST POSITION LEFT		Changes to POSITION menu or moves the specimen holder to the left in the X-direction in steps of 100 µm. Keep the key pressed to increase the speed.
③ STOP		Stops the cutting process. The recirculation water, if selected, is turned off.	⑨ FAST POSITION RIGHT		Changes to POSITION menu or moves the specimen holder to the right in the X-direction in steps of 100 µm. Keep the key pressed to increase the speed.
④ FUNCTION KEYS		Controls for various purposes. See the bottom of the individual screens.	⑩ POSITION LEFT		Changes to POSITION menu or moves the specimen holder slowly to the left in the X-direction in steps of 5 µm. Keep the key pressed to increase the speed.
⑤ MENU		Scrolls up (▲) or down (▼) in the menu tree structure of Accutom-5. When setting a parameter the value is increased (▲) or decreased (▼).	⑪ POSITION RIGHT		Changes to POSITION menu or moves the specimen holder slowly to the right in the X-direction in steps of 5 µm. Keep the key pressed to increase the speed.
⑥ ENTER		Selects a marked parameter value or chooses a menu.	⑫ POSITION UP/DOWN		Changes to POSITION menu or moves the specimen holder up- or downwards in the Y-direction in steps of 100 µm. Keep the key pressed to increase the speed.

Display



- A Heading
- B Inverted text. Cursor
- C Function key options
- D Arrow indicates, that there are more lines in the picture




**Please Note**  
The examples of display screens in this Instruction Manual show a number of possible texts. The actual display screen may differ from the examples in the Instruction Manual.

## Reading the Display



The display can show various kinds of information, for example about the cutting method or about the sample position. A screen for a cutting method could look as the following example:

Pressing one of the X-axis positioning keys will change the screen to the following:

Pressing one of the Y-axis positioning keys will change the screen to the following:

Press Esc, to return to the cutting method screen.

## Changing/Editing Values

Depending on the type of value, there are two different ways of editing.

### Numeric Values

**Cutting Method: 1. EMPTY METHOD**

Wheel : **B0D15** Speed : 2700 rpm  
Feed : 0.100 mm/s Force : MEDIUM  
Rotation: OFF  
Water : ON

↓


 Press MENU UP/DOWN  to select the numeric value you want to change, e.g. Feed:

↓

**Cutting Method: 1. EMPTY METHOD**

Wheel : B0D15 Speed : 2700 rpm  
Feed : **0.100** mm/s Force : MEDIUM  
Rotation: OFF  
Water : ON

↓

 Press ENTER ↵, to edit the value.

↓ Two square brackets [ ] appear around the value.

**Cutting Method: 1. EMPTY METHOD**

Wheel : B0D15 Speed : 2700 rpm  
Feed [0.100]mm/s Force : MEDIUM  
Rotation: OFF  
Water : ON

↓

 Press MENU UP/DOWN  to increase or decrease the numeric value.




↓

 Press ENTER ↵, to accept the new value.



↓ or

 Press Esc, to keep the original value.


Alphanumeric Values

A submenu appears.

or

## **Positioning the Sample**

After clamping the sample in the specimen holder, the holder is placed in the specimen holder head and fixed there. To move the specimen holder and the sample, the POSITION keys ▼▲▶▶◀◀ are used. The display shows the position in either X- or Y- direction.

### *Reference Position*

Every time the power is switched on while the cover is closed, Accutom-5 checks its own reference position. The sample holder head will be moved back as far as possible, to the reference position (X=0.000 Y=0.0), and thereafter it will return to the position where it was before the power was switched on.

### *Absolute Position*

The absolute position shows the total distance the specimen holder has travelled from the reference position.

### *Relative Position*

The relative position equals that of the absolute position until it is set to zero at a desired point. By setting it to zero, calculation of the sample movement close to the cut-off wheel is made easier. The screen value always relates to the distance the specimen holder has travelled since being set to zero.

### *Relative Zero*

The relative zero position is the point where the relative position in either the X- or Y- position was set to zero. Having completed the cutting process, the sample holder automatically returns to this point. Pressing F3 when in the X-or Y- position screen also returns the sample to the relative zero position.

### *Stop Position*


A stop position can be set to halt the cutting process at a precise point. After reaching this point the sample will be retracted and returned to the relative zero position.

### **Compensating for Wheel Wear**

Please make sure when setting the stop position to compensate for possible wear of the cut-off wheel. This is especially important when using Al<sub>2</sub>O<sub>3</sub> or SiC wheels.

## Cutting


### Changing the Cut-off Wheel

- Open the cover and swing the right cooling tube in the air.
- Insert the stop pin in the hole of the inner flange.
- Use the spanner (17 mm) to loosen the flange screw.
- Remove the outer flange and the wheel.

#### **IMPORTANT**

The tolerance between the spindle and inner flange is very small which means that the two surfaces must be absolutely clean. Never try to squeeze the cut-off wheel on as this may damage the spindle. If there are any small burrs remove them with grinding paper (grit size 1200).

- Mount the new cut-off wheel and remount the outer flange, with the machined face towards the inner flange.
- Insert the locking pin in the hole in the inner flange.
- Fasten the flange screw gently using the spanner.
- Swing the right cooling tube back in its place.

### Clamping the Sample and Specimen Holder


- Clamp the sample securely in the specimen holder using the appropriate Allen key.
- Fasten the specimen holder in the specimen holder head by pushing the specimen holder into the dovetail fixture and tightening the screw with an Allen key.
- When cutting with rotation or oscillation, the sample and the specimen holder should be clamped so that they rotate evenly around the centre of the sample. This way the fastest cutting is obtained as the cut-off wheel will be cutting most of the time and the possibility of damaging the cut-off wheel is limited.

### Positioning the Sample

- Move the sample into the correct start position, close to the cut-off wheel by using the POSITION keys.
- Reset the relative position in both X- and Y-direction by pressing F1 when in the respective screens.
- Set the stop position to define the length of the cut:

There are two ways of setting the stop position:


1. If you know the size of your sample:
  - Press ENTER ↵ and use MENU UP/DOWN ⏮▲⏭▼ to set the stop value.
  - Press ENTER ↵ again to accept the value or Esc to cancel.
2. If you do not know the exact size of your sample:
  - Position the sample using the POSITION key ► so you can move it along the cut-off wheel.
  - Move the sample along the cut-off wheel to the required stop position using the POSITION key ▲.
  - Press F2:SET STOP to record the stop position.
  - Press F3:MOVE TO REL.ZERO to move the sample back into the initial position.
  - Position the sample correctly in front of the cut-off wheel using the POSITION key ◀. If possible use F3:MOVE TO REL.ZERO.

**Setting the Cutting Parameters** All cutting parameters can be edited independently of each other.

See section Changing/Editing values for details on how to change values.

**Method Screen**

*Wheel* All Struers cut-off wheels are listed with their default values for force limit, wheel thickness and speed. Also eight undefined wheels are included for your own choice of cut-off wheels. (See also configuration of cut-off wheels).

*Speed* The speed of the cut-off wheel can be set between 300 and 3000 rpm in steps of 100 rpm. The default values for all Struers cut-off wheels are already saved together with the wheel definitions. Whenever a different cut-off wheel is selected the recommended speed will automatically be inserted.

*Feed* The feed speed can be set between 0.005 mm/s and 3.000 mm/s. (See the section Cutting Parameters in chapter 3. Consumables for recommended settings).

*Force* The force limit is a protection for the cut-off wheel and/or the samples to be cut. Depending on the thickness or strength of the cut-off wheels the force limit can be set to three different levels, LOW, MEDIUM and HIGH. The default values for all Struers cut-off wheels are already saved together with the wheel definitions. Whenever a different cut-off wheel is selected, the recommended force limit will automatically be inserted.

Sample rotation can be an advantage when cutting large, very hard, coated or very long samples. On Accutom-5, three different settings are possible.


### Rotation

Only for Accutom-5 with rotation

Rotation: Off


The sample does not rotate.

Rotation: Continuous

The sample is rotating around its centre.  
The speed can be set to three different levels, 1, 2 or 3.  
The direction of rotation can be set to either counter- or co-rotation. Counter-rotation is recommended.  
Press F3 to test sample rotation with the set parameters.

Rotation: Oscillating

The sample is rocking/oscillating around its centre.  
The speed can be set to three different levels, 1, 2 or 3.  
The angle can be adjusted from 10° to 400°.  
Press F3 to test sample rotation with the set parameters.

#### Note

When the sample is moved in the x- or y- direction, with Rotation set to Cont. or Osc. a warning message, "Rotation Mode!", will briefly flash on the screen as a reminder that rotation mode is selected.

Water

The cooling water can be set to either on or off. For all normal cutting operations the setting should be on.


### Starting the Cutting

Before you start cutting make sure that the cutting chamber is clean to ensure a free flow of cooling water.

- Position the sample correctly.
- Set the correct cutting parameters.
- Make sure that the correct cut-off wheel is mounted.
- Close the cover of the machine.
- Press START .

### During Cutting

After Start has been pressed, the display changes to the following:

Both the pre-set and the actual feed speed are displayed. The horizontal bar is used to display the feed speed graphically.


### Changing the Feed Speed

During the cutting process the feed speed can be changed. Simply press Enter ↵ and change the feed speed to the correct value. Press Enter ↵ again to confirm the change or Esc to cancel.

### Retracting the Sample

The only possibility of moving the sample after the cutting process has been started is to retract the sample from the cut-off wheel. This can be done by pressing the POSITION DOWN ▼ key. The forward movement of the sample is stopped and the sample moved backwards until the key is released again. Then the sample will start again to move forward with the pre-set feed speed. If the pre-set feed speed can not be achieved because the force limit is reached, Accutom-5 automatically reduces the feed speed to the maximum possible speed. This value is displayed, and after the cut is finished, the following message is shown on the display:

### Force Limit Reached

For similar samples to be cut afterwards, the feed speed should be reduced to the new value or below.

**Stopping the Cutting**  
*Automatic Stop*

- Accutom-5 automatically stops the cutting process at the pre-set stop position. (Make sure to compensate for possible wheel wear when setting the stop position).
- The sample is then retracted and the cut-off wheel is stopped.

*Manual Stop*

- The cutting process can be stopped at any time during operation by pressing the STOP ⏹ key. The sample remains in its actual position and the process stops there.
- To avoid stopping the process while the cut-off wheel still is in the sample, press the POSITION key ▼ for the Y-direction and then press F2:SET STOP. The sample is retracted immediately and the process will be stopped when the relative zero position is reached.
- The display changes back to what it was before cutting was started.

### 3. Maintenance

#### Daily Service

- Clean the cutting chamber with a damp cloth. Do not use tap water as you risk overflow in the recirculation cooling water tank. Remember to remove all dirt particles from the grate.
- Clean the specimen holder head and the clamps for the dovetail feed.
- Clean the flanges.
- Clean the transparent cover with a damp cloth.

**WARNING!**

Do not use alcohol, acetone or similar solvents.

#### Checking the Recirculation Unit

The cooling unit should be checked for cooling water after 8 hours use or at least every week. Refill the cooling unit if the flushing pump cannot reach the cooling water.

Remember to add Struers Additive for Cooling Fluid: One part of Additive for 33 parts of water.

Use a refractometer to check the concentration of additive.

Concentration = Brix value. The concentration of additive should lie between 2.7 and 3.3 %. Add Struers Additive for Cooling Fluid if the concentration is too low.

**IMPORTANT**


Always maintain the correct concentration of Struers Additive in the cooling water (percentage stated on the container of the Additive). Remember to add Struers Additive each time you refill with water.

## Weekly Service

- Clean the specimen holder(s): movable parts, dovetail feeds and screws. Lubricate with acid free oil.
- Clean the cutting chamber and cover thoroughly. Remember to remove all dirt particles from the grate.

## Refilling the Cooling Water Tank

### Emptying and Cleaning the Tank

Replace the cooling water in the Recirculation Cooling Unit at least once a month.

- Carefully pull out the drawer and lower the hinged front plate.
- Carefully pull out the flat tank.
- Remove the cover plate and empty the tank into a drain approved for waste chemicals.
- Clean the tank, the dividers and the cover with tap water.
- Place the dividers correctly in the tank and replace the cover plate.
- Press the cover plate into the tank and make sure it is seated firmly.

#### **IMPORTANT**

The container and the cover plate have to be placed correctly, otherwise the pump cannot suck up the water.

#### **IMPORTANT**

Flush the recirculation system with clean water if Accutom is not to be used over longer periods of time. This will prevent any dried residue of cutting material from damaging the inside of the pump.

### Refilling the Tank

- Fill the tank with 3.88 l water and 120 ml Struers Additive. The water level should be 5 mm below the edge of the front hole in the tank cover.
- Push the drawer with the recirculation tank back into place.

#### **IMPORTANT**

Always ensure that there is sufficient water in the tank as the recirculation pump will be damaged if it is run dry.

#### **IMPORTANT**

Always maintain the correct concentration of Struers Additive in the cooling water (percentage stated on the container of the Additive). Remember to add Struers Additive each time you refill with water.

**Maintenance of Cut-off Wheels**

*Maintenance of Diamond and  
CBN Cut-off Wheels*

The precision of diamond and CBN cut-off wheels and thus the cut depends on how carefully the following instructions are observed:

- Never expose the cut-off wheel to overload, such as heavy mechanical load, or heat.
- Store the cut-off wheel in a dry place, horizontally on a plane support, preferably under light pressure.
- A clean and dry cut-off wheel does not corrode. Therefore, clean and dry the cut-off wheel before storing. If possible, use ordinary detergents for the cleaning.
- Regular dressing of the cut-off wheel is also part of the general maintenance (see ADVANCED OPERATIONS).

*Storing of*

*Abrasive Cut-off Wheels*

These cut-off wheels are sensitive to humidity. Therefore, do not mix new, dry cut-off wheels and used, humid ones. Store the cut-off wheels in a dry place, horizontally on a plane support.





# Reference Guide

Table of Contents	Page
<b>1. Advanced Operations</b>	
Configuration Menu .....	26
Cut-off Wheels .....	27
Options.....	28
Manual Functions.....	29
Sample Rotation .....	30
Recirculation Pump.....	31
Dressing.....	32
Method Options.....	34
Saving a Method .....	34
Copying a Method.....	35
Inserting a Method .....	36
Resetting a Method.....	37
Editing Names.....	38
Name Editing Principles.....	39
<b>2. Accessories</b> .....	40
<b>3. Consumables</b>	
Cut-off Wheels .....	41
Optimising the Cutting Results.....	41
Cutting Parameters .....	41
Common Cutting Questions.....	42
Consumables .....	42
<b>4. Trouble-Shooting</b> .....	43
<b>5. Technical Data</b> .....	45
<b>6. Menu Overview</b> .....	46

## 1. Advanced Operations

### Configuration Menu

Press Esc until you reach the Main Menu. (A long beep can be heard).

## Cut-off Wheels

If you want to use other than Struers cut-off wheels, these must be configured to contain the correct values for force limit, wheel thickness and speed.

To configure your own cut-off wheels,


 Press MENU UP/DOWN ▲▼ to select CUT-OFF WHEELS.

↓  
 Press ENTER ↵ to activate the CUT-OFF WHEELS MENU.


↓  

CONFIGURATION OF CUT-OFF WHEELS	
Cut-off Wheel :	USER2
Force Limit :	<b>MEDIUM</b>
Wheel Thickness:	0.50 mm
Speed :	3000 rpm
F1:PREV. WHEEL F2:NEXT WHEEL	
F4:RENAME	

↓  
 Press MENU UP/DOWN ▲▼ to select the value to be changed.

↓  
 Press ENTER ↵ to edit the value

After changing the value,

↓  
 Press ENTER ↵ to accept the new value and to return to the previous screen.

↓  
**F4** Press F4 to rename the cut-off wheel.  
(see Editing Names for details).

In the options menu various items can be configured.

*Options*

Language

See how to change the language in the Getting Started section of this Instruction Manual.

X-retract 

After the stop position is reached, the sample is normally moved straight back to the relative zero position. To avoid any contact with the cut-off wheel during that movement, the sample can be retracted from the wheel before it is repositioned by setting X-retract to ON. After reaching the relative zero position, the sample is then moved back into the original position.

***IMPORTANT***

Make sure that the Y-stop position is set correctly. If the sample is not cut through before the sample is retracted the cut-off wheel might be damaged.

Contrast

F1: DEFAULT VALUE

The contrast setting of the display can be adjusted.



Pressing F1 will reset the inverted option to the factory adjusted default value.

## Manual Functions

Several functions on Accutom-5 can also be carried out manually. All of these functions are available from the Manual Functions Menu.

Press Esc until you reach the Main Menu.

Sample Rotation

Press MENU UP/DOWN ▲▼ to select SAMPLE ROTATION.

Press ENTER ↵ to activate the SAMPLE ROTATION menu.

SAMPLE ROTATION	
Sample Rotation:	<b>CONTINUOUS</b>
Speed Level :	1
Direction :	> <
F1:START ROTATION	

Press MENU UP/DOWN ▲▼ to select the value to be changed.

Press ENTER ↵ to edit the value.

After changing the value,


Press ENTER ↵ to accept the new value and to return to the previous screen.

Press F1 to start rotation and check that everything has been aligned correctly. Press F1 again to stop rotation.


Recirculation Pump

 From the MANUAL FUNCTIONS MENU press MENU UP/DOWN  to select RECIRCULATION PUMP.

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
*Dressing*

Diamond and CBN cut-off wheels may have to be dressed to open up the wheel again and achieve maximum effect after cutting of ductile materials. This can be done using the following routine:

Dressing of Cut-off Wheels

- Mount the cut-off wheel to be dressed.
- Clamp the dressing stick in the sample holder and position it correctly.
  - The dressing stick has to be in front of the wheel so that the wheel is cutting slices of the dressing stick.

 ▲ From the MANUAL FUNCTIONS MENU Press MENU  
 ▼ UP/DOWN  ▲ ▼ to select DRESSING OF CUT-OFF WHEEL.



 Press ENTER ↵ to activate the dressing option.

DRESSING OF CUT-OFF WHEEL	
Feed Speed :	<b>0.5</b> mm/s
Wheel Speed :	3000 rpm
Cooling Water:	ON
Y-Distance $\blacktriangle\blacktriangledown$ :	10.0 mm
F1:START DRESSING	

↓  
 Press MENU UP/DOWN  to select the value to be changed.

↓  
 Press ENTER  to edit the value.

↓  
 Press ENTER  to accept the new value and to return to the previous screen.

↓  
 Press F1 to start the dressing operation.

## Method Options

### Saving a Method

While working with a specific method, you can save the changes you have made in the database.

- As soon as a parameter is changed in a cutting method, F4:SAVE will be shown on the bottom line of the display.

Cutting Method: 1. EMPTY METHOD	
Wheel :	<b>B0D15</b> Speed : 2700 rpm
Feed :	0.200 mm/s Force : MEDIUM
Rotation:	OFF
Water :	ON
F4:SAVE	

**F4** Press F4:SAVE to save the changes after all necessary parameters have been altered.

SAVE METHOD	
Save changes in method 1. EMPTY METHOD ?	
ENTER: YES	ESC: NO

 Press ENTER ↵ to save the modified method.

Unless you have already named the method, the name will automatically change from EMPTY METHOD to UNNAMED METHOD. This will always show that at least one parameter has been changed compared to a default method.







### **Important**

When saving changes, the original method will be overwritten. If you want to keep the original method, you should make a copy of the method with a new name, thus making the changes in a copy rather than changing the original method. See Copying a Method.

### Copying a Method


Copying a method is a shortcut to creating a new cutting method on the basis of an existing one.

- Select the cutting method you want to copy from.

### Inserting a Method

If you want to insert the method in a different method:

 Press MENU UP/DOWN ▲▼ to select the method, in which you want to insert the method.

If it is not an EMPTY METHOD:

 Press ENTER  to see the method.

**Esc** Press Esc to return to the above display.

**F2** Press F2:INSERT to insert the method from the buffer.

**INSERT METHOD**

Do you want to copy from buffer to method 4. EMPTY METHOD ?

ENTER: YES                      ESC: NO

 Press ENTER  to accept insertion of the method.







### Resetting a Method

If a method is not in use any more, it should be reset.

The parameters will change to default values, which can easily be changed to a new method later on.

The method name will change to: EMPTY METHOD, telling you, that you can copy to or modify in this method without replacing any valuable data.






- Select the method you want to reset.

Editing Names

The names for the methods and cut-off wheels can be edited and changed to suit your preference.

- Select the method or cut-off wheel you want to rename.

### Name Editing Principles

- Place the main cursor on the character you want to change, using F1:LEFT or F4:RIGHT. Use MENU DOWN  to move to the character set in the next line. An auxiliary cursor in the text line shows the position in the method name.
  - Write the new name using the following keys:
    - F1** Moves the main cursor to the left
    - F2** Inserts a space in the text
    - F3** Deletes one character to the left in the text
    - F4** Moves the main cursor to the right
- ↵ ENTER places the new character in the method name and moves the auxiliary cursor in the name to the right. Repeat the procedure for each character.

↓ Write the new name using the above name editing procedures.

 Press Esc to leave the editor again.

↓

**NAME CHANGED**

The name is changed, do you want to keep the old name, or accept the new name ?

ENTER:ACCEPT NEW NAME ESC:KEEP OLD NAME

↓

 Press ENTER ↵ to accept the new name.

## 2. Accessories

### Flange Sets

Specification	Cat. No:
For cut-off wheels 65 mm dia. For use on Accutom-5 to achieve highest precision.	04946902

### Specimen Holders

Specification	Cat. No:
For general use. Vice type with max. opening 60 mm	04946903
For round or square specimens. Teardrop type with max. opening $\varnothing$ 40 mm / $\varnothing$ 1 1/2"	04946904
For irregular specimens, with 7 screws. Max. width 40 mm / 1 1/2"	04946905
With goniometer	04276911
For adhering specimens	04276912
With ceramic vacuum chuck for thin sections	04276913
For small specimens. Vice type	04276915
Double parallel vice	04946909
<i>Joints to be mounted between the dovetail and the specimen holder</i>	
Tilting joint with max. angle $\pm 10^\circ$	04946906
Angling joint with max. angle $+30/-90^\circ$	04946908
<i>Base Plate</i>	
With dovetail. For mounting of other types of specimen holders	04276914

### 3. Consumables

#### Cut-off Wheels

Please refer to the Selection Guide in the [Struers Cut-off Wheels brochure](#).

Accutom-5 accepts cut-off wheels with a diameter of 75 mm (smallest flange is  $\varnothing 42$  mm) up to 152 mm. Hole 12.7 mm.

**IMPORTANT!**  
Always use large flanges (65 mm) for the highest possible accuracy.

#### Optimising the Cutting Results Cutting Parameters

Use the following table to select proper wheel and cutting parameters according to the sample material.

Recommended Cutting Parameters				
Material	Hardness	Force limit	Feed speed [mm/s]	Wheel speed [rpm]*)
Ceramics, minerals and crystals	> HV 800	LOW	0.005-0.15	3000
		LOW	0.005-0.20	3000
		HIGH	0.005-0.30	3000
		HIGH	0.005-0.30	2700
Sintered carbides and hard ceramics	> HV 800	MEDIUM	0.005-0.25	3000
		MEDIUM	0.005-0.25	2700
Extremely hard ferrous metals	> HV 500	MEDIUM	0.005-0.25	3000
Hard and very hard ferrous metals	HV 350-800	MEDIUM	0.05-0.30	1000-3000
		MEDIUM	0.05-0.30	1000-3000
Hard and very hard ferrous metals with larger dimensions	HV 350-800	MEDIUM	0.05-0.30	1000-3000
Soft and medium soft metals	HV 30-350	MEDIUM	0.05-0.30	1000-3000
		MEDIUM	0.05-0.30	1000-3000
Soft and ductile non ferrous metals	HV 70-400	MEDIUM	0.05-0.30	1000-3000
Plastics and very soft metals	< HV 100	MEDIUM	0.05-0.30	max. 1200

\*) To obtain low wheel wear and better surface quality, always use highest recommended wheel speed.

Common Cutting Questions

The following table shows the possible answers to a number of common cutting questions:

<b>Optimising the Cutting Results</b>	
<b>Objective</b>	<b>How to achieve it</b>
Better surface quality	Use lowest recommended feed speed and no specimen holder rotation.
Lower wheel wear	Use the lowest recommended feed speed and no specimen holder rotation. This is especially important when using resin bonded wheels and all abrasive cut-off wheels.
Problems with abrasive cut-off wheels?	Abrasive cut-off wheels should not be used outside their recommended feed speed range. At lower than recommended feed speeds they will produce irregularly cut surfaces. At higher feed speeds excessive wheel wear will occur, along with increased risk of wheel breakage.
Flatter samples	Use primarily low feed speeds, highest recommended wheel speed, largest possible flanges and no specimen holder rotation.  The initial kerf is especially critical. If the initial feed speed is too high the wheel will bend and start cutting at an angle. Such a cut will never end up flat.
Better parallelism	Use the lowest recommended feed speed.
Faster cutting	Orientate the sample so that the wheel will cut the smallest possible cross-section and then use maximum recommended feed speed.

Consumables

<b>Specification</b>	<b>Cat. No:</b>
<i>Corrozip</i> <i>Additive for Cooling Fluid</i> Environment friendly. To protect the machine from corrosion and to improve cutting and cooling qualities. 1 l 5 l	    49900045 49900046
<i>Cutting Fluid</i> Water free Cutting Fluid for cutting of water-sensitive materials 5 l	  49900030

## 4. Trouble-Shooting

Error Messages		
Display Message	Explanation	Action
Searching for XY-reference position	The machine searches for the initial reference positions	Wait for the search process to finish
Reference pos. not found, close cover	Accutom-5 was started with the cover open and did not search for the reference position	Close the cover, wait for the reference search to finish and press START $\diamond$ again
Cover open!	You are trying to start the cutting process or a dressing while the cover is open	Close the cover and press START $\diamond$ again
Process stopped by open cover	You have opened the cover during a cutting process or a dressing process	Close the cover and restart the interrupted process
Process running	You are trying to start two different processes at the same time	Wait for the previous process to finish and try again
Y-stop position outside cutting range	The set stop position is outside the maximum range in the Y-direction. (max. 105.0 mm)	Adjust Y-stop position.
Y-stop position less than Y-relative	You try to start a cutting process when the relative Y-position is larger than or equal to the Y-stop	Position the specimen holder correctly and/or change the Y-stop position.
Feed less than 90% of pre-set value	The force has exceeded the force limit (LOW, MEDIUM, HIGH) and the feed speed is automatically reduced	Decrease the feed speed or increase the force limit next time you cut the same material
Cutting stopped, feed too low	The feed speed has been reduced to less than 10% due to the force exceeding the force limit	Decrease the feed speed, increase the force limit or change the cut-off wheel and re-start the cutting process
Cutting method not selected	START $\diamond$ was pressed before a cutting method was selected	Select a method and press START $\diamond$ again
Dressing range exceeded	The parameter Y-distance is set to a value larger than the possible movement of the table	Adjust the Y-distance and press START $\diamond$ again

*Accutom-5  
Instruction Manual*

<b>Error Messages</b>		
<b>Display Message</b>	<b>Explanation</b>	<b>Action</b>
X-MIN sensor not activated	The inductive sensor cannot detect the end position of the X-table	Check for mechanical blockage of the table. If that is not the case, call a Struers service technician
X-MIN sensor not deactivated	The inductive sensor for the X-position cannot be deactivated	Check for mechanical blockage of the table. If that is not the case, call a Struers service technician
Encoder error	The Y-table does not move, or the encoder is defective	Check for mechanical blockage of the table. If that is not the case, call a Struers service technician
Cutting motor overloaded	The main motor has been exposed to a high load for a long time	Wait until the motor has cooled down "Cutting motor ready after overload" will be shown on the display. Adjust process parameters to avoid repetitions
<b>Machine Problems</b>	<b>Explanation</b>	<b>Action</b>
The recirculation pump is not working correctly	The water level in the recirculation tank is too low	Check that there is sufficient water in the recirculation tank
	The water outlet is clogged	Remove the cooling tubes and run the recirculation pump. Flush the cooling tubes with clean water before replacing
	The recirculation pump is contaminated with cutting residue	Flush the pump with clean water by gently pressing fresh water into the inlet tube


## 5. Technical Data

Subject	Specifications	
Cutting	<p><i>Cutting Speed:</i> 300-3000 rpm (adjustable in steps of 100 rpm)  <i>Feed Speed:</i> 0.005-3.000 mm/s (adjustable in steps of 0.005 mm/s)  <i>Max. Positioning Speed:</i> Y = 13 mm/s, X = 13 mm/s  <i>Force Limits:</i> LOW: approx. 20N; MEDIUM: approx. 40N; HIGH: approx. 60N</p>	
Positioning Range	<p>Y direction: 105 mm (precision 0.1 mm)  X direction: 60 mm (precision 0.005 mm)</p>	
Sample Size	<p><i>Max. Length of Cut-off Sample:</i> 30 mm, 140 mm at <math>\varnothing 20</math> mm  <i>Max. length of Sample to be clamped:</i> 225 mm  <i>Max. Sample Cross Section:</i>  <math>\varnothing 127</math>mm cut-off wheel and <math>\varnothing 42</math>mm flange: <math>\varnothing 40</math>mm (without rotation)  <math>\varnothing 80</math>mm (with rotation)  <math>\varnothing 152</math>mm cut-off wheel and <math>\varnothing 42</math>mm flange: <math>\varnothing 50</math>mm (without rotation)  <math>\varnothing 100</math>mm (with rotation)</p>	
Cut-off Wheels	<p><i>Diameter:</i> <math>\varnothing 75</math> to <math>\varnothing 152</math> mm / 3 to 6"  <i>Max. Thickness:</i> 2 mm / 0.08"  <i>Hole:</i> <math>\varnothing 12.7</math> mm / 0.5"</p>	
Recirculation Cooling Unit	<p><i>Contents:</i> 4 l  <i>Flow:</i> 800 ml/min</p>	
Motor	370 W at 3000 rpm                      Continuous Torque 0.8 Nm	
Noise Level	Approx. 65 dB (A) measured at idle running, at a distance of 1.0 m / 39.4" from the machine.	
Altitude	Min. 50 m. below sea level. Max 3000 m above sea level	
Surrounding temperature	5-40°C/41-104°F	
Humidity	0-95% RH non condensing	
Software and Electronics	<p><i>Display:</i> 8 x 40 characters  <i>Controls:</i> touch pad  <i>Database:</i> 20 cutting methods  <i>PCB fuses:</i> 3.15A , 1,6A SUB mini fuses</p>	
Power Supply	50-60Hz (max. load 2.7A)                      1/3 x 220 - 240V	
Dimensions and Weight	Width	510 mm / 20"
	Depth	700 mm / 27.6"
	Height	270 mm / 10.6"
	Weight	45 kg / 99 lbs
Safety Standard	Please refer to the Declaration of Conformity	

# 6. Menu Overview

Accutom-5  
Menu Structure


# Quick Reference Guide

## Loading a Cutting Method

- Press Esc until the Main Menu appears.
- Select Cutting Methods and press Enter ↵.
- Select the method you want to use and press Enter ↵.

## Positioning the Sample

The “heavy” arrows on the display indicate the direction in which the sample holder is moved, for example:

- Press one of the positioning keys ▼ ▲ to change to the Y-Positions screen.
- Place the specimen holder with the sample in the specimen holder head and clamp it.

*If you already know the cutting length (i.e. the diameter of the sample):*

- Position the sample correctly in front of the cut-off wheel using the positioning keys.
- Press one of the positioning keys ▼ ▲ to change to the Y-Positions screen.
- Press F1:RESET to reset the relative Y-position.
- Press Enter ↵ and use the MENU arrows ⏪⏩⏴⏵ to adjust the Y-stop value. Press Enter ↵ again to accept the new value, or Esc. to cancel.

*If you do not know the cutting length in advance:*

- Position the sample correctly in front of the cut-off wheel using the positioning keys.
- Press F1:RESET to reset the relative X-position.
- Press one of the positioning keys ▼ ▲ to change to the Y-Positions screen.
- Press F1:RESET to reset the relative Y-position.
- Move the sample to the right using the positioning key ► so that it can pass along the cut-off wheel.
- Move the sample along the cut-off wheel to the desired stop position using the positioning key ▲. Press F2:SET STOP to record the stop position.
- Press F3:MOVE TO REL: ZERO to move the sample back to the initial Y-position.
- Press one of the positioning keys ◀▶ or ◀ to change to the X-Positions screen.
- Press F3:MOVE TO REL: ZERO to move the sample back to the initial X-position.

## Starting the Process

- Press START ⬡ to start the cutting process.

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
```
INSERT METHOD
Do you want to copy from buffer to
method 4. EMPTY METHOD      ?
ENTER: YES                  ESC: NO
```

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```
SELECT CUTTING METHOD ↓
1.UNNAMED METHOD      2.EMPTY METHOD
3.EMPTY METHOD        4.UNNAMED METHOD
5.EMPTY METHOD        6.EMPTY METHOD
7.EMPTY METHOD        8.EMPTY METHOD
9.EMPTY METHOD        10.EMPTY METHOD
11.EMPTY METHOD       12.EMPTY METHOD
F1: COPY  F2: INSERT F3: RESET F4: RENAME
```

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```
RESET METHOD
Do you want to reset this method ?
4. UNNAMED METHOD
ENTER: YES          ESC: NO
```

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SELECT CUTTING METHOD ↓	
1.UNNAMED METHOD	2.EMPTY METHOD
3.EMPTY METHOD	4. <b>UNNAMED METHOD</b>
5.EMPTY METHOD	6.EMPTY METHOD
7.EMPTY METHOD	8.EMPTY METHOD
9.EMPTY METHOD	10.EMPTY METHOD
11.EMPTY METHOD	12.EMPTY METHOD

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UNNAMED METHOD	EDIT NAME
Text: UNNAMED METHOD	
Press ENTER to accept the text, or press "↓" to select new character.	

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
↓  

```
NAME CHANGED
The name is changed, do you want to keep
the old name, or accept the new name ?
ENTER:ACCEPT NEW NAME  ESC:KEEP OLD NAME
```

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*Accutom-5*  
*Spare Parts and Diagrams*

---

**Always state *Serial No* and *Voltage/frequency*  
if you have technical questions or when ordering spare parts.**

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The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2012.

**Struers A/S**  
Pederstrupvej 84  
DK-2750 Ballerup  
Denmark  
Telephone +45 44 600 800  
Telefax: +45 44 600 801

---

## **Spare Parts and Diagrams**

	<b>Table of contents</b>	<b>Drawing</b>
<b>Drawings</b>	Cabinet .....	14940021D
	Water Cock.....	14940044B
	Recirculation Tank .....	14940041G
	Bottom .....	14940051U
	Back plate.....	14940054F
	Various Parts	
<b>Diagrams</b>	Wiring Diagram .....	15013452E
	Block Diagrams.....	15013052E

The drawings are not to scale.  
Some of the drawings may contain position numbers  
not used in connection with this manual.


*Accutom-5*  
*Spare Parts and Diagrams*

## Spare Part List for Accutom-5

<b>Drawing</b>	<b>Pos.</b>	<b>Spare Part</b>	<b>Cat no.</b>
14940021		<b>Cabinet</b>	
	15	Seal for cover	R4940218
	25	Safety switch	2SS00012
	70	Cover	14940240
14940044		<b>Water cock</b>	
	10	Pipe branch	14940420
	15	O-ring 7.65-1.78 (10 pcs)	RIO17812
	20	O-ring 10.82-1.78 (10 pcs)	130MP057
	30	Socket for water nozzle	R4940422
	40	End piece for water nozzle	14940428
	60	Water nozzle with internal thread for stopper	R4940421
	70	Stopper for water nozzle with thread	R4940429
14940041		<b>Recirculation Tank</b>	
	10	Bottom plate	14949005
	20	Fitting for Drawer	2JX04230
	50	Front plate for Recirculation tank	14949006
	60	Foil for recirculation drawer	15010315
	90	Tension spring	R4940416
14940051		<b>Bottom</b>	
	20	Rubber foot, with distance (pos. 30) and screw	R4910140
	910	Silicone suction hose $\varnothing 8/\varnothing 12$ (3 pcs)	14942908
14940054		<b>Back Plate</b>	
	50	Fuse glass tube 6.3A (10 pcs)	315MP134

## **Spare Part List for Accutom-5**

<b>Drawing</b>	<b>Pos.</b>	<b>Spare Part</b>	<b>Cat no.</b>
		<b>Various Parts</b>	
		Flange ø42, back	14940527
		Flange ø42, front	14940529
		Flange screw M5	427MP012
		Backstop with 1 hole	14278017
		Backstop with 2 holes	14278016
		Pressure spring	R2600718
		Grate	14940220
		Grate, large	R4940222
		Backstop pin for cut-off wheel flange	14940701
		2 pin spanner for cut-off wheel mounting	427MP085
		Allen key 2 (10 pcs)	RGR00020
		Allen key 2,5 (10 pcs)	RGR00025
		Allen key 3 (10 pcs)	RGR00030
		Allen key 4 (5 pcs.)	RGR00040
		Allen key 5 (5 pcs)	RGR00050
		Set of screws for Specimen Holder for general use (Cat. no. 04946903)	14942906
		Debris collection tray	14940221

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













Rev B: Motor Control. S1, S2, S3 added.  
Rev C: Motor Control. Connections changed.  
Rev D: Z1 changed to Z2.  
Rev E: Z3 and X22 removed.

WALDORF & WALTERS ALLE 176  
D-30900 WEDDINGEN  
TELEFON +49 5130 3000

BLOCK AND WIRING DIAGRAM: ACCUTOM 5/50  
230V BLOCK- / WIRINGDIAGRAM

FILE NAME: 5013052.DSN  
Rev: A1  
15013052  
E

Wednesday, November 21, 2012  
KFC/KFC

English

## Declaration of Conformity

**Manufacturer** Struers A/S  
Pederstrupvej 84  
DK-2750 Ballerup, Denmark  
Telephone +45 44 600 800

Herewith declares that

<i>Product Name:</i>	Accutom-5
<i>Type No:</i>	494
<i>Machine Type:</i>	Precision cut-off machine

is in conformity with the provisions of the following directives:

**Safety of Machinery** 2006/42/EC according to the following standard(s):  
EN ISO 12100:2011, EN60204-1:2006/AC:2010.

**Low Voltage Directive** 2006/95/EC according to the following standard(s):  
EN60204-1:2006/AC:2010.

**EMC-Directive** 2004/108/EC according to the following standard(s):  
EN61000-6-1:2007, EN61000-6-3:2007/A1:2011.

**RoHS** 2011/64/EC according to the following standard(s):  
EN 50581:2012.

**Supplementary Information** The equipment complies with the American standards:  
FCC part 15, subpart B and UL508.

**The above has been declared according to the global method, module A**

Date: 14.01.2014

  
Christian Skjold Heyde,  
Vice President, R & D and Production, Struers A/S

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... and grinding machine

### Common features of Accutom-5/-50

- Positioning accuracy of 5  $\mu\text{m}$
- Pre-set, constant feed speed (0.005 to 3 mm/sec)
- Adjustable force limit
- Sample rotation or oscillation
- Variable cut-off wheel speed up to 3,000 rpm
- Menu control and digital read-out
- Built-in cut-off wheel list
- Database with user-definable cutting methods
- Built-in recirculation cooling unit
- Enclosed cutting chamber with safety switch

### Common application areas:

Precise and deformation-free cutting of:

- Metals
- Electronic components
- Ceramics
- Composites
- Crystals
- Sintered carbides
- Minerals
- Fibre reinforced materials
- Biomaterials (teeth and bones)

### Additional features of Accutom-50

- Serial cutting of slices (Multi-cut)
- Grinding of thin sections
- Target grinding in components
- Variable cut-off wheel speed up to 5,000 rpm
- Database for both cutting and grinding methods
- Improved menu control

### Additional application areas:

- Precise serial cutting of all materials
- Thin sections of minerals
- Thin sections of ceramics
- Target grinding in components

*Oscillating movement is used for cutting oblong samples, thus reducing the cutting time considerably*

Accutom-50 is based on the known and well tested technology of Accutom-5. Sample positioning within 5 µm combined with constant feed speed and sample rotation or oscillation represents the utmost in precision cutting. The addition of fully automatic serial cutting and precision grinding facilities, together with even higher wheel speeds, makes Accutom-50 the most universal tool in precision manufacturing of samples.

### Working principle

Both Accutom-5 and Accutom-50 are built on the same, precisely manufactured, cast frame. The use of step motors to control the movement of the sample in both X- and Y-directions enables very precise positioning.

### Positioning accuracy of 5 µm

The high resolution of the step motor makes it possible to position the sample with an extraordinary precision of 5 µm. The total movement of the sample in the X-direction is 60 mm.

### Pre-set, constant feed speed

The feed speed can be pre-set and controlled in the range of 0.005 to 3 mm/sec. This wide range allows even the most delicate materials to be cut. Both the pre-set and the actual speeds are shown on the display. If the preset force limit is reached, the feed speed will be reduced and the process stopped if the speed has fallen below a specified limit.

### Adjustable force limit

To protect thin cut-off wheels or fragile samples, the force limit can be set to three different levels: Low, Medium and High.

### Sample rotation or oscillation

Sample rotation or oscillation reduces the contact area between cut-off wheel and sample to a minimum. This results in a decrease in frictional heat and allows fast cutting of extremely hard materials.

### Variable cut-off wheel speed up to 3000 rpm


The motor of Accutom-5 has variable cut-off wheel speeds from 300 up to 3000 rpm in steps of 100 rpm. This allows both high and low speed cutting and enables a more accurate choice of wheel speed for different cut-off wheel diameters.

### Menu control and digital readout

Accutom-5 is operated from a touch pad on the ergonomically designed, slanting front. All parameters, both for positioning and cutting, are presented on the LC-Display.

### Built-in cut-off wheel list

A list of all Struers cut-off wheels is stored in an internal database together with their name, respective thickness and force limit.

*Accutom-50 can be programmed to automatically cut series of slices from the same samples*

### Database with user-definable cutting methods

After establishing a method for a specific material, all settings can be saved in a database and recalled at any time. Accutom-5 has a storage capacity of ten methods.

### Built-in recirculation cooling unit

With a capacity of four litres, and a very high through-put, the recirculation cooling unit ensures a constant low temperature during the cut, and an efficient removal of cutting debris.


### Enclosed cutting chamber with safety switch

The fully transparent cover is equipped with a safety switch to prevent the cut-off wheel from starting whilst the cover is open. The sample can be inserted and positioned but the cut-off wheel itself can only be activated when the cover has been closed properly, thus providing maximum safety for the operator.

### Accutom-50 provides the following additional features:


#### Serial cutting of slices (Multi-cut)

Accutom-50 can be programmed to automatically cut series of slices from the same sample. The advanced software will also compensate for the wheel thickness and for slight drift of the cut-off wheel. This results in very precise sections with an identical thickness, as required in failure analysis or the examination of biomaterials.

*After the process has been started, the display changes and the relevant parameters are shown*

*Grinding on Accutom-50 can be carried out in two different ways: The grinding operation can be stopped at a specified position, or a well-defined amount of material can be removed*

### **Grinding of thin sections**

The thickness ground during each different operation associated with a thin section can be easily adjusted or even saved in different grinding methods. After grinding the glass plate to the precise thickness, the position of the sample holder is set to zero. Thereafter the thickness of the thin section is entered directly.

### **Target grinding in components**

When the precise position of a component within a sample is known, Accutom-50 can be programmed to grind down to precisely that position. The position can be adjusted in steps of 5 µm.

### **Variable cut-off wheel speed up to 5000 rpm**

The variable cut-off wheel speed can be adjusted from 300 up to 5000 rpm in steps of 100 rpm. This grants even higher precision as the drift of the cut-off wheel is reduced because of the higher speed. Also, the feed speed can be increased without experiencing any deterioration in the high quality of the cut surface.

### **Database for both cutting and grinding methods**

Accutom-50 is equipped with a database capable of storing 20 cutting and 10 grinding methods.

### **Improved menu control**

Subsequent menus for the selection of cutting or grinding parameters greatly facilitate the operation. Together with the larger display the new menu control improves the user-friendliness of Accutom-50.

### **Cut-off wheels**

As Accutom-50 works with programmable constant feed speeds, all types of cut-off wheels, including abrasive wheels, can be used.

Struers offer a large variety of Diamond, CBN, Al<sub>2</sub>O<sub>3</sub> and SiC cut-off wheels covering all possible applications. Also, saw blades for the cutting of very soft metals or plastics can be used on Accutom-5/50 due to the cover with the safety switch.

### **Cup wheels**

Two different types of cup wheels, resin and metal bonded, are available for grinding on Accutom-50.

### **Specimen holders**

All specimen holders for Accutom-5/50 are provided with a dovetail plate. This means that the samples can be clamped into the sample holder outside the cutting chamber, and then easily inserted into the dovetail fixture of Accutom-5/50.

*A wide range of specimen holders can be used on Accutom machines*

## Technical data

Cutting	Accutom-5	Accutom-50
	Cutting speed: 300-3000 rpm (adjustable in steps of 100 rpm) Feed speed: 0.005-3.000 mm/s (adjustable in steps of 0.005 mm/s) Max.pos.speed: Y=13 mm/s X=10 mm/s Force limits: low: approx. 20 N medium: approx. 40 N; high: approx. 60 N	Cutting speed: 300-5000 rpm (adjustable in steps of 100 rpm) Feed speed: 0.005-3.000 mm/s (adjustable in steps of 0.005 mm/s) Max.pos.speed: Y=13 mm/s X=10 mm/s Force limits: low: approx. 20 N medium: approx. 40 N; high: approx. 60 N
Pos. range	Y direction: X direction:	105 mm (precision 0.1 mm), 60 mm (precision 0.005 mm)
Sample size	Max length of cut-off sample: Max length of sample to be clamped: Max sample cross section: ø127 mm cut-off wheel and ø42 mm flange: ø152 mm cut-off wheel and ø42 mm flange:	30 mm, 140 mm at ø20 mm 225 mm ø40 mm (without rotation) ø80 mm (with rotation) ø50 mm (without rotation) ø100 mm (with rotation)
Recirculation cooling unit	Contents: 4 l Flow: 800 ml/min	Contents: 4 l Flow: 800 ml/min
Motor	370 W at 3000 rpm	550 W at 5000 rpm
Software and electronics	Display: 4 x 40 characters Controls: touch pad Database: 10 cutting methods	Display: 4 x 40 characters Controls: touch pad Database: 20 cutting and 10 grinding methods
Power supply	1-3 x 200-240 V / 50/60 Hz	
Dimensions and weight	Width: 510 mm / 20", Depth: 700 mm / 27.6", Height: 270 mm / 10.6", Weight: 45 kg / 99 lbs	

## Specifications

	Cat. no.
<b>Accutom-5</b> Precision cut-off machine with variable speed. Automatic feeding with electronic control of feed speed. Motorised positioning system and digital read-out, with built-in sample rotation and oscillation. Complete with recirculation cooling unit, flange sets 42 mm dia. and specimen holder (04946903). Cut-off wheels, additive for cooling fluid and additional specimen holders are ordered separately.	04946227
<b>Accutom-5</b> Precision cut-off machine with variable speed. Automatic feeding with electronic control of feed speed. Motorised positioning system and digital read-out. Complete with recirculation cooling unit, flange sets 42 mm dia. and specimen holder (04946903). Without built-in sample rotation and oscillation. Cut-off wheels, additive for cooling fluid and additional specimen holders are ordered separately.	04946127
<b>Accutom-50</b> Precision cut-off and grinding machine with variable speed. Equipped with Multi-cut feature. Automatic feeding with electronic control of feed speed. Motorised positioning system and digital read-out, with built-in sample rotation and oscillation. Complete with recirculation cooling unit, flange sets 42 mm dia. and specimen holder (04946903). Cut-off wheels, additive for cooling fluid, cup wheels, additional specimen holders, and the flange set for cup-wheels are ordered separately.	05016227
<b>Accessories</b>	
<b>Specimen holders</b>	
Three jaw-chuck	05036919
For general use. Vice type with max. opening 60 mm	04946903
For round or square specimens. Teardrop type with max. opening ø 40 mm / ø 1½"	04946904
For irregular specimens, with 7 screws. Max. width 40 mm / 1½"	04946905
With goniometer	04276911
For adhering specimens	04276912
With ceramic vacuum chuck for thin sections	04276913
For small specimens. Vice type	04276915
For long specimens. With double parallel vice	04946909
For cylindrical specimens, max. dia. 40 mm	04946907
<b>Joints to be mounted between the dovetail and the specimen holder</b>	
Tilting joint with max. angle ±10°	04946906
Angling joint with max. angle +30/-90°	04946908
<b>Base Plate</b>	
With dovetail. For mounting of other types of specimen holders	04276914
<b>Flange Sets</b>	
For cut-off wheels, 65 mm dia., to achieve highest precision	04946902
For cup wheels, only to be used on Accutom-50	05016901
<b>Additive for Cooling Fluid</b>	
For Recirculation Cooling Unit. To protect the machine from corrosion, and to improve cutting and cooling qualities	1 l 49900045 5 l 49900046
See special brochure for more information about cut-off wheels.	

Struers' equipment is in conformity with the provisions of the applicable International Directives and their appurtenant Standards. (Please contact your local supplier for details)

Struers' products are subject to constant product development. Therefore, we reserve the right to introduce changes in our products without notice.

**Struers A/S**  
Pederstrupvej 84  
DK-2750 Ballerup, Denmark  
Phone +45 44 600 800  
Fax +45 44 600 801  
struers@struers.dk  
www.struers.com

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### CANADA

**Struers Ltd.**  
7275 West Credit Avenue  
Mississauga, Ontario L5N 5M9  
Phone +1 905-814-8855  
Fax +1 905-814-1440  
info@struers.com

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**Struers GmbH**