Unitom-5

Instruction Manual

Manual No.: 15337001

Date of Release 15.11.2005

Unitom-5
Instruction Manual

Table of Contents	Page
User's Guide	1
Reference Guide.....	19
Quick Reference Guide.....	43

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the Serial No. and Voltage on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2004.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Fax +45 44 600 801

Unitom-5 Safety Precaution Sheet

To be read carefully before use

1. The operator(s) should be fully instructed in the use of the machine and its cut-off wheels according to the Instruction Manual and the instructions for the cut-off wheels.
2. The machine must be placed on a safe and stable support table. All safety functions and guards of the machine must be in working order.
3. Use only intact cut-off wheels. The cut-off wheels must be approved for min. 50 m/s.
4. Observe the current safety regulations for handling, mixing, filling, emptying and disposal of the additive for cooling fluid.
5. The workpiece must be securely fixed in the quick-clamping device or the like. Large or sharp workpieces must be handled in a safe way.
6. Do not touch anything inside the cutting chamber while positioning the cut-off wheel with the position button.
7. Struers recommends the use of exhaustion as the materials to be cut may emit harmful gasses or dust.

If optional rubber curtains are needed in order to cut long, protruding workpieces:

1. Check with your local authorities whether rubber curtains are approved as protection in your country.
2. Do not put the hands through the rubber curtain during operation of the machine.
3. Protruding workpieces should be shielded or marked.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

User's Guide

Table of Contents	Page
1. Getting Started	
Checking the Contents of the Packing Box	2
Recirculation Cooling Unit.....	2
Placing Unitom-5	2
Getting Acquainted with Unitom-5	3
Noise Level.....	4
Supplying Power	4
Direction of the Cut-off Wheel	4
Connecting the Recirculation Unit	5
Software Settings.....	6
Configuration Menu	6
Setting the Language	7
2. Basic Operations	
Using the Controls.....	8
Front Panel Controls of Unitom-5	8
Groups of Keys.....	8
Acoustic Signals.....	8
Front Panel Controls	9
Cooling Valve Positions.....	10
Reading the Display	11
Reading the Display during the Cutting Process	11
Sleep Mode.....	11
Changing/Editing Values	12
Numeric Values	12
Alphanumeric Values.....	13
Changing the Cut-off Wheel.....	14
Clamping the Workpiece	14
Positioning the Cut-off Wheel	14
Setting the Feed Rate before Cutting.....	15
Starting the Cutting	16
Fast down.....	16
Setting the Feed Rate during Cutting.....	16
Stopping the Cutting	17
Autostop	17
Position Stop.....	17
STOP  Key	17
Flushing the Cutting Chamber.....	18

1. Getting Started

Checking the Contents of the Packing Box

In the packing box you should find the following parts:

- 1 Fork spanner (30 mm), for changing the cut-off wheel
- 1 Grease gun, for lubrication of tilting frame bearings
- 2 Cover plugs for threads for cutting tables
- 1 Set of Instruction Manuals

Recirculation Cooling Unit

- 1 PVC hose 1/2", 2 m
- 1 Outlet hose, 2 m
- 1 Drain elbow pipe
- 5 Hose clamps

Placing Unitom-5


Unitom-5 should be placed on the optional Struers table (LABUL), or a table which is able to carry a weight of min. 200 kg.

Adjust the table height using the screws in the table legs. It is important that the table and Unitom-5 are level. Place Unitom-5 on the table. If the machine is not horizontal, you must adjust the height for each corner, as follows:

- Loosen the upper 13 mm nut under the corner screw.
- Adjust the lower 8 mm nut with a spanner.
- Fasten the upper 13 mm nut again.
- Repeat the adjustment for each of the other three corners until Unitom-5 is levelled.

**Getting Acquainted with
Unitom-5**

Take a moment to familiarise yourself with the location and names of the Unitom-5 components.

- ① Front panel controls
- ② Cutting Table
- ③ Quick-clamping device
- ④ Table LABUL (optional)
- ⑤ Recirculation cooling unit
- ⑥ Removable plate. Can be replaced by a rubber curtain (optional)
- ⑦ Main switch
- ⑧ Protection guard with window
- ⑪ Cooling valve

Noise Level

Approx. 76 dB (A) measured at idle running, at a distance of 1.0 m/39.4" from the machine.

Supplying Power

IMPORTANT

Check that the mains voltage corresponds to the voltage stated on the type plate on the side of the machine.

- Unitom-5 is factory mounted with an electric cable. Mount a plug on the cable according to the following:

50 hz		60 hz	
Phase	Black and brown	L ₁	Black
Earth	Yellow/green	L ₂	Red
		L ₃	Orange
		Earth	Green

Direction of the Cut-off Wheel

Check that the cut-off wheel rotates in the direction indicated on the label attached to the guard of Unitom-5. If the direction of rotation is incorrect:

50 hz	60 hz
Switch two of the phases	Switch phases L ₁ and L ₂

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
Software Settings

Switch on the power at the main switch. The following display will appear briefly:

Afterwards, the display will change to the screen for cutting.

 Press ESC to come to the MAIN MENU.

Configuration Menu

The MAIN MENU is the highest level in the menu structure. From here you can go to cutting or to the configuration menu.

 Press MENU DOWN  to select CONFIGURATION.

*Unitom-5
Instruction Manual*

 Press ENTER ↵ to activate the CONFIGURATION menu.

CONFIGURATION		POSITION	
Contrast	26	Wheel	100.0mm
Language	ENGLISH	Stop	223.0mm
Unit	mm		
Autostop	On		
Return to	Start		
Feed reduction	On		
F1:Default			


Setting the Language

 Press MENU UP/DOWN  to select Language.

 Press ENTER ↵ to activate the Language pop-up menu.

 Press MENU UP/DOWN  to select the language you prefer.


 Press ENTER ↵ to accept the language.

The CONFIGURATION Menu now appears in the language you have chosen.

2. Basic Operations

Using the Controls

Front Panel Controls of Unitom-5


Groups of Keys

- A** START /STOP  of Unitom-5
- B** Position of the cut-off wheel
- C** Programming and monitoring functions

Acoustic Signals

-  A short beep when a key is pressed indicates that the command has been accepted.
-  A long beep indicates that the key is inactive at the moment.
-  Three long beeps indicate that the cutting process is finished.











Location of the Main Switch

Please Note...




The contents of the program memory are not lost when the main switch is turned off.

Front Panel Controls


Name	Key	Function	Name	Key	Function
① MAIN SWITCH		The main switch is located to the right hand side of Unitom-5. Turn clockwise to switch on the power.	⑧ ESC		Leaves the present menu or aborts functions/changes.
② POWER		Lights when main power is switched on.	⑨ ENTER		Accepts a marked parameter value or chooses a menu.
③ START		Starts the cutting process.	⑩ FUNCTION KEY		Controls for various purposes. See the bottom of the individual screens.
④ STOP		Stops the cutting process.	⑩ FUNCTION KEY		Controls for various purposes. See the bottom of the individual screens.
⑤ ⑥ POSITION		Moves the cut-off wheel up or down. Keep the key pressed to move the cut-off wheel.	⑩ FUNCTION KEY		Controls for various purposes. See the bottom of the individual screens.
⑦ MENU UP		Moves the cursor up or increases the parameter value.	⑩ FUNCTION KEY		Controls for various purposes. See the bottom of the individual screens.
⑦ MENU DOWN		Moves the cursor down or decreases the parameter value.	⑫ EMERGENCY STOP		- Push the red button to stop. - Pull the red button to release.

⑪ COOLING VALVE: See next page!

Cooling Valve Positions


<p>⑪ COOLING VALVE</p>	<p>The cooling valve for cooling/cleaning is located in the back of the cutting chamber to the left.</p>
	<p>Turn the knob counter-clockwise to vertical position for cooling of cut-off wheel.</p>
	<p>Turn the knob to regulate the flow while flushing the cutting chamber.</p>
	<p>Turn the knob clockwise to vertical position for flushing the cutting chamber with maximum flow.</p>

Location of Cooling Valve for cooling/cleaning

Reading the Display

The display can show various kinds of information, for example parameters for an ongoing process. During the cutting process the screen could look as the following example:

- A Heading.
- B Inverted text: cursor position.
- C Function key options.

Reading the Display during the Cutting Process

Feed The feed rate of the cut-off wheel pre-set before starting the cutting.

Load Load indicator of the cutting motor. The moving black column shows the status of the cutting motor:

Inside the white area: The cutting motor can be loaded continuously.

In the grey overload area: The cutting motor can only be intermittently loaded. The built-in thermo fuse might shut off the cutting motor.

Wheel Position. The actual position of the cut-off wheel counted from the top position of the cut-off wheel.

Stop Position The pre-set stop position counted from the top position of the cut-off wheel.

Sleep Mode

To increase the lifetime of the display, the backlight is automatically switched off if Unitom-5 has not been used for 15 min. Press any key to re-activate the backlight.


Please Note

The sample screens in this Instruction Manual show a number of possible texts. The actual display screen may differ from the examples in the Instruction Manual.

Changing/Editing Values

Depending on the type of value, there are two different ways of editing.

Numeric Values

Press ENTER ↵, to edit the value.

Two square brackets [] appear around the value.

Press MENU UP/DOWN ⏮ ▲ ▼
⏭ to increase or decrease the numeric value.



Press ENTER ↵, to accept the new value.

Pressing Esc, aborts the changes, preserving the original value.

Alphanumeric Values


CONFIGURATION		POSITION	
Contrast	26	Wheel	100.0mm
Language	ENGLISH	Stop	223.0mm
Unit	mm		
Autostop	On		
Return to	start		
Feed reduction	On		
F1:Default			


↓
 Press MENU UP/DOWN  to select the alphanumeric value you want to change

↓
 Press ENTER , to edit the value.

A pop-up menu appears.


↓
 Press MENU UP/DOWN  to select the correct value.

↓
 Press ENTER , to accept the new value and to return to the previous screen.

Pressing Esc, preserves the original value, returning to the previous screen.

Changing the Cut-off Wheel

- Open the protection guard and swing the right cooling tube to the right.
- Press the black locking knob on the left-hand side of the cut-off wheel, turning the cut-off wheel until the spindle lock clicks.
- Remove the nut with a fork spanner. Remove the flange and the old cut-off wheel.

IMPORTANT
The spindle in Unitom-5 is right-hand threaded.

- Mount the new cut-off wheel.

IMPORTANT
Conventional cut-off wheels based on Al₂O₃/SiC abrasives should be placed between two cardboard washers, to protect the cut-off wheel.
For maximum precision with diamond or CBN cut-off wheels, do not use cardboard washers.

- Mount the flange and the nut. Tighten carefully.
- Swing the right cooling tube back into position.

Clamping the Workpiece

- Place the workpiece between the clamp and the back stop at the right hand cutting table.
- Push the clamp towards the workpiece and lock the quick-clamping device with the locking handle.

Positioning the Cut-off Wheel

- After clamping, you may lower the cut-off wheel to a position at least 1 mm above the workpiece, to test the correct position of the cut on the workpiece.
- Press POSITION DOWN ▼ to lower the cut-off wheel.
- The cut-off wheel stops when you stop pressing the key.
- You can follow the position of the cut-off wheel in the display:

CUTTING		Load	POSITION	
Feed	0.50mm/s	█	Wheel	100.0mm
		█	Stop	223.0mm
SET STOP				
F3:ACT. F4:MAX				

- Leave the cut-off wheel in position before you press START ◊.


Setting the Feed Rate before Cutting


CUTTING		POSITION	
Feed	0.50mm/s	Wheel	100.0mm
		Stop	223.0mm
		SET STOP	
		F3:HCT. F4:MAX	

- ↓
-  Press ENTER ↵, to edit the value.
- Two square brackets [] appear around the value.
- ↓
-  Press MENU UP/DOWN ⏮ ▲ ▼ to increase or decrease the numeric value.
- ↓
-  Press ENTER ↵, to accept the new value.

Pressing Esc, aborts the changes, preserving the original value.

Starting the Cutting

CUTTING		POSITION
Feed	0.50mm/s	Wheel 106.3mm
	Load 	Stop 223.0mm
F1:Fast down		SET STOP
		F3:ACT. F4:MAX

- Before starting, lower the cut-off wheel into correct position above the workpiece.
- Carefully close the protecting guard.
- Press START . The cut-off wheel starts rotating and the cooling water starts running. The cut-off wheel moves down with the pre-set feed rate.



Fast down

You may increase the feed rate to twice the pre-set value between pressing start and the cut-off wheel physically touching the workpiece.

- Press the function key F1. The cut-off wheel moves down fast until it touches the workpiece. From that minute, Unitom-5 automatically changes the feed rate back to the pre set value.

Setting the Feed Rate during Cutting

CUTTING		POSITION
Feed	0.50mm/s	Wheel 141.6mm
	Load 	Stop 223.0mm
		SET STOP
		F3:ACT. F4:MAX

- ↓
-  Press ENTER ↵, to edit the value.
- Two square brackets [] appear around the value.
- ↓
-  Press MENU UP/DOWN   to increase or decrease the numeric value.
- ↓
-  Press ENTER ↵, to accept the new value.
- Pressing Esc, aborts the changes, preserving the original value.

Stopping the Cutting

Autostop

- Unitom-5 automatically stops cutting when the workpiece has been cut through (default setting). The cut-off wheel moves to start position (default setting). The cut-off wheel stops rotating and the cooling water stops.


Position Stop

- To be used if you want to stop the cutting at a certain depth (please see Advanced Operations, for further information).

STOP Key

- You can always stop the cutting process by pressing the STOP  key. The cut-off wheel moves to its start position (default setting). The cut-off wheel stops rotating and the cooling water stops.

Flushing the Cutting Chamber


To ensure a longer lifetime of your Unitom-5, Struers strongly recommends that you clean the cutting chamber with the cleaning hose every day.

- Point the cleaning hose towards the bottom of the cutting chamber.
- Turn the cooling valve to cleaning position (see sketch in table) and adjust the flow as needed.
- Clean the cutting chamber thoroughly.


NB! Do not clean the protection guard with the cleaning hose; use a damp cloth instead.

Cooling Valve


The cooling valve for cooling/cleaning is located in the back of the cutting chamber to the left.

Turn the knob counter-clockwise to vertical position for cooling of cut-off wheel.

Turn the knob to regulate the flow while flushing the cutting chamber.


Turn the knob clockwise to vertical position for flushing the cutting chamber with maximum flow.

- Turn the cooling valve to cooling position.
- Leave the protection guard open to let the cutting chamber dry completely.

AVOID RUST!

Leave the protection guard open to let the cutting chamber dry completely.

Reference Guide

Table of Contents	Page
1. Advanced Operations	
Configuration Menu	21
Display contrast.....	22
Language.....	22
Unit	22
Autostop	22
Return position	22
Feed reduction	22
F1 Default value	22
Adjusting the Back Stop.....	23
Clamping Irregular Workpieces.....	23
Using the Continue Function.....	23
Opening the Protection Guard before using the Continue Function	23
Using Autostop	24
Using Fixed Stop.....	25
Using Fixed Stop together with Autostop	25
Setting the Stop Position	25
The Feed Reduction Feature.....	26
Optimising the Cutting Results	27
Connection to External Exhaust System.....	28
Mounting the Rubber Curtain	28
Explanation of Safety Factors.....	28
2. Accessories	29
3. Consumables	
Other Consumables	30
4. Trouble-Shooting	31

5. Maintenance

Daily Service 35
 Checking the Recirculation Unit 35
Weekly Service..... 35
Monthly Service 36
 Replacing the Cooling Water 36
 Lubrication of Bearings for Tilting Frame..... 37
Maintenance of Cutting Tables..... 38
Maintenance of Cut-off Wheels 38
 Storing of Bakelite Bonded Al₂O₃ Cut-off Wheels 38
 Maintenance of Diamond and CBN Cut-off Wheels 38
Maintenance of Clamping Devices..... 38

6. Technical Data 39

7. Menu Overview 41

1. Advanced Operations


Configuration Menu


- Press Esc until you reach the Main Menu.

MAIN MENU	POSITION
CUTTING	Wheel 100.0mm
CONFIGURATION	Stop 223.0mm
	SET STOP
	F3:ACT. F4:MAX

Press MENU DOWN  to select CONFIGURATION

Press ENTER  to activate the CONFIGURATION menu

↓

CONFIGURATION	POSITION
Contrast 26	Wheel 100.0mm
Language ENGLISH	Stop 223.0mm
Unit mm	
Autostop On	
Return to Start	
Feed reduction On	SET STOP
F1:Default	F3:ACT. F4:MAX

<i>Display contrast</i>	The contrast settings of the display can be adjusted to suit the individual preferences.
<i>Language</i>	The language can be set to either English, German or French.
<i>Unit</i>	The position of the cut-off wheel and the Automatic Cutting Table (optional) can be displayed in either mm or inches.
<i>Autostop</i>	<p>In the setting: On, Unitom-5 automatically stops cutting, when the workpiece has been cut through.</p> <p>In the setting: Off, Unitom-5 will only stop when it reaches the pre-set stop position or the STOP Ⓢ key is pressed.</p>
<i>Return position</i>	<p>After cutting or after pressing STOP Ⓢ, the return movement of the cut-off wheel can be set to three different functions:</p> <p><i>Top:</i> Unitom-5 automatically retracts the cut-off wheel to the top position.</p> <p><i>Start:</i> Unitom-5 automatically retracts the cut-off wheel to the original position of the cut-off wheel, at the time you pressed START Ⓞ.</p> <p><i>Stay:</i> The cut-off wheel stays down. Use the Stay function for bakelite bonded diamond or CBN cut-off wheels, as retraction might destroy the rim of the cut-off wheel.</p>
<i>Feed reduction</i>	<p>In the setting: On, Unitom-5 will reduce the pre-set feed speed, if it is too high in relation to the cutting situation.</p> <p>In the setting: Off, Unitom-5 will maintain the pre-set feed speed throughout the cutting process. However, this might result in a motor stopping due to overload.</p>
F1 <i>Default value</i>	The factory setting of every single value can be restored by pressing the function key F1 when the appropriate value is highlighted.

Adjusting the Back Stop

Adjust the back stop of the quick-clamping device to align the centre of the workpiece beneath the centre of the cut-off wheel:

- Loosen the back stop.
- Insert the workpiece and adjust the back stop.
- Tighten the back stop.
- If you also need to adjust the quick-clamping device, place it parallel to the back stop and tighten it.

Clamping Irregular Workpieces

Irregular workpieces without plane clamping surfaces must be clamped using special clamping tools, as the workpieces must not move during the cutting. This could result in damage to the cut-off wheel or to the sample itself. Use the T-slots to mount the special clamping tools. Struers offers a kit of Clamping Tools (See Accessories).

To achieve faster cutting, orientate the workpiece so that the wheel will cut the smallest possible cross-section.

Using the Continue Function

If the cutting is temporarily stopped by pressing STOP  during cutting, the cut-off wheel moves up to the position from which you started the cutting (the cut-off wheel does not, however, retract when Unitom-5 is in STAY function). If you want to continue from the position at which the cutting was stopped, you may press F1: Continue in the CUTTING menu:

CUTTING		Load	POSITION	
Feed	0.50mm/s		Wheel	97.9mm
			Stop	223.0mm
F1:Continue			SET STOP	
			F3:ACT. F4:MAX	

After pressing F1: Continue, the cut-off wheel moves down quickly to the position at which the cutting was stopped, and starts cutting at the pre-set feed speed.

Opening the Protection Guard before using the Continue Function

In case the protection guard has been opened before using the continue function, a message appears, reminding you that the workpiece may have been moved and the cut-off wheel may have been changed.

Using Autostop

When you want to use Autostop only, Autostop is set to On in the Configuration Menu and the stop position is set to MAX. Press **F4**.

When cutting pipes or other workpieces with changing cross sections, the cut-off wheel may retract before the workpiece has been cut. It can then be necessary to switch Autostop off. In the setting: Off, Unitom-5 will only stop when it reaches the pre-set stop position or the STOP  key is pressed.



Using Fixed Stop

*Using Fixed Stop
together with Autostop*

Fixed stop can be used together with the *Autostop On* function. In this situation, the cutting stops at the first occurring stop signal.

If you want to be sure that the cutting is not stopped before a certain position, Autostop should be set to *Autostop Off* in the configuration, using fixed stop only.

Setting the Stop Position

- Clamp the workpiece so that you can move the cut-off wheel along the workpiece.

CUTTING		POSITION
Feed	0.50mm/s	Wheel 100.0mm
	Load	Stop 223.0mm
SET STOP		
F3:ACT. F4:MAX		

- Press POSITION DOWN ▼ to move the cut-off wheel down, along the workpiece, to the desired stop position. Remember to take the wear of the cut-off wheel into consideration.

CUTTING		POSITION
Feed	0.50mm/s	Wheel 200.3mm
	Load	Stop 223.0mm
SET STOP		
F3:ACT. F4:MAX		

- Press F3: ACT. to record the actual stop position. The stop position is now identical to the desired position.

CUTTING		POSITION
Feed	0.50mm/s	Wheel 200.3mm
	Load	Stop 200.3mm
SET STOP		
F3:ACT. F4:MAX		


- Move the cut-off wheel back into top position, and clamp the workpiece correctly.

The Feed Reduction Feature

On Unitom-5, the desired Feed Speed is pre-set and kept constant throughout the cutting process.

If the pre-set Feed Speed is too high in relation to the hardness of the material or size of the sample, Unitom-5 will reduce the Feed Speed automatically (default setting).


The display shows the following text:

Unitom-5 continues to cut at the reduced Feed Speed and reduces the Feed Speed repeatedly, if necessary.

Next time the same type of sample has to be cut, the reduced value for Feed Speed can be used.

After the Feed Reduction Feature has been in operation, the display shows the following, next time START \diamond is pressed:

If the same type of sample has to be cut, press F1.

If Feed Reduction is not desired, feed reduction is set to Off in the configuration. This can be used to force a dressing of the cut-off wheel. However, this might result in a motor stopping due to overload.

Optimising the Cutting Results

The following table shows possible answers to a number of common questions:

Optimising the Cutting Results	
Question	Answer
How can I avoid discoloration or burning of the sample?	Use a lower Feed Speed
	Change the cut-off wheel as the hardness of the present cut-off wheel may be inappropriate for the hardness of the sample*)
How can I avoid burrs?	Use a softer cut-off wheel*)
	Clamp the workpiece securely at the right hand clamping device Tighten the left hand clamping device just enough to prevent the workpiece to shifting on being cut
How can I avoid the cut-off wheels wearing too quickly?	Use a lower Feed Speed or a harder cut-off wheel*)
How can I achieve faster cutting?	Place the workpiece in a direction which allows the cut-off wheel to cut the smallest possible cross-section. Use a high Feed Speed

*) Please refer to the Selection Guide in the [Struers Cut-off Wheels brochure](#).

Connection to External Exhaust System

Struers recommends the use of an exhaust as the materials to be cut may emit harmful gasses or dust.

On the back of Unitom-5 you will find a joint for a 50 mm dia. exhaust hose.

- Remove the red cap.
- Mount an exhaustion hose from your local exhaust system onto the joint.

Mounting the Rubber Curtain

If you need to cut long, protruding workpieces, it is possible to replace one or both of the plastic side plates with a rubber curtain (optional). When you only need one opening in the guard, we recommend that you use the left side to minimize the risk of water spills outside the machine. Please check in advance with your local authorities if this type of protection is approved in your country.

- Remove the 6 nuts holding the plastic side plate and splice plate on the protection guard. Remove the plate.
- Mount the rubber curtain and the splice-plate on the inside of the protection guard. Fasten the rubber curtain with the 6 nuts.
- If necessary, cut a few slots in the curtain to minimise the water spills outside the machine.

Explanation of Safety Factors

The protection guard has a safety switch to prevent the cut-off wheel from starting while the guard is open. Furthermore, a locking mechanism prevents the opening of the protection guard before the cut-off wheel is standing absolutely still.

2. Accessories

Specification	Code
<p>Automatic Cutting Table For Unitom-5. Motorized table for cutting plane parallel sections. Positions the workpiece at a right angle to the cut-off wheel. With 10 mm T-slots and exchangeable steel bands. The table is controlled through the Unitom-5.</p>	UNIEL
<p>Table Unit For Unitom machines. With compartment for Recirculation Cooling Unit.</p>	LABUL
<p>Rubber Curtains To be mounted on Unitom-5 and Unitom-2 for cutting long protruding workpieces. Set of two.</p>	UNIRU
<p>Flange Set 119 mm diameter For use with cut-off wheels 350 mm diameter on Unitom machines.</p>	UNIFL
<p>Cutting Table with 10 mm T-slots Additional table for support of large workpieces. With exchangeable stainless steel bands. 140 x 315 mm (5.5" x 12.4"). For Exotom and Unitom machines.</p>	UNIFT
<p>Replacement Stainless Steel Bands For Automatic Cutting Table (UNIEL) and Cutting Table with 10 mm T-slots (UNIFT). Set of three.</p>	UNIBA
<p>Quick-clamping Device for 10 mm T-slots, Left For securing the workpiece. Complete with back stop. To be mounted on a cutting table placed at the left hand side of the cut-off wheel.</p>	UNILE
<p>Quick-clamping Device for 10 mm T-slots, Right For securing the workpiece. Complete with back stop. To be mounted on a cutting table placed at the right hand side of the cut-off wheel.</p>	UNIRI
<p>Clamping Tools for 10 mm T-slots For clamping complicated workpieces to the cutting table on Unitom-5 and Unitom-2. Complete with clamps, supports and bolts.</p>	UNIOF

3. Consumables

Cut-off Wheel

Please refer to the Selection Guide in the [Struers Cut-off Wheels brochure](#).

IMPORTANT

The max. rotational speed of Unitom-5 is 2775 rpm when running idle. The Struers EXO range of cut-off wheels is marked max. 2720 rpm. However, the rotational speed will always be within an acceptable $\pm 5\%$ tolerance.

Note:

In certain countries, a flange of min. 119 mm (UNIFL) must be used in connection with $\varnothing 350$ mm cut-off wheels.

Other Consumables

Specification	Code
<i>Additive for Cooling Fluid</i> 1 l 5 l	ADDUN ADDFI
<i>Bacterial Cleaning Liquid</i> 0.5 l	BACTE
<i>Disposable Plastic Insert</i> for collecting waste material in the recirculation cooling unit. Package with 10 pcs.	EXOSP

4. Trouble-Shooting

Error	Explanation	Action
Error Messages in the Display		
Motor overload.	The cutting motor has been overloaded for a period and the thermo switch has stopped the motor.	Open the protection guard and wait until the cutting motor has cooled down. The machine can be re-started when the error message disappears in the display.
Position error.	The feed speed has been reduced to the lowest speed (0.05mm/s), by the user or automatically. In spite of this, the feed speed could not be kept up. The cutting process stops.	Choose a cut-off wheel which is more suitable for the material to be cut.
	The cutting arm has been blocked physically.	Check whether the cutting arm is really blocked physically, e.g. by the flange or by the workpiece.
Guard open.	The desired cutting operation cannot start before the guard has been closed.	Close the guard and try again.
No cooling.	The cooling valve is not in cooling position.	Check cooling valve setting.
Display is blank.	The backlight is automatically switched off.	Press any key to reactivate the backlight.
Problems that require a Struers Service Technician		
Unable to lock guard.	The safety lock in the guard is out of order and it is possible to open the guard during the cutting process.	Press the protection guard completely down. If this does not help: Call a Struers service technician.
	The safety switch still stops the cutting motor and the cutting process.	The mechanical safety switch still functions, but the safety lock must be repaired as soon as possible.
Pump overload.	The motor of the recirculation pump has been overloaded. The cutting process stops.	The pump has probably been blocked with cutting debris. Clean the pump, and call a Struers service technician to release the thermo relay.

*Unitom-5
Instruction Manual*

Error	Explanation	Action
Machine Problems		
Chamber light does not work.	Fuse F2 or the bulb may have blown.	Replace fuse F2 or the bulb (NB 230V). The F2 fuse is placed behind the back plate.
Water leaking.	Leak in a hose of Recirculation Cooling Unit.	Check the hose and tighten the hose clamp.
	Water overflow in the cooling water tank.	Remove the excess water in the tank.
Samples or cutting chamber rusty.	Insufficient additive for cooling fluid.	Add Struers Additive for cooling fluid to the cooling water, using the correct concentration. Check with a refractometer. Follow the instructions in the Section Maintenance.
	The machine is left with closed protection guard.	Leave the protection guard open to let the cutting chamber dry.
Protection guard window blurred.	Insufficient cleaning.	Clean with mild soapy water (do not use the cleaning hose).
The protection guard will not open.	The safety lock has stuck.	Open the lock with a thin flat screw driver. Clean and lubricate the lock. Follow the Weekly service instructions in Section Maintenance to avoid the situation.
Quick-clamping device unable to hold the workpiece.	The quick-clamping device is not in balance.	Adjust the screw underneath the clamping column. Use a 3 mm hexagon key.
	Clamping heart worn.	Call Struers service technician.
Cutting Problems		
Discoloration or burning of the sample.	The hardness of the cut-off wheel is inappropriate for the hardness / dimensions of the sample.	See Section Consumables, Cut-off Wheels.
	Inadequate cooling.	-Check that there is enough water in the recirculation cooling unit. -Check the cooling water hoses.
	Feed speed too high.	Reduce the Feed Speed.
Unwanted burrs.	Wheel too hard.	See Section Consumables, Cut-off Wheels.
	Lack of support.	Further support of the workpiece.

*Unitom-5
Instruction Manual*

Error	Explanation	Action
The cutting quality differs.	Cooling water hose clogged.	Clean the cooling water hose and the cooling tube. Check the water flow by turning the cooling valve to cleaning position.
	Insufficient cooling water.	Refill tank with water. Remember Struers Additive.
The cut bends to a side.	The Feed Speed is too high.	Reduce the Feed Speed.
The cut-off wheel breaks.	Incorrect mounting of the cut-off wheel.	- Check that the hole has the correct diameter. - Check cardboard washer on both sides of the cut-off wheel. The nut must be tightened properly.
	Incorrect clamping of the workpiece.	Make sure that only one of the quick-clamping devices is tight. The other device should only press lightly. Use support tools if the geometry of the workpiece makes support necessary.
	Wheel too hard.	See Section Consumables, Cut-off Wheels.
	Feed speed too high.	Decrease the Feed Speed.
	Inadequate cooling.	- Check that there is enough water in the recirculation cooling unit. - Check the cooling water hoses.
The cut-off wheel wears down too quickly.	Feed speed too high.	Reduce the Feed Speed.
	Insufficient cooling.	- Check that there is enough water in the recirculation cooling unit. - Check the cooling water hoses.
	The cut-off wheel is too soft for the task.	See Section Consumables, Cut-off Wheels.
	Unitom-5 vibrates (worn bearings).	Call Struers service technician.

*Unitom-5
Instruction Manual*

Error	Explanation	Action
The cut-off wheel does not cut through the sample.	Incorrect choice of cut-off wheel.	See Section Consumables, Cut-off Wheels.
	Cut-off wheel worn.	Replace the cut-off wheel.
	Diameter of cut-off wheel too small.	Change to a cut-off wheel with another diameter, e.g. ø350 mm.
	The cut-off wheel gets caught in the workpiece.	Support the workpiece and clamp it on both sides of the cut-off wheel in such a way to allow the cut to stay open.
	The cutting is stopped by the Autostop function.	Set Autostop to Off in the configuration.
	The stop position is set too low.	See Section Basic Operation, Stopping the cutting.
The workpiece breaks when clamped.	The workpiece is brittle.	Place the workpiece between two polystyrene plates. NB! Always cut brittle workpieces very carefully.
The sample is corroded.	The sample is not resistant to water.	Use a neutral liquid as cooling fluid or cut without using cooling fluid at all. DO NOT USE AN INFLAMMABLE LIQUID (You may lift the pump carefully out of the recirculation cooling unit).
	The sample has been left in the cutting chamber for too long.	Leave the protection guard open, when you leave the machine.
	Insufficient additive for cooling fluid.	Add Struers Additive for cooling fluid to the cooling water in the correct concentration. Check with a refractometer. See Section Maintenance.

5. Maintenance

Service and greasing of the cutting table and quick-clamping devices is part of a regular Struers Service Schedule. Please contact your Struers Representative for details.

Daily Service

- Clean the cutting chamber, especially the cutting table with the T-slots.
- Clean the protection guard window with a damp cloth. Do not use the cleaning hose.

IMPORTANT

Do **not** clean the lamp glass with alcohol.
Use a damp cloth.

Checking the Recirculation Unit

- Empty and clean the sieve in the Recirculation Cooling Unit.

The cooling unit should be checked for cooling water after 8 hours use or at least every week. Refill the cooling unit if the flushing pump cannot reach the cooling water.

Remember to add Struers Additive for Cooling Fluid: One part of Additive for 33 parts of water.

Use a refractometer to check the concentration of additive.

Concentration = Brix value.

The concentration of additive should lie between 2.7 and 3.3 %.

Add Struers Additive for Cooling Fluid if the concentration is too low.

IMPORTANT

Always maintain the correct concentration of Struers Additive in the cooling water (percentage stated on the container of the Additive). Remember to add Struers Additive for cooling fluid each time you refill with water.

Weekly Service

- Clean the cutting chamber thoroughly.
- Remove the pump and clean the sieve.
- Clean and lubricate the safety lock.

Monthly Service

Replacing the Cooling Water

- Replace the cooling water in the Recirculation Cooling Unit at least once a month.
- Place a disposable plastic insert in the tank and fold it over the edge.

Please Note

Too high a level of coolant in the tank might damage the pump. To avoid this: place the disposable container so that the pre-punched hole is in front of the overflow aperture in the tank.

- Fill the tank with 63 l of water and 2 l Struers Additive for cooling fluid. The water level should be 8-10 cm below the upper edge of the tank.

IMPORTANT

We strongly warn against refilling without removing the cooling tank from the cabinet. This could lead to overflow in the cooling tank, as you cannot control the water level when the cooling tank is inside the cabinet.

- Put the lid containing the sieve on top of the tank.
- Place the pump in the lid (fig. 1).
- Wheel the trolley into the table. Place the return hose in the sieve (fig. 2).


Fig. 1



Fig. 2

IMPORTANT

Always maintain the correct concentration of Struers Additive in the cooling water (percentage stated on the container of the Additive). Remember to add Struers Additive for cooling fluid each time you refill with water.

*Lubrication of Bearings for
Tilting Frame*

- Press emergency stop.
- Open the protection guard.
- Remove the red protection caps from each of the two grease nipples, located as indicated with a white circle in the photo.
- Lubricate in the two grease nipples with the grease gun and move the spindle/cutting motor up and down by hand, until grease comes out at the shaft. Use water-repelling grease, e.g. Shell Albida Grease LX.
- Clean superfluous grease away and mount the protection caps at the grease nipples.
- Release the emergency stop.

Maintenance of Cutting Tables

The iron rails (available as spare parts) should be replaced if they have become worn or damaged.

Maintenance of Cut-off Wheels

Storing of Bakelite Bonded Al_2O_3 Cut-off Wheels

These cut-off wheels are sensitive to humidity. Therefore, do not mix new, dry cut-off wheels with used humid ones. Store the cut-off wheels in a dry place, horizontally on a plane support.

Maintenance of Diamond and CBN Cut-off Wheels

The precision of diamond and CBN cut-off wheels (and thus the cut) depends on how carefully the following instructions are observed:

- Never expose the cut-off wheel to overload, such as heavy mechanical load, or heat.
- Store the cut-off wheel in a dry place, horizontally on a plane support, preferably under light pressure.
- A clean and dry cut-off wheel does not corrode. Therefore, clean and dry the cut-off wheel before storing. If possible, use ordinary detergents for the cleaning.
- Regular dressing of the cut-off wheel is also part of the general maintenance.

Maintenance of Clamping Devices

Important

It is recommended to thoroughly clean and lubricate the Quick Clamping Device and Vertical Clamping Device at regular intervals.


6. Technical Data

Subject		Specifications	
		Metric/International	US
Cutting	Three-phased cutting motor Rotational speed (running idle) Feed Speed range (adjustable in steps of:) Max. positioning speed Cut-off wheel position (Adjustable in steps of) Cutting force	4.7kW 2775rpm 0.05 - 3mm/s (0.05mm/s) 30mm/s 0-223mm (1mm) Max. 350N	6.4HP 2775rpm 0.002-0.118"/s (0.002"/s) 1.2"/s 0-8.75" (0.04") Max. 350N
Cut-off Wheel	Diameter x thickness	ø350x2.5x32mm ø300x2x32mm	13.8x0.1x1.25" 11.8x0.08x1.25"
Cutting Chamber Geometry	<i>Max. workpiece dimensions:</i> Thickness Width Length <i>Max cutting capacity:</i> Max. cutting thickness Max. length of cut (with a sample thickness of:) (The max. cutting capacity might be less, depending on the sample material and the cut-off wheel).	150mm 415mm 630mm 120mm 200mm (90mm)	5.9" 16.3" 24.8" 4.7" 7.9" (3.5")
Cutting Table	An optional table may be mounted (EXOBA)	150x300mm	5.9x11.8"
Recirculation Cooling Unit	Approx. flow: Tank volume	12 l/min 65 l	3.2 gallons/min 17 gallons
Dimensions and Weight	<i>Unitom-5</i> Width Depth Height Weight <i>Optional Table Unit (LABUL)</i> Width: Depth: Height:	700mm 750mm 660mm 160kg 900mm 750mm 800mm	27.5" 29.5" 26" 353lbs 35.4" 29.5" 31.5"

Unitom-5
Instruction Manual

Subject		Specifications
Software and Electronics	LC display Controls Memory	64 x 128 dots (8 x 40 characters) <i>Touch Pad</i> EEPROM/RAM
Environments	Safety standard Dynamic Balancing Noise level	IEC 204-1/EN 60204-1(VDE 0113), 89/392/EEC, 89/339/EEC Max. admissible unbalance U_{per} according to ISO 1940/1, Balance Quality Grade G6.03: $U_{per} = 2$ gmm/kg = 110 gmm. Approx. 76dB(A) during cutting, at a distance of 1.0 m/39.4" from the machine.
Electrical Data	<i>Voltage/frequency:</i> 3 x 200V/50Hz 3 x 200-210V/60Hz 3 x 220-230V/50Hz 3 x 220-240V/60Hz 3 x 380-415V/50Hz 3 x 460-480V/60Hz	<i>Max. Load:</i> 21.3A 21.3A 19.8A 19.8A 11.8A 11.8A

7. Menu Overview

Quick Reference Card

Clamping the Workpiece

- Place the workpiece between the clamp and the back stop at the right hand cutting table.
- Push the clamp towards the workpiece and lock the quick-clamping device with the locking handle.

Starting the Cutting

- Lower the cut-off wheel into correct position above the workpiece by pressing POSITION DOWN ▼. The cut-off wheel should be min. 1 mm above the workpiece.
- Carefully close the protecting guard.
- Press START ◇. The cut-off wheel starts rotating and the cooling water starts running.

Stopping the Cutting

Automatic Stop

- Unitom-5 automatically stops cutting when the sample has been cut through (default setting). The cut-off wheel moves to start position (default setting). The cut-off wheel stops rotating and the cooling water stops.

Manual Stop


- Press STOP ⊕ and the cutting process stops. The cut-off wheel moves to start position (default setting). The cut-off wheel stops rotating and the cooling water stops.

Change Feed Speed

The feed speed can be set before the cutting and changed during the cutting.

- Press ENTER ↵ to edit the value.
- Press MENU UP/DOWN ⏮▲▼ to change the value.
- Press ENTER ↵ to accept the new value.

Changing the Cut-off Wheel

- Open the protection guard and swing the right cooling tube to the right.
- Press the black locking knob and turn the cut-off wheel until the spindle lock clicks.
- Remove the nut with a spanner. Remove the flange and the old cut-off wheel.
- Mount the new cut-off wheel. Mount the flange and the nut. Tighten carefully.
- Swing the right cooling tube back into position.

Flushing the Cutting Chamber

- Point the cleaning hose towards the bottom of the cutting chamber.
- Turn the cooling valve clockwise to flushing position.
- Clean the cutting chamber.
- Turn the cooling valve counter-clockwise to cooling position.
- Always leave the protection guard open to let the cutting chamber dry.

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www.struers.com

Unitom-2/-5/-50
Spare Parts and Diagrams

**Always state *Serial No* and *Voltage/frequency*
if you have technical questions or when ordering spare parts.**

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Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to change without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2001 .

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Telefax: +45 44 600 801

Spare Parts and Diagrams

Table of contents

Drawing

Unitom-2

Drawings

Unitom-2, complete.....	15370001U
Case bed	15370010J
Spindle unit, complete	15370015E
High power plate, rear side	15370030E

Diagrams

Main Circuit Diagram	15373100D
Power Circuit Diagram.....	15373110C
Control In/Out Circuit Diagram	15373120C
Circuit Diagram, Logic control circuit	15373200D
Wiring of Power connection diagram.....	15373420E
Wiring diagram, Back-Panel, Control	15373410E

Unitom-5

Drawings

Unitom-5, complete.....	15330002G
Case Bed, assembly	15330010D
Cutting Arm, assembly.....	15330020D
Spindle unit, complete	15330015G
Stepper motor/gear unit, complete	15330060F
High power plate, rear side.....	15330030D
UNIEL, Automatic Cutting table complete.....	15330080F
Automatic Cutting table.....	15330081A
Stepper Motor with pulley, complete.....	15330082
Threaded Spindle, assembly	15330083

Unitom-50

Drawings

Unitom-50, complete.....	15310001J
Case Bed, assembly	15310010G
Cutting Arm, assembly.....	15310020H
Spindle unit, complete	15330015G
Stepper motor/gear unit, complete	15330060F
High power plate, rear side.....	15310050C
Cutting table w. excicut mechanism, complete.....	15310030B
Excicut drive mechanism	15310060E
Excicut cutting table, cross section.....	15310070

Common Parts to Unitom-2/-5 and Unitom-50

Drawings

Cutting Table, complete	15330003B
Locking system, complete.....	15330005C
Protection guard, complete.....	15330042F
Lamp, complete.....	15330044B
Water Valve, complete.....	15330025D
Recirculation tank	14980068A
Quick clamping device, left hand side	15330006B
Quick clamping device, right hand side.....	15330007B
Outlet for Unitom-5/-50	15317049

Diagrams for Unitom-5 and Unitom-50*

Block/Wiring diagram, Unitom-5	15333100E
Block/Wiring diagram, Unitom-50	15313100C
Power circuit diagram, Unitom-5	15333110A
Power circuit diagram, Unitom-5/-50	15313110A
Control circuit diagram, Unitom-5	15333120C
Control circuit diagram, Unitom-50	15313120C
Circuit diagram, Unitom-5/-50 (2 pages).....	15313200C
Wiring diagram, back panel, control, Unitom-5	15333450D
Wiring diagram, back panel, control, Unitom-5/-50	15313455D
Wiring diagram, power connection, Unitom-5	15333460B
Wiring diagram, power connection, Unitom-5/-50	15313465A
Wiring diagram, front panel, control PCB	15313475B
Wiring diagram, front panel, control PCB, Unitom-5/-50	15333470D
X-Table Mover, UNIEL	15333480

Some of the drawings may contain position numbers
not used in connection with this manual.

*Unitom-5's chassis and back panel will gradually be replaced by the
Unitom-50 chassis and back panel.

Unitom-2/-5/-50
Spare Parts and Diagrams

Spare Parts list, Unitom-2

Drawing	Pos.	Spare Part	Cat no:
15370001		Unitom-2, complete	
	170	50Hz timing belt 600-8MP-30mm	2JT20600
	170	60Hz timing belt 560-8MP-30mm	2JT20560
	190	Gas spring 700N Lift_O-Mat S100 with bellows	RYS00190
	310	1.5 AT Fuse Mini TR5	2FU32950
	310	Unitom-2 PCB, complete	15373001
	382	Gas Spring 400 N CAMLOC SL=150 (2 pcs.)	RYS00370
	660	Fixed flange with key	R4981210
	661	Cardboard blotter ø32/105 f. xxUNI (25 pcs)	15332902
	665	Key A 4x4x10	2ZF10410
	670	Loose flange	R4981209
	680	Nut M20	2TA10200
	690	Disc ø21/37/3 A2 DIN 125	2ZA10020
	770	Transparent silicon sealant for window glass	2LL00330
	800	Cooling pump (3x220-240V/50Hz)	381MP408A
	800	Cooling pump (3x380-440V/50Hz)	381MP408B
	800	Cooling pump (3x200-220V/60Hz)	381MP409A
	800	Cooling pump (3x380-480V/60Hz)	381MP409B
	801	Sieve for 2AB4 (ø99)	14980074
	15370010		Case Bed
20		Machine supports	14980015
40		Gasket for flange	15330150
45		Flange for exhaust (not shown)	13811000
200		Water Valve, complete without knob	15330025
210		Knob for Cooling Water valve with grease	15332904
220		Safety Switch, T3C 236-O2Z	15330120
260		Sealing ring JS16	2IF03713
290		Flushing hose, complete	15330075
15370015		Spindle Unit, Complete	
	20	Ball bearing, 6007 2RS1 ø35/ø68	2BK00110
	50	Sealing Ring w. grease	14982901
	120	Stop pin	14981231
	130	Spring for stop pin	14981232
	140	Scraper ring	2IT00006
	150	Ball handle	2GH00112
15370030		High Power Plate, Rear Side	
	90	0.1AT fuse glass tube	2FU13300
	100	2.50 AT fuse glass tube (10 pcs.)	RFU14000

Unitom-2/-5/-50
Spare Parts and Diagrams

Spare Parts list, Unitom-5

Drawing	Pos.	Spare Part	Cat no:	
15330002		Unitom-5, complete		
	170	50Hz timing belt 600-8MP-30mm	2JT20600	
	170	60Hz timing belt 560-8MP-30mm	2JT20560	
	190	Gas spring 600N Lift-O-Mat S100 with bellows	RYS00182	
	375	Bellows,V6-26G	RIB00007	
	410	Gas Spring 400 N CAMLOC SL=150 (2 pcs)	RYS00370	
	600	Back Jaws with screws for 10 mm T-slots	R5330310	
	660	Fixed flange with key	R4981210	
	661	Cardboard blotter ø32/105 f. xxUNI (25 pcs)	15332902	
	665	Key A 4x4x10	2ZF10410	
	670	Loose flange	14981209	
	680	Nut M20 A2	2TA10200	
	690	Disc ø21/37/3 A2 DIN 125	2ZA10020	
	770	Transparent silicon sealant for window glass	2LL00330	
	800	Cooling pump (3x220-240V/50Hz)	381MP408A	
	800	Cooling pump (3x380-440V/50Hz)	381MP408B	
	800	Cooling pump (3x200-220V/60Hz)	381MP409A	
	800	Cooling pump (3x380-480V/60Hz)	381MP409B	
	801	Sieve for 2AB4 (ø99)	14980074	
	15330010		Case bed, assembly	
		20	Machine supports	14980015
40		Gasket for flange	15330150	
45		Flange for exhaust	13811000	
200		Water Valve, complete without knob	15330025	
210		Knob for Cooling Water valve with grease	15332904	
220		Safety Switch, T3C 236-O2Z	15330120	
290		Flushing hose, complete	15330075	
15330020		Cutting arm, assembly		
	20	Bearing housing, complete	14990061	
	40	Oil Seal Ring SD 25x35x4	2IS30025	
	50	Grease nipple H1 M6x1 Straight Male	2GN20004	
	51	Grease Gun, Shell Albida LX female	R5332901	
	55	Protection cap with strap	2GN20004	
	135	Timing Belt PU6 T2.5 L=245	2JT50245	
	180	Emergency stop	260MP007	

Unitom-2/-5/-50
Spare Parts and Diagrams

Spare Parts list, Unitom-5

Drawing	Pos.	Spare Part	Cat no:
	200	Switch element	381MP440
	205	Protective seal for PCB	15331123
	210	PCB Unitom-5, w. Display & Prom	15333001
	211	PROM, UNITOM-5, latest version	R5333901
	213	1.60 AT Mini Fuse TR5	2FU32800
	214	3.15 AT Mini Fuse TR5	2FU32950
	270	Locking system, complete	15330005
15330015		Spindle unit, complete	
	120	Stop pin	14981231
	130	Spring for stop pin	14981232
	140	Scraper ring	2IT00006
	150	Ball handle	2GH00112
15330060		Stepper motor/gear unit, complete	
	20	Resolver, complete	15332903
	120	Stepper motor with toothed wheel	15330060
	130	Gear for stepper motor	15332000
15330030		High power plate, rear side	
	20	Transformer 200-480V/50-60Hz	2MT70711
	40	Contactora CI-12 for 24V/50-60Hz	2KM10232
	50	Contactora K2-09A01 24V/50-60Hz	2KM09024
	60	Thermal relay U12/16em 0.6-0.9 A	2FS10024
	60	Thermal relay U12/16em 0.4-0.6 A	2FS10023
	60	Thermal relay U12/16em 0.18-0.27 A	2FS10021
	90	0.25 AT Fuse, Glass Tube	2FU13500
	100	4.00 AT Fuse, Glass Tube	2FU14200
15330080		UNIEL, Automatic Cutting table complete	
	10	Threaded Spindle, complete	15330083
	30	Bellows	15330517
	40	Flange for Bellows	15330516
	70	Inductive sensor, DJ 2G	2HQ00023
	110	Sealing band (1 m)	2IP10202
	170	Automatic cutting table	15330081
	180	Guide rail	15330505
	190	O-ring, ø20.2x 3	2IO30007
	220	Timing Belt, HTD 180-3M-9 ZG 60	2JT20180
	230	Stepper Motor with pulley, complete	15330082
	235	Vibration damper	15330521

Unitom-2/-5/-50
Spare Parts and Diagrams

Spare Parts list, Unitom-5

Drawing	Pos.	Spare Part	Cat no:
	310	Gasket	15330511
	350	Shell Helix 10-40W oil (1L)	2LS11040
	360	DBI Dut no. 8	2GT20008
15330081		Automatic Cutting table	
	20	Radial Bearing $\varnothing 25/35 \times 25$	2BG00020
	30	Sealing ring SD 25x35x4	2IS30025
	40	Scraper ring AS 25-35-7-10	2IT00025
	50	Spindle Bushing	15330510
	60	Lock Ring, J36 DIN 472 (10 pcs)	RZL20360
15330082		Stepper Motor with pulley, complete	
	20	Pulley for Stepper motor	15330514
	30	Lock pin, $\varnothing 2 \times 12$ DIN 1481	2ZS04242
15330083		Threaded Spindle, assembly	
	10	Threaded spindle	15330083
	20	Spherical ball bearing 2201 E-2RS1	2BK20012
	30	Key 4x4x10	2ZF10410
	40	Pulley for spindle	15330513

Unitom-2/-5/-50
Spare Parts and Diagrams

Spare Parts list, Unitom-50

Drawing	Pos.	Spare Part	Cat no:	
15310001		Unitom-50, complete		
	170	50Hz timing belt 600-8MP-30mm	2JT20600	
	170	60Hz timing belt 560-8MP-30mm	2JT20560	
	190	Gas spring 600N Lift-O-Mat S100	RYS00182	
	410	Gas Spring 400 N CAMLOC SL=150	RYS00370	
	660	Fixed flange with key	R4981210	
	661	Cardboard blotter ø32/105 f. xxUNI (25 pcs)	15332902	
	665	Key A 4x4x10	2ZF10410	
	670	Loose flange with nut M20 and disc	R4981209	
	680	Nut M20 A2	2TA10200	
	690	Disc ø21/37/3 A2 DIN 125	2ZA10020	
	770	Transparent silicon sealant for window glass	2LL00330	
	800	Cooling pump (3x220-240V/50Hz)	381MP408A	
	800	Cooling pump (3x380-440V/50Hz)	381MP408B	
	800	Cooling pump (3x200-220V/60Hz)	381MP409A	
	800	Cooling pump (3x380-480V/60Hz)	381MP409B	
	801	Sieve for 2AB4 (ø99)	14980074	
	805	Reset Push button	2SA00301	
	15310010		Case bed, assembly	
		20	Machine supports	14980015
40		Gasket for flange	15330150	
45		Flange for exhaust	13811000	
200		Water Valve, complete without knob	15330025	
210		Knob for Cooling Water valve with grease	15332904	
220		Safety Switch, T3C 236-O2Z	15330120	
290		Flushing hose, complete	15330075	
15310020		Cutting arm, assembly		
	20	Bearing housing, complete	14990061	
	40	Oil Seal Ring SD 25x35x4	2IS30025	
	50	Grease nipple H1 M6x1 Straight Male	2GN20004	
	51	Grease Gun, Shell Albida LX female	R5332901	
	55	Protection cap with strap	2GN20004	
	90	Stepper motor/gear unit complete	15330060	
	135	Timing Belt PU6 T2.5 L=245	2JT50245	
	180	Emergency stop	260MP007	

Unitom-2/-5/-50
Spare Parts and Diagrams

Spare Parts list, Unitom-50

Drawing	Pos.	Spare Part	Cat no:
	200	Switch element	381MP440
	205	Protective seal for PCB	15331123
	210	PCB Unitom-5/-50 w. display	15313001
	211	PROM, UNITOM-5, latest version	R5333901
	213	1.60 AT Sub Mini Fuse TR5	2FU32800
	214	3.15 AT Sub Mini Fuse TR5	2FU32950
	270	Locking system, complete	15330005
15330015		Spindle unit, complete	
	120	Stop pin	14981231
	130	Spring for stop pin	14981232
	140	Scraper ring	2IT00006
	150	Ball handle	2GH00112
15330060		Stepper motor/gear unit, complete	
	20	Resolver, complete	15332903
	120	Stepper motor with gear wheel	15330060
	130	Gear for stepper motor	15332000
15310050		High power plate, rear side	
	20	Transformer 200-480V/50-60Hz	2MT70711
	110	Contactora K2-09A01 24V/50-60Hz	2KM09024
	120	Thermal relay U12/16em 0.6-0.9 A	2FS10024
	120	Thermal relay U12/16em 0.4-0.6 A	2FS10023
	120	Thermal relay U12/16em 0.18-0.27 A	2FS10021
	130-	Contactora CI-12 for 24V/50-60Hz	2KM10232
	140		
	150	Mechanical blocking, Danfoss	2KH10100
	160	Contactora CI-12 for 24V/50-60Hz	2KM10232
	190	4.00 AT Fuse, Glass Tube (10 pcs)	RFU14200
	200	0.25 AT Fuse, Glass Tube	2FU13500
15310030		Cutting table with ExciCut mechanism, complete	
	10	ExciCut drive mechanism, complete	15312902
	40	Bellows V6-307	2IB00307
	50	Sparks Guard	15310260
	130	Spark cover for bellows	15310130
	140	Bearing house, complete	15310267
	150	Grease nipple	2GN20004
	151	Grease Gun, Shell Albida LX female	R5332901
	170	Cutting table	15310070
	235	T-rail steel bands	05336905

Unitom-2/-5/-50
Spare Parts and Diagrams

Spare Parts list, Unitom-50

Drawing	Pos.	Spare Part	Cat no:
15310060		Excicut drive mechanism	
	10	Planet gear PLE 40-M2	15312903
	40-120	Crankshaft with connection rod, complete	15312901
	130-160	Excicut Motor, pre-adjusted	15312902
15310070		Excicut Cutting table, cross section	
	10-35	Excicut Cutting table w. bearings	15312904
	20	Grease nipple	2GN20004
	21	Grease Gun, Shell Albida LX female	R5332901
	40	Shaft for excicut table	15310161
	50	Sleeves $\varnothing 20/\varnothing 30 \times 15$	15310157

Unitom-2/-5/-50
Spare Parts and Diagrams

Common Parts to Unitom-2/-5/-50




Drawing	Pos.	Spare Part	Cat no:
15330003		Cutting table, complete	
	10	Cutting table, complete	03816901
	20	T-rail plate, complete	15330509
15330005		Locking System	
	80	Actuator, MK1 24VDC	2MK70524
15330042		Protection Guard without panes	
	20	Rubber side curtain I+II	05336902
	20	Side Guard w. rubber sealings	15330045
	90	Rubber profile, 1038-02 EPDM SHORE A 1m)	2IP21038
	100	Sealing profile 209-0203	RIP00203
	110	Rubber curtain, front	14981450
14980510		Lamp, complete	
	20	O-ring	2IO18065
	90	Bulb E27	2HG11741
	110	End cap	15330044
15330025		Water Valve, complete	
	20	Inductive Sensor, DJ 2G	2HQ00023
	30	O-ring 25.00-2.00 72 NBR 872	2IO20055
	90	PVC Hose ø12/ø18, armed (1m)	2NU21218
14980068		Recirculation tank	
	90	Cooling Tank, 65 l	381MP140
	100	Lid	381MP142
	110	Sieve for tank, complete	315MP052
	115	Red Plastic Cover for tank (ø117)	2GT20117
	117	Drain hose ø38	2NU30238
15330006/7		Quick clamping device (right and left)	
	60	Clamping Rings for left & right (2 pcs)	315MP005
	70	Clamping Spring for left & right (2 pcs)	12150155
	100	Grease Gun (60 ml), Shell Albida LX male	381MP415
	130	Screw head protection knobs	2GT10005
	140	Jaw for quick clamping device, left	15330311
	141	Jaw for quick clamping device, right	15330312
	170	T-Bolt for 10 mm rails	RTF40810

Unitom-2/-5/-50
Spare Parts and Diagrams


Common Parts to Unitom-2/-5/-50

Drawing

Pos.	Spare Part	Cat no:
	Various	
	Fork Spanner, 30 mm	2GR00230
	Angle Pipe wrench, 13 mm	2GR00313
	Pipe Wrench 17 mm	2GR00007
	Allen Key 6	2GR00060
	Grease Gun, Shell Albida LX female	R5332901
	Tray (for Unitom-50 only)	15310301

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



The visible part of the axle, must be coated with tectyl.

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






The back side of the sealing must be cut


E	22.11.2007	Pos.50 2KM09024->2KM10640	SPE	22.11.2007	CJE
D	07.06.2007	Translated to english. Pos. 220+225+230->Pos.250 Coolistik	BMJ	07.06.2007	CJE
C	03.04.1998	F1->F2;F2->F3. F2 rykket så diag.+tegn. stemmer overens	BMJ		
B					
A	10.06.1997		BMJ	02.07.1997	JTV
Revision	Crea. date dd-mm-yy	Revision description	Draw. Init	Appr. date dd-mm-yy	Appr. Init
		Material:	Scale: 1:2	Format: A3	Tolerance: DS/ISO 2768- Surface treat.: mK None
		ID: Description:	15370030 Electrical back plate, assembled		Rev: E

Pedersstrupvej 84
 DK-2750 Ballerup/Copenhagen
 Denmark
 Phone: +45 44 600 800
 Fax: +45 44 600 904

M1
CUTTING MOTOR

M2
RECIRCULATION MOTOR

CUTTING MOTOR M1,
NOT IN OVERLOAD.

RECIRCULATING PUMP
NOT IN OVERLOAD

RECIRCULATING PUMP ON.

CUTTING MOTOR RUN

SAFETY GUARD LOCK

K1
CONTACTOR

K2
CONTACTOR

YM1
AKTUATOR

START
PROCES ON

STOP
PROCES OFF

WATER VALVE IN
COOLING POSITION.


WATER VALVE IN
SLICING POSITION.

EMERGENCY STOP
AND
OVERLOAD AT

RECIRCULATION -
OR
CUTTING MOTOR

SAFETY GUARD
CLOSED.

SAFETY GUARD
RELEASE

15373110 POWER CIRCUIT DIAGRAM


15373200 PCB CIRCUIT DIAGRAM

15373410 WIRING DIAGRAM, CONTROL, BACK-PANEL


15373420 WIRING DIAGRAM, POWER, BACK-PANEL

15373430 WIRING DIAGRAM, CONTROL, FRONT-PANEL

Rev. C: Test change 1 YW1: 10V → 24V F1 → F2, G.R.40.3 F2 → F3, G.R.40.2 B1 → B3, G.R.40.2 S2 → S4, G.R.40.2/40.2 G.R.40.40.40.40.40: B Rev. D: (P111 13-12-04) Ref. G1 changed to QS1 Colors on BQ14802 changed G.R. = Gnd Referenced Friday, December 17, 2004	Rev. A: COLOR CHANGE: YW1: 10V → 24V YW1: RD → BK, W1 → BN, BN → VT, YW1: BN → BU, BU → RD, X5.1 → -5, -5 → -1, X5.4 → BU → X5.5 X5.5 RD → X5.4	STROEBERG AG MAINFABRIK ALLE 178 D-62303 FRIEDRICHSDORF FRIEDRICHSDORF PHONE: +49 69 3900 3900 FAX: +49 69 3900 3900 E-MAIL: info@stroebert.de WWW: www.stroebert.de MAIN CIRCUIT DIAGRAM: UNITOM - 2: Size A1 SCALE 1:2 mg 97 / Dan Molding DWG NO 15373100 Sheet 1 of 2
--	---	---

Rev.: A. Text change	STRUERS A/S VALHOESALLE 176 DK-2610 ROEDOVRE DENMARK PHONE: +45 3670 3500		
Rev.: B.: Text change ! F1 -> F2, G.R.->J,4 F2 -> F3, G.R.->J,3	UNITOM -2: POWER CIRCUIT DIAGRAM		
Rev.: C. Ref. Q1 changed to QS1	Size A2	CAGE Code	DWG NO 15373110
Monday, December 13, 2004	Scale	D 16.06.97 / DEM	Sheet 1 of 1

Rev.: B.: Text change 1
S1 -> S3, G.R.=E.2
S2 -> S4, G.R.=F.2/G.2
HL** -> Wiped off.

Rev.: C.: Colors on
BQ1 & BQ2 changed

Rev.: A.: COLOR
CHANGE:
YM1: 12V -> 24V,
BN -> VT,
YM1: GN -> BU; BU -> RD,
X5:1 -> -5, -5 -> -1.


Friday, December 17, 2004

STRUJERS A/S
VALHOESALLE 176
DK-2810 ROEDDovre
DENMARK
PHONE: +45 3670 3500


Control In / Out Circuit Diagram UNITOM - 2:

Size	A2	CAGE Code	DWG NO	Rev
			15373120	I C
Scale		D.16.06.97 / DEM	Sheet	1 of 1


X1: Power Supply


Chassis


Signal In:


Signal In:


Signal In:

Signal In:

Signal In:

X2:

24V AC, TO CONTROL CIRCUIT I.
 IN I 24V AC, FROM CONTROL CIRCUIT.
 CUTTING MOTOR ON.
 IN II 24V AC, FROM CONTROL CIRCUIT.
 RECIRCULATING PUMP ON.
 24V



X3:

ALARM FOR OVERLOAD OR EMERGENCY STOP
 SLICING ON
 +24V DC OUT.
 "SAFETY GUARD" RELEASE
 "SAFETY GUARD" LOCK
 GND-P









X4:

START PROCES ON
 STOP PROCES OFF
 SAFETY GUARD CLOSED
 ACTUATOR: "SAFETY GUARD" IN POS. LOCK
 OVERLOAD, EMERG. CUT. MOTOR REC. MOTOR EMERGENCY STOP
 COOLING VALVE READY
 SLICING VALVE READY

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


For connection of Unitom-5/-50 to supporting table, use drain elbow pipe*(55mm diameter). A 65mm drain hose throughhole in the table is recommended. Use hose clamp* for mounting drain hose* onto drain elbow pipe*.

*) supplied with Unitom-5/-50

Matr.:	Overfl.beh.:	Målförhold:	Ikke ang. tol. efter DS/ISO 2768-	
	Projektionsmetode	1:5	Date	Sign.
			Tegn:	200300 BMJ
			Kontr.:	220300 HAB
Outlet for Unitom-5/-50		Erst.:	15317049	

M1 CUTTING MOTOR

M2 RECIRCULATION MOTOR

CUTTING MOTOR M1, NOT IN OVERLOAD.

RECIRCULATING PUMP NOT IN OVERLOAD

RECIRCULATING PUMP ON.

CUTTING MOTOR RUN

CUTTING MOTOR DELTA.

CUTTING MOTOR STAR

SAFETY GUARD LOCK

STEPMOTOR, MOVING TABLE

STEPMOTOR, ARM DRIVE.

GRAPHICS LCD DISPLAY MODULE (240 X 64 DOTS)

TASTATUR 6 IN. 2 OUT.

CUTTING ARM POSITION.

WATER VALVE IN COOLING POSITION.

WATER VALVE IN SLICING POSITION.

OVERLOAD FROM M1 THERMOSWITCH

RECIRCULATION IN OVERLOAD.

SAFETY GUARD OPEN.

EMERGENCY STOP ON


CUTTING MOTOR M1, CURRENT ALT.

SAFETY GUARD RELEASE

X-TABEL REFERENCE POSITION

SEE DRAWING NO. 1533200

Rev. 4: MATHON, J. 28 Rev. 5: COLOR CHANGE New: Z1 FILTERS YMH: ZV = 5V YMH: RD = BK, WH = BN, BN = VT. YMH: GN = BU, BU = RD. Rev. 3: Test removed Rev. 2: Ref. Q1 changed to Q51, BQ1 & BQ2 wires changed.	STROBERE Z1E RECIRCULATING PUMP M4 CNC-21753 MALLERUP DINNAMAK TROMSØ - 44 5008 BQ1 BLOCKWIRUNG DIAGRAM Mark 10 - 0404 - 0206-1 Rev. 03: Test removed Rev. 02: Ref. Q1 changed to Q51, BQ1 & BQ2 wires changed.	UNITOM - 5:
Size A1	CAGE CODE	DWG NO 15333100
SCALE	21 March 97 (Date Making)	Sheet 1 of 1
December 16, 2004		

TASTATUR
6 IN. 2 OUT.

CUTTING ARM
POSITION.

WATER VALVE IN
COOLING
POSITION.

WATER VALVE IN
FLUSHING POSITION.

OVERLOAD
FROM M1
THERMOSWITCH

WIRE STRAP:
CODE SIGN FOR
"EXCICUT TABLE CONNECTION"

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
STEPMOTOR, ARM DRIVE.

GRAPHICS LCD
DISPLAY MODULE
(240 X 64 DOTS)


- Rev.C: (03-01-2025 CH)
- removed YM1, X5
 - removed wires BU13, BU14, BU15, BK, BN, VT
 - changed part S2 safety
 - removed wires BU11 + RD, BU12 + BU
 - new wires: X43 to S2A2 <BU12>, X41 to S2A1 <BU11>
 - removed X23 and wire RD5
 - new wire: X62 to S221 <RD105>
 - removed X12, X3.1, and wires RD1, RD2, BU2
 - new wire: X83 to S211 <BU102>
 - removed X11, RD1, BU1 changed to BU101
 - new connection: S2:12 to X3:2
 - removed RD6
 - new connection: S2:22 to X2:4 to M1 to X8:1 <RD106>
 - changed label FT2 to FT1

Rev. A: Test removed		STRUENEGE AS	
Rev. B: (FF11-13-2024)		PROJECT: UNITOM V1E1_04	
Rev. C: changed in 03/01/2025		DINNAMIK	
Rev. D: changed in 03/01/2025		S21MS22 color changed	
BLOCK: WIRING DIAGRAM		UNITOM - 50:	
With Excicut Table.			
Size	EDGE Code	DWG NO	Rev
A1		15313100	C
Date	03 mch 09 / 2025	DEM	1 of 1

RECIRCULATION MOTOR


CUTTING MOTOR


SEE DRAWING NO. 15333120

SEE DRAWING No. 15113470

Rev: A New: Z1 FILTERS Wire connection change M1 V2 BK4 -> BK3, M1 V2 BK5 -> BK4, M1 L2 BK5 -> BK3.		SYRDIERS AVS VALDREISALLE 178 DK-2610 RØDEDØVRE DENMARK PHONE +45 2670 3500	
- CIRCUIT DIAGRAM -		UNIFORMS POWER CIRCUIT DIAGRAM	
Size: A3	CASE Code:	DWG NO:	Rev:
Scale:	Date: Tuesday, August 01, 2000	153331 10	1 A
Scale:		29 JAN 97 / JTV	Sheet 1 of 1

Rev. A: (FTH 13-12-2004) Ref. G1 changed to QS1		STRUERS A/S VALHOEJS ALLE 176 DK-2610 ROEDOVRE DENMARK PHONE: +45 3670 3500	
CIRCUIT DIAGRAM:		POWER CIRCUIT DIAGRAM UNITOM-5, - 50:	
Wednesday, January 29, 1997	Size A2	CAGE Code	DWG NO 15313110
Monday, December 13, 2004	Scale	08.03.99 / DEM	D.E.M. Sheet 1 of 1


SEE REMAINING 153131410


Rev.C: (03-01-2005 OH)
 - removed YM1, X5
 - removed wires BU13, BU14, BU15, BK, BN, VT
 - changed: S2:21 to 12, S2:22 to 11, S2:11 to 22, S2:21 to 21
 - removed wires BU11 + RD, BU12 + BU
 - new wires: X4:3 to S2:A2 <BU112>, X4:1 to S2:A1 <BU111>
 - RD6 changed to RD106
 - removed X2:3 and wire RD5
 - new wire: X6:2 to S2:21 <RD105>

Rev. B. Colors on BQ1 & BQ2 changed		Sivners A/S Vilhøjs Alle 176 DK 2610 Roedovre Denmark Phone +45 3670 3500	
		Control Circuit Diagram	
		Unitom - 50: With Excit. Table.	
Size A2	CAGE Code EK	DWG NO 15313120	Rev C
Tuesday, January 03, 2006	Scale	08.03.99 / DEM D.E.M.	Sheet 1 of 1


POWER SUPPLY: +5V, +12V AND +32V


TASTATUR KEYBOARD


DISPLAY INTERFACE


CONTROL INPUT FOR ARM POS.


CUT-MOTOR CURRENT:


CPU MEMORY


DELAY 50ns







AD CONVERTER


ANALOG TO DIGITAL INPUT:

Page 1	Rev. A Change R48 150K -> 120K	Siemens AG	Order No.
Rev. B	R49 200K -> 180K; (24-09-01 DEM)	Order No. 6ES7 312-1EG00-0AB0	175
Rev. C	(14-10-2004)	Order No. 6ES7 312-1EG00-0AB0	175
Rev. D	R77-R80 changed from 1k to 150R	Order No. 6ES7 312-1EG00-0AB0	175
UNITOM - 5 - 50:			
STEP-MOTOR CONTROL AND DRIVER:			
Tuesday, May 25, 1999	Rev. AIP	15313200	C
Tuesday, August 10, 2004	Rev. B	2.000.36.100V.DEM.	2 of 2

SEE CONNECTIONS
DIAGRAM 15113470


SEE DIAGRAM
15333460

SEE DIAGRAM
15333460




TO X-TABLE MOVER.

See Drawing
15333460

SLICING VALVE

CUTTING ARM


ALL WIRINGS 0.75mm², EXCEPT OTHERWISE MARKED.
 COLOR CODES:
 BK = BLACK
 BRN = BROWN
 RD = RED
 ORG = ORANGE
 YE = YELLOW
 GN = GREEN
 BU = BLUE
 VT = VIOLET
 GY = GREY
 WH = WHITE

TO THE FRONT PANEL:


Rev. A WHS YE-GN-DW-VE RPT-PHNS-VBECOR Rev. G. 24.01.02 Mod. 02 Rev. C. Reengineering Rev. SC03.1/3.2/4.0/4.1		STRUERS JOB VALUOR ISALE 175 BK-QVAL-ROE-DOVE EPSCONE-ES-3030 3000 WIRING DIAGRAM	
SHEET A2	CODE CODE 15333450	DWG NO 15333450	Rev. D
TUESDAY, AUGUST 01, 2000 Scale	25 JAN 27, 2000 Sheet 1 of 1	BACK-PANEL, CONTROL	

ALL WIRINGS 0.75mm² EXCEPT OTHERWISE MARKED.
 COLOR CODES:
 BK = BLACK
 BR = BROWN
 RD = RED
 OG = ORANGE
 YE = YELLOW
 GN = GREEN
 BU = BLUE
 VT = VIOLET
 GR = GREY
 WH = WHITE



SEE CONNECTIONS DIAGRAM 15113470


SEE DIAGRAM 15313465

SEE DIAGRAM 15313465


ONLY AT UNIT UNITOM-50:

TO X-TABLE MOVER.


ONLY AT UNIT UNITOM-5:


COOLING VALVE

CUTTING ARM

W4

T2
ITEMS
1 Amp. M1

W2

K4

K3


K2

K1

FT1

X1:1

J4 REAR CHASSIS


TO CONTROL BOARD FRONT PANEL

TO CONTROL BOARD FRONT PANEL

TO THE FRONT PANEL:
 SEE DRAWING NO.
 15313475

- Rev. D: (03-01-2005 CH)
- removed S2
- removed RD01, RD02, RD03, RD04, RD05
- removed segment of wire RD06 (J2-3-J2-5)
- changed label of RD06 to RD106
- changed wiring of BU01
- changed label of BU01 to BU101

Rev. B: Test Change BU4 → WH, BU5 → BK (25-01-01 BEM) Rev. C: BU1 & BQ2 colors changed. (25-12-04 KFC)		SINERG A/S Pallevej 14 DK-2750 Ballerup Copenhagen N, DENMARK	
Wiring Diagram			
Unitom - 5 & 50: BACK-PANEL, CONTROL.			
Sheet A1	EGG Code EL	Proj. No. 15313455	Rev. D
Tuesday, January 03, 2006		Scale 20 mA/J. 99 / DEM. BEM.	Sheet 1 of 1



SEE CONNECTIONS
DIAGRAM 151347.0

COLOR CODES: BK=BLACK, BN=BRNWN, RD=RED, OG=ORANGE, YE=YELLOW, GN=GREEN, BL=BLUE, VT=VIOLET, GY=GREY, WH=WHITE

Must be shown circuit protected with external fuses according to local regulations.
I_{th} max. = 145 A/m. In 5 sec.
Fuses max.:
200 - 240 V_o k 25 Amp. type gl.
380 - 480 V_o k 25 Amp. type gl.

WIRING AWG
200-230V - 3,3mm² 12AWG
380-480V - 2,6mm² 13AWG

Rev. A: 14.08.11.100 Rev. B: Wire connection change M1 U2 BK4 -> BK5, M1 U3 BK5 -> BK4, M1 U2 BK5 -> BK3, M1 Thermal Switch Symbol, NO -> NC.		STYRERS AS VALANCE SCALE 570 DK-25.10 RICE DOVRE SEWMAK PHONE +45 3070 3500	
WIRING DIAGRAM: UNITOM - 5			
Wiring of Power connection.			
Sheet: AZ	ENIG Code:	ENIG NO:	15333460
Title: Tuesday, August 01, 2000		Scale:	17.02.97 Dan Modig

SEE CONNECTIONS DIAGRAM 1513465

Must be short circuit protected with external fuses according to local regulations.
 In plus: + 100 Amp. In 5 sec.
 Fuses max:
 200 - 240 Vol.: 25 Amp. 1500 gL
 380 - 480 Vol.: 20 Amp. 1500 gL

WIRE GAUGE
 200-220V ~ 3.3mm GAWS
 380-480V ~ 2.1mm 14AWG

- COLOR CODES:
- BK = BLACK
 - BN = BROWN
 - RD = RED
 - OG = ORANGE
 - YE = YELLOW
 - GN = GREEN
 - BL = BLUE
 - VT = VIOLET
 - GY = GREY
 - WH = WHITE



M1
 CLIPPING MOTOR
 3 AC
 4500W

M2
 RECIRCULATION

Rev: A	Wire Gauge 2,5/13 -> 2mm/14AWG	Siemens ACS	Upholster Able 175	Rev: A
		Siemens	DK 2010 Protection	
		Phoenix	451 3870 3000	
Circuit Diagram		Unitom -5 & 50:		
		BACK-PANEL, POWER CONNECTION.		
Rev: A2	SK	DRG NO	15313465	Rev: A
DRG	SK	10.24.01 DEM	2.5.01	1 of 1
Wednesday, November 29, 2000				


TO THE BACK PANEL:
SEE DRAWING NO. 15313455

S2
Safety Interlock Switch
AZM 170-02zrk 24V

* ONLY AT UNIT UNITOM-50:
WIRE STRAP:
CODE SIGN. FOR
"EXCICUT TABLE CONNECTION"

EMERGENCY STOP
S1





ALL WIRINGS 0,75mm². EXCEPT OTHERWISE MARKED.
COLOR CODES:
BK = BLACK
BN = BROWN
RD = RED
OG = ORANGE
YE = YELLOW
GN = GREEN
BU = BLUE
VT = VIOLET
GY = GREY
WH = WHITE

Ps.: The drawing show the wiring seen from the front side.

Rev.B: (03-01-2005 OH)
- removed YM1, X5, added S2
- changed wiring of X8, X6, X4
- removed wires BU13, BU14, BU15, BK, BN, VT
- removed BU, RD, BU11 changed to BU111, BU12 to BU112
- BU02 to BU102, RD06 to RD106, BU01 to BU101, RD05 to RD105
- changed labels of wires: X3:2 BU1 to BU101, X2:4 RD6 to RD106

Rev.A: (07-01-2005 FTH) - Term.no. on S1 changed, - Additional text on Z1		Struers A/S Vejhoj Alle 176 DK 2610 Rødovre Denmark Phone +45 3670 3500	
Wiring Diagram		Unitom - 5 & 50: Front Panel; Control PCB.	
Size A2	CAGE Code EL	DWG NO 15313475	Rev 1 B
Tuesday, January 10, 2006	Scale	08.12.98 / DEM D.E.M.	Sheet 1 of 1

YM1
SAFETY GUARD LOCK

ALL WIRINGS 0.75mm², EXCEPT OTHERWISE MARKED.


- COLOR CODES:
- BK = BLACK
 - BN = BROWN
 - RD = RED
 - OG = ORANGE
 - YE = YELLOW
 - GN = GREEN
 - BU = BLUE
 - VT = VIOLET
 - GY = GREY
 - WH = WHITE

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Herewith declares that

<i>Product Name:</i>	Unitom-5
<i>Type No:</i>	533
<i>Machine Type:</i>	Cut-off machine with automatic cutting table

is in conformity with the provisions of the following directives:

Safety of Machinery 98/37/EEC according to the following standard(s):
EN ISO 12100-1:2005, EN ISO 12100-2:2005, EN ISO 13850:2006, EN 60204-1:2006.

EMC-Directive 2004/108/EEC according to the following standard(s):
EN 61000-6-1:2002, EN 61000-6-3:2002.

Low Voltage Directive 2006/95/EEC according to the following standard(s):
EN 60204-1:2006.

Supplementary Information The equipment complies with the American standards:
FCC part 15, subpart B and UL508.

The above has been declared according to the global method, module A

Date: 27.10.2008

Christian Skjold Heyde,
Vice President, R& D and Production, Struers A/S

Hersteller
Struers A/S
Pederstrupvej 84
DK-2750 Ballerup, Danmark
Telefon +45 44 600 800

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Christian Skjold Heyde,
Vice- President, R& D et Production, Struers A/S

UNIEL

Automatic Cutting Table for Unitom-5
Instruction Manual

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Automatic Cutting Table Safety Precaution Sheet

To be read carefully before use

1. The operator(s) should be fully instructed in the use of Unitom-5 and the automatic cutting table according to the Instruction Manuals.
2. The automatic cutting table must only be installed by a qualified technician.
3. The workpiece must be securely fixed in the quick-clamping device or the like.
4. Do not touch anything inside the cutting chamber while positioning the automatic cutting table.

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User's Guide

Table of Contents	Page
1. Getting Started	
Introduction	2
Checking the Contents	2
Installation of Automatic Cutting Table	2
Software Settings.....	3
Reading the Display	3
Reading the Display during the Cutting Process.....	4
Programming the Automatic Cutting Table	5
Positioning the Automatic Cutting Table.....	5
Clamping the Workpiece	6
Positioning the Cut-off Wheel	6
Using the Automatic Cutting Table as a Stationary Cutting Table	6

1. Getting Started

Introduction

The automatic cutting table is an optional accessory for Unitom-5. The Instruction Manual of the automatic cutting table is an addition to the Unitom-5 manual. The Instruction Manual explains the operation of a Unitom-5 with an automatic cutting table.

<p>Before reading this manual, please read the User's Guide section in the Unitom-5 manual.</p>

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- 1 Automatic cutting table (UNIEL)
- 2 Shims 0.1 mm and 0.2 mm
- 4 Cover plugs for threads in the bottom of the Automatic cutting table

Installation of Automatic Cutting Table

The automatic cutting table is mounted in the cutting chamber of Unitom-5, replacing the left cutting table.

- Turn off the main switch.
- Remove the left-hand side cutting table. Watch out for the shims placed under the cutting table in the bottom cabinet. The shims are to be re-used later for levelling the automatic cutting table.
- Place the automatic cutting table in the cutting chamber without fastening it.
- Remove the cover from the socket for connection of the automatic cutting table to Unitom-5. The socket is placed in the cutting chamber behind the cutting motor, near the back plate.
- Lead the cable for the automatic cutting table under the cutting motor, parallel to the back of Unitom-5. Mount the plug in the socket. Make sure that the cable is not jammed when the cutting motor moves up and down.
- Turn the main switch on again. The display on Unitom-5 now includes the positions for the automatic cutting table, to be controlled from the front panel on Unitom-5.
- Check with a ruler that the top side of the automatic cutting table is parallel to the fixed cutting table. The top side of the automatic cutting table should be 0.1 - 0.2 mm above the surface of the fixed cutting table. Level the automatic cutting table with the shims. The legs of the automatic cutting table must be firmly planted on the contact face of the bottom cabinet.


- Place the four screws used for the left-hand side cutting table loosely in the automatic cutting table.
- Align the edge of the automatic cutting table to be parallel to the fixed cutting table within ± 0.1 mm.
- Move the automatic cutting table using the controls on the Unitom-5 front panel, for free access to the screws in the automatic cutting table.
- Fasten the screws and check that the alignment is still within the limits.
- Mount cover plugs in any unused threads in the bottom.
- Mount the back stop at a distance of 190 mm from the front of the cutting table, at a right angle to the T-slot.
- Mount the quick clamping device at a distance of 25 mm in front of the back stop.

Software Settings

Unitom-5 automatically detects the presence of the automatic cutting table when the plug is connected to Unitom-5. The parameters for the automatic cutting table are integrated in the display on Unitom-5. The software settings for Unitom-5 will also be valid for the automatic cutting table.

Reading the Display

The display shows various kinds of information, for example parameters for an ongoing process. During the cutting process, the screen may look like this:

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Feed: The pre-set feed rate of the cut-off wheel.

No of cuts: the pre-set number of cuts. Please note that the first sample cut from the workpiece is not necessarily plane parallel or conforming to the pre-set thickness.

Sample width: the desired thickness of the sample. (Unitom-5 automatically compensates for the pre-set cut-off wheel thickness.)

Wheel thick.: the actual thickness of the cut-off wheel plus possible side run-out.

Load: load indicator of the cutting motor. The moving black column shows the status of the cutting motor:

Inside the white area: the cutting motor can be loaded continuously.

In the grey overload area: the cutting motor can only be intermittently loaded. The built-in load sensor might shut off the cutting motor.

Wheel Position: the actual position of the cut-off wheel counted from the top position of the cut-off wheel.

Stop Position: the pre-set stop position counted from the top position of the cut-off wheel.

Table Position: the actual position of the automatic cutting table. The table may travel for a distance of 0-90 mm, with the 90 mm position facing the cut-off wheel.


Cut no.: the actual number of completed cuts in a cutting sequence.

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
Programming the Automatic Cutting Table

The desired number of cuts, the sample width and the wheel thickness are pre-set before starting the cutting.

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
Positioning the Automatic Cutting Table

CUTTING		Load	POSITION	
Feed	0.50mm/s		Wheel	100.0mm
No of cuts	8		Stop	223.0mm
Sample width	3.0mm		Table	90.0mm
Wheel thick.	2.0mm		Cut no.	0
			SET STOP	
F2:Table			F3:ACT. F4:MAX	

↓
F2 Press F2, to toggle between SET STOP and AUT. TABLE in the display. Set to AUT. TABLE.

↓
 Use F3 <</F4>> to move the automatic cutting table.

The position of the automatic cutting table can be read on the display.

Place the automatic cutting table as close to the cut-off wheel as possible, yet leaving sufficient room for the desired number of samples to be cut.

CUTTING		Load	POSITION	
Feed	0.50mm/s	[Progress Bar]	Wheel	100.0mm
No of cuts	8		Stop	223.0mm
Sample width	3.0mm		Table	40.0mm
Wheel thick.	2.0mm		Cut no.	0
			HUI TABLE	
F2: Stop			F3: <<	F4: >>

Clamping the Workpiece

- Place the workpiece between the clamp and the back stop of the automatic cutting table. The protruding part of the workpiece should pass slightly across the fixed cutting table on the right-hand side, if possible. Make sure that the sample does not run into the fixed cutting table when it advances.
- Push the clamp against the workpiece and lock the quick-clamping device with the locking handle.

Do not clamp the sample with the fixed cutting table.

Positioning the Cut-off Wheel

- After clamping the workpiece, you may lower the cut-off wheel to a position at least 1 mm above the workpiece, to test the correct position of the cut on the workpiece.
- Press POSITION DOWN ▼ to lower the cut-off wheel.
- The cut-off wheel stops when you cease pressing the key.
- You can watch the position of the cut-off wheel in the display.

CUTTING		Load	POSITION	
Feed	0.50mm/s	[Progress Bar]	Wheel	150.3mm
No of cuts	8		Stop	223.0mm
Sample width	3.0mm		Table	40.0mm
Wheel thick.	2.0mm		Cut no.	0
			HUI TABLE	
F2: Stop			F3: <<	F4: >>

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Using the Automatic Cutting Table as a Stationary Cutting Table

If you do not want the automatic cutting table to move, set the number of cuts to 1 and move the automatic cutting table to the right-hand position.

CUTTING		Load	POSITION	
Feed	0.50mm/s	[Progress Bar]	Wheel	100.0mm
No of cuts	1		Stop	223.0mm
Sample width	3.0mm		Table	90.0mm
Wheel thick.	2.0mm			
			SET STOP	
F2: Table			F3: ACT.	F4: MAX

Reference Guide

Table of Contents	Page
1. Trouble-Shooting	8
2. Technical Data	8
3. Accessories	8

1. Trouble-Shooting

Error	Explanation	Action
Samples not cut to set thickness +/- 0.2 mm.	The sample moves during cutting.	Clamp the sample more solidly.
	The cut-off wheel bends.	Use a lower feed speed.
	Width of cut-off wheel not taken into consideration.	Measure the width of the cut-off wheel and compensate for the actual width.
	The movement of the automatic cutting table is blocked.	Remove possible obstructions underneath the automatic cutting table and clean around the table. Check that the workpiece is not clamped with the fixed cutting table

2. Technical Data

Subject		Specifications	
		Metric/International	US
Sample	Number of cuts:	1-20	
	Thickness: adjustable in steps of:	1-90 mm 0.1 mm	0.04"-3.54" 0.004"
Positioning	Positioning range at right angles to the cut-off wheel- adjustable in steps of:	90 mm 0.5 mm	3.6" 0.02"
	Max. positioning speed	4 mm/s	0.16 "/s

3. Accessories

Specification	Code
Replacement Stainless Steel Bands For Automatic Cutting Table (UNIEL) and Cutting Table with 10 mm T-slots (UNIFT). Set of three.	UNIBA
Quick-clamping Device for 10 mm T-slots, Left For securing the workpiece. Complete with back stop. To be mounted on a cutting table placed at the left hand side of the cut-off wheel.	UNILE
Clamping Tools for 10 mm T-slots For clamping complicated workpieces to the cutting table on Unitom-5 and Unitom-2. Complete with clamps, supports and bolts.	UNIOF

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Christian Skjold Heyde,
Vice- President, R& D et Production, Struers A/S

Unitom

Automatic and Manual Cut-off machines for large and difficult workpieces up to 120 mm (4.7") in diameter

*The Unitom range:
Choose the model to suit
your requirements*

Unitom is the family of Struers materialographic cut-off machines sharing a platform with a 350 mm cut-off wheel and a powerful 6.4 HP cutting motor. From the manually operated Unitom-2 to the top Unitom-50, with the unique ExciCut feature, the Unitom range will match all cutting requirements in the mid-size.

3 models

Unitom-2: Manual operation
Unitom-5: Automatic operation
Unitom-50: Automatic operation and ExciCut

Specifications

- 4.7 kW (6.4 HP) cutting power
- Cut-off wheel: Max diameter: 350 mm (13.7")
- Cutting capacity: Max cutting diameter: 120 mm (4.7")

The Unitom family is based on

proven, superior cutting technology. The Unitom concept of high quality and high capacity combined with ease-of-use has been extremely successful.

Choose the Unitom that fits your requirements:


The manually operated Unitom-2 for the limited volume or the Unitom-5 for high volume automated cutting.

If you are dealing with hardened materials, our Unitom-50 with its unique ExciCut feature, is the right choice.

Unique features

- *ExciCut:* Cutting of workpieces with ultrahigh hardness (Unitom-50)
- *Automatic Feed Reduction:* Fastest cutting without damaging sample and cut-off wheel
- *AutoStop:* Permits unattended cutting
- *Constant Feed Speed:* Patented cut-off wheel movement
- *Automatic Table:* For cutting multiple identical slices (Unitom-5)

Built-in flush hose with adjustable nozzle for convenient cleaning

ExciCut


-for cutting of the hardest materials

Unitom-50 comes standard with the unique ExciCut feature. During cutting, the cutting table moves slightly forwards and backwards. This reduces the contact area between cut-off wheel and workpiece, and enables troublefree cutting of very hard materials in very large sizes. An added benefit is better access of the coolant into the cut, eliminating the risk of thermal damage.

The ExciCut technology was invented by Struers more than a decade ago, and appears on Unitom-50 in a highly refined version.

Direct Cut

-the classic cutting mode

For cutting of less hard materials, the Direct Cut mode is ideal. This classic cutting mode moves

the cut-off wheel into the workpiece in a linear movement.

Automatic Feed Reduction

-fastest cutting without damaging sample and cut-off wheel

On an ordinary automatic cut-off machine, feed speed or force is preset and kept constant throughout cutting. Struers have proven that maintaining a constant feed speed offers far better results than constant force. However, if the preset feed speed is too high in relation to the workpiece hardness or size, there is a risk of thermal damage. To overcome this, Unitom-

Manual cutting on the Unitom-2

5 and -50 have been fitted with an ingenious feature, Automatic Feed Reduction, which will automatically reduce the feed speed and proceed to cut with the reduced speed.

In this way, cutting is performed as fast as possible, without risking damage to sample or cut-off wheel.

AutoStop

-for unattended cutting

Conduct your cutting unattended, and let Unitom-5 and -50 stop automatically.

Beyond being able to work with preset stop positions, Unitom has an AutoStop feature, where it automatically detects when the workpiece has been cut through. Thereafter the cutting motor and coolant are stopped and the cut-off wheel retracted. This allows completely unattended cutting.

Automatic Cutting Table

-for cutting of multiple, uniform slices

Unitom-5 can be equipped with a motorized table for cutting plane parallel sections. The table makes it possible to cut a number of samples without any manual operation. Programming and operation of the automatic table is an integrated part of the menu control of Unitom-5. Thickness and number of cuts can be programmed.

Constant Feed Speed

-patented cut-off wheel movement

A unique, patented method of moving the cut-off wheel dramatically increases the quality of a cut. Struers has proven that this method effectively counters the extreme forces applied to a workpiece when a cut-off wheel enters and exits a workpiece.

Rapid Advance

-fast contact with the workpiece

Having started Unitom-5 or -50, the cut-off wheel can be advanced with increased speed until it contacts with the workpiece. Unitom then automatically changes to the preset speed.


Continue Function

-stop and resume as you please

If cutting has been temporarily stopped, it is possible to press "Continue", and the cut-off wheel will move quickly to the position where the cutting was

Unitom-5 is available with automatic cutting table for cutting of multiple, uniform slices

Unlimited access to the spacious cutting chamber

stopped, and continue cutting at the preset feed speed.

Programmable Return Position

-automatic cut-off wheel positioning

The return position of the cut-off wheel can be set to three different positions: Unitom-5 and -50 retract the cut-off wheel to the top position, to the position where the cutting was started or let the cut-off wheel stay in the cut. In this way it is very easy to get convenient access to the workpiece or a quick turn-around when cutting many uniform workpieces.

Cutting Large Workpieces

If the user wants to cut large or protruding workpieces, it is possible to replace one or both of the side walls on the machine with optional rubber curtains.

User-Friendly Design

Extra large windows and fully, indirect lightning, give an excellent, non-glare view of the cutting chamber. A balanced mechanism ensures effortless opening and closing of the protection guard.

Corrosion-protected Cutting Chamber

All metal parts in the cutting chamber are made of either stainless steel or surface protected materials, practically eliminating corrosion. The top surface of the cutting table is exchangeable and made of stainless steel. A powerful lamp lights up the cutting chamber,

Unitom on optional LABUN table

and a flush hose with adjustable flow takes care of fast and efficient daily cleaning.

Simple Installation

Installation and levelling of the machine are made simple, by virtue of four height adjustable vibration dampers.

Safety Features

Unitom complies with international safety standards, and safety features include emergency stop and cut-off switch for cutting motor and all machine functions. The protection

Simple operation via large display

guard cannot be opened before the cut-off wheel has stopped. The windows of the protection guard are made of impact resistant material. Unitom is prepared for connection to an exhaust system, as some cutting materials may emit harmful gases or irritating fumes.

Special Table Unit Available

An optional strong table unit, matching the look of Unitom, has been specially designed to hold the recirculation cooling unit, and to position the machine at the right working height.

A Sturdy Workhorse

Unitom will enable the user to cut a large number of samples, hour after hour. All parts of the machine subject to stress have been put through a careful life-

Cutting of very long, protruding workpieces using optional rubber curtains

time test. The rotating parts of the machine have been dynamically balanced to ensure stability and to maximize the working life of your Unitom.

Accessories

The Unitom family comes with an extensive range of accessories. A combined table and cabinet is available for bringing the cut-off machine into an ideal working height. Quick-clamping devices and vertical clamping devices ensure full flexibility in mounting.

Please refer to complete list in this brochure.


Cut-off Wheels

Unitom is designed for cut-off wheels of 300 and 350 mm diameter. For normal use, the 300 mm wheels (with a thickness of 2 mm) are recommended. For large workpieces, the 350 mm wheels (with a thickness of 2.5 mm) are used, allowing Unitom to cut samples up to 120 mm diameter.

Unitom-50 with its unique ExciCut, will cut through even large workpieces with 350 mm cut-off wheels.

A range of Struers cut-off wheels is available for cutting any material.

Please refer to table (right).



Even very bulky workpieces can be cut - here using quick-clamping devices (optional)

Cut-off wheels

A range of cut-off wheels for all ferrous and non-ferrous metals has been specially developed for Struers cut-off machines. Diamond cut-off wheels are available for cutting of ceramic and mineralogical materials and sintered carbides.

Application	Dimensions	Grinding agent	Code
For ExciCut and direct cutting of extremely hard steels (HV 350-800)	300 x 2 x 32 mm	Al ₂ O ₃	51UNI
	350 x 2.5 x 32 mm	Al ₂ O ₃	81EXO
For ExciCut and direct cutting of steels (< 500 HV). Extremely hard (white cast iron) or ductile (18/8 or St60) steel types are cut with ExciCut. Ensures good wheel economy	300 x 2 x 32 mm	Al ₂ O ₃	53UNI
	350 x 2.5 x 32 mm	Al ₂ O ₃	83EXO
Universal cut off-wheel for steel-alloys (< 700 HV). For ExciCut only	350 x 2.5 x 32 mm	Al ₂ O ₃	8UEXO
For ExciCut and direct cutting of steels (< 500 HV). Ensures very good wheel economy	300 x 2 x 32 mm	Al ₂ O ₃	54UNI
	350 x 2.5 x 32 mm	Al ₂ O ₃	84EXO
For ExciCut and direct cutting of non-ferrous metals (< 300 HV). Primarily for ExciCut	300 x 2 x 32 mm	SiC	56UNI
	350 x 2.5 x 32 mm	SiC	86EXO
For direct cutting of titanium and its alloys, as well as other non-ferrous metals (< 400 HV)	300 x 2 x 32 mm	SiC	56UNI
	350 x 2.5 x 32 mm	SiC	90EXO
For direct cutting of sintered carbides and ceramics	305 x 1.8 x 32 mm	Diamond	25EXO
	350 x 1.8 x 32 mm	Diamond	26EXO
Ceramics and minerals	305 x 1.5 x 32 mm	Diamond	55UNI
	356 x 1.5 x 32 mm	Diamond	85EXO
Plastics and resins, mounted materials	300 x 1.5 x 32 mm	Diamond	59UNI
	356 x 1.5 x 32 mm	Diamond	89UNI
Extremely hard ferrous metals	305 x 1.8 x 32 mm	CBN	58UNI
	356 x 1.8 x 32 mm	CBN	88EXO

Cutting Capacity

Technical Data

Subject	Specification	Metric/International	US
CUTTING SPECIFICATIONS			
Workpiece Dimensions (Max)	<i>Workpiece inside cutting chamber:</i>		
	Height	120 mm	5.9"
	Width	415 mm	16.3"
	Depth	630 mm	24.8"
	<i>Workpiece protruding cutting chamber:</i>		
	Height	120 mm	5.9"
Depth	160 mm	6.3"	
Cutting Capacity*	Max. cutting diameter	120 mm	4.7"
	<i>Rectangular sample**:</i>		
	Height	90 mm	3.5"
	Depth	200 mm	7.9"
PHYSICAL SPECIFICATIONS			
Cutting Motor	Cutting power	4.7 kW	6.4 HP
	Max power (depending on voltage)	5.7 - 6.7 kW	7.6 - 9.0 HP
Cut-off Wheel	Diameter x Thickness x Centre-hole Rotational speed (running idle)	350 x 2.5 x 32 mm 2775 rpm	13.7 x 0.09 x 1.25" 2775 rpm
Positioning & Feed	Positioning range (of cut-off wheel) (adjustable in steps of)	0 – 223 mm 1 mm	0 – 8.75" 0.04"
	Max. positioning speed	30 mm/s	1.2"/s
	Feed Speed range	0.05 – 3 mm/s	0.002 – 0.118"/s
	(adjustable in steps of)	(0.05mm/s)	(0.002"/s)
	Cutting force	Max. 350 N	Max. 88 lbf
Cutting table	Width	286 mm	11.26"
	Depth	316 mm	12.44"
	T-slots	10 mm	0.39"
Dimensions and Weight	Height	660 mm	26"
	Width	700 mm	27.5"
	Depth	750 mm	29.5"
	Weight	162 kg	353 lbs
	<i>Optional Table Unit (LABUN):</i>		
	Height	795 mm	31.3"
	Width	700 mm	27.5"
Depth	750 mm	29.5"	
Recirculation Cooling Unit	Tank volume	65 l	17 gallons
	Approx. flow	12 l/min	3.2 gallons/min
ENVIRONMENTS			
Directives	Machinery: 98/37/EEC Low Voltage: 73/23/EEC EMC: 89/336/EEC		
Dynamic Balancing	Max. admissible unbalance U_{per} according to ISO 1940/1, Balance Quality Grade G6.03: $U_{per} = 2 \text{ gmm/kg} = 110 \text{ gmm}$.		
Noise level	Approx. 73 dB(A) measured at idle running, at a distance of 1.0 m/39.4" from the machine.		

*Projected cutting capacity with new cut-off wheel. Actual cutting capacity depends on sample material, cut-off wheel and clamping technique

**Example of cutting capacity. Please refer to Cutting Diagram above

Struers A/S
 Pederstrupvej 84
 DK-2750 Ballerup, Denmark
 Phone +45 44 600 800
 Fax +45 44 600 801
 struers@struers.dk
 www.struers.com

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Specification	Code
Unitom-2 Manual Cut-off machine. 4.7 kW (6.4 HP) motor. For 350 mm cut-off wheels. With recirculation cooling unit (capacity 65 l) and two cutting tables with 10 mm T-slots. Table unit and clamping devices are ordered separately.	UNIW2
Unitom-5 Automatic Cut-off machine. Automatic feeding with electronic control of feed speed. 4.7 kW (6.4 HP) motor. For 350 mm cut-off wheels. With recirculating cooling unit (capacity 65 l) and two cutting tables with 10 mm T-slots. Table unit and clamping devices are ordered separately.	UNIW5
Unitom-50 Automatic Cut-off machine with ExciCut. Automatic feeding with electronic control of feed speed. 4.7 kW (6.4 HP) motor. For 350 mm cut-off wheels. With recirculating cooling unit (capacity 65 l) and two cutting tables with 10 mm T-slots. Table unit and clamping devices are ordered separately.	UNIWL
Accessorie	Code
Automatic Cutting Table (for Unitom-5) Motorized table for cutting plane parallel sections. Positions the workpiece at a right angle to the cut-off wheel. The table is controlled through Unitom-5	UNIEL
Table Unit For support of Unitom-machines. With compartment for recirculation cooling unit. Width: 900 mm, depth: 750 mm, height 800 mm. When recirculation cooling unit is used, an optional rolling pallet (TREJU) is ordered separately	LABUL
Rolling Pallet For Unitom recirculation cooling unit	TREJU
Cutting table for 10 mm T-slots (for Unitom-2/-5) Additional table for support of large workpieces. With exchangeable stainless steel bands. 140 x 315 mm (5.5" x 12.4")	UNIFT
Replacement Stainless Steel Bands For Cutting Table with 10 mm T-slots (UNIFT) and Automatic Cutting Table (UNIEL). Set of three	UNIBA
Quick-clamping Device for 10 mm T-slots, Left For securing the workpiece. Complete with back stop. To be mounted on a cutting table placed at the left hand side of the cut-off wheel.	UNILE
Quick-clamping Device for 10 mm T-slots, Right For securing the workpiece. Complete with back stop. To be mounted on a cutting table placed at the right hand side of the cut-off wheel.	UNIRI
Vertical Clamping System for 10 mm T-slots For clamping irregularly shaped workpieces. Clamping height adjustable up to 58 mm. Complete with operating key and one flat swivel shoe.	TREVS
Riser Block for 10 mm T-slots For elevating Vertical Clamping System (TREVS) when clamping high workpieces. To increase clamping height by 60 mm.	TREKS
Swivel Shoes for Vertical Clamping System Set of 4 multi-shaped swivel shoes for (TREVS).	TREVI
Rubber Curtains To be mounted on Unitom machines for cutting of long, protruding workpieces. Set of two	UNIRU
Flange Set, 119 mm dia. Alternative to 100 mm standard flange. For use with 350 mm cut-off wheels on Unitom cut-off machines.	UNIFL

Struers' products are subject to constant product development. Therefore, we reserve the right to introduce changes in our products without prior notice.

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Struers Ltd.
Struers Inc.