Metason

Instruction Manual
Gebrauchsanweisung
Mode d'emploi

Manual No.: 15147001

Date of Release FF.01 .201F

*Metason
Instruction Manual*

Instruction Manual

Table of Contents	Page
User's Guide	1
Reference Guide.....	7

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the Serial No. and Voltage on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

Original instructions. The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2011.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Fax +45 44 600 801

Metason Safety Precaution Sheet

To be read carefully before use

1. The operator should be fully aware of the use of the machine according to the Instruction Manual. The machine must be placed in an adequate working position.
2. Be sure that the actual voltage corresponds to the voltage stated on the back of the machine. Do not operate the cleaner without proper electrical grounding.
3. Do keep the control panel and the area around the cleaner clean and wipe up solution which spills over the tank brim with a dry cloth.
4. Do not use inflammable liquids in the tank. Do not use solvents in the tank.
5. Do not use mineral acids or bleaches. These could damage the tank.
6. Do not allow the level of liquid to drop below 25 mm (1") from the top of the tank.
7. Do not immerse the ultrasonic cleaner in water.
8. Do not place your fingers in the ultrasonic cleaner when it is in operation. Do not touch the stainless steel tank or cleaning solution - they may be hot.
9. Do not operate the ultrasonic cleaner with the cover open or with a solution temperature higher than 70 °C (160 °F).
10. Do not allow unauthorised personnel to disassemble the cleaner. High voltage is present in the electronic components creating a possible shock hazard.

The equipment should only be used for its intended purpose and as detailed in the Instruction Manual.

The Metason ultrasonic cleaner is designed to operate with aqueous solutions and for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician.

User's Guide

Table of Contents	Page
1. Installation	2
2. Operation	2
Introduction	2
Construction	2
The Ultrasonic Cleaning Principle	3
Cleaning	4
Direct Cleaning	4
Indirect Cleaning	4
Start Cleaning	4
Cleaning with Heat (Metason 200HT)	4
Safety Precautions	5
Hints	5
Temperature	6
Cleaning Liquids	6


1. Installation

- Start with a clean tank and clean accessories.
- Check that the voltage printed on the model plate is identical to your supply voltage.
- Fill the tank with warm water. The liquid should come to 25 mm (1") from the top of the tank when the tray/beakers and items are submerged.
- Turn the TIMER to 5-10 min and set Metason to run, to allow the solution to degas.

2. Operation

Introduction

3 standard models are available: Metason 60T, Metason 120T, and Metason 200HT.

Construction

The cleaners consist of a generator and a stainless steel tank, all in one unit. In the generator a high frequency voltage is generated, which is transmitted to one or several transducers placed under the tank. These transducers are piezoelectric crystals, which contract and expand, i.e. oscillate, with the frequency of the applied voltage. This oscillation is then conveyed to the tank and thus to the liquid.

A 60 min mechanical timer is standard on each model. Metason 200HT is also equipped with heating, which increases the chemical activity of the cleaning solution. The heating is adjusted to 60 °C and automatically switches off if the tank runs out of cleaning solution.

Metason 60T and Metason 120T are designed for single specimens and Metason 200HT is designed for both single specimens and specimen holders.

**The Ultrasonic
Cleaning Principle**

The ultrasound which is transferred to the liquid causes the propagation of pressure waves in the tank. These pressure waves involve a continuous change of the pressure in the liquid so that, as a result of the sub-pressure, innumerable microscopic evacuated bubbles are formed which will soon expand and implode (cavitation).

These implosions give the cleaning liquid a high speed locally so that dirt and loose particles are removed from the specimens which are immersed into the liquid.

Cleaning

Direct Cleaning

The direct method is satisfactory for many cleaning jobs. The cleaning solution is filled directly into the tank and the items are placed in a perforated tray and lowered into the tank (specimen holders are placed in the integrated grip in the cover). The items can also be suspended on a wire or a hook and submerged in the tank. This method is simple and very effective. However, this method implies that all removed material remains in the tank, and that only one solution can be used at a time.

Indirect Cleaning

When using this method the tank is filled with a liquid driving medium such as water with a wetting agent (liquid soap) to improve cavitation.

The cleaning solutions are filled into beakers or solid insert trays. These beakers or trays are then suspended in the tank.

The indirect method enables the simultaneous use of several completely different cleaning solutions.

Removed material stays in the beaker/tray and can easily be examined, filtered or discarded.

The indirect method implies a reduction of the ultrasonic power compared to the direct method.

Start Cleaning

Cleaning with Heat (Metason 200HT)

- Set the TIMER for the amount of time you wish the specimen to be cleaned.
- Turn the HEAT on.
- Set the TIMER for the amount of time you wish the specimen to be cleaned.

NOTE

To stop Metason at any time turn the TIMER to zero.
When the specimens are cleaned turn the HEAT OFF.

Safety Precautions

Your Metason ultrasonic cleaner is designed to operate with aqueous solutions. Use tap water or distilled water and add a cleaning agent such as Struers Concentrated Soap Solution.

Solvents and other organic liquids are traditionally used in cleaning operations. However, such liquids should not be used since they imply an occupational hazard and may damage the cleaning unit.

The level of the cleaning liquid should come to 25 mm (1") from the top of the tank. Do not operate the unit with an empty tank.

Read and follow the instructions listed in the yellow front page Safety Precaution Sheet.

Do not use toxic or inflammable cleaning liquids.

Hints

- The best cleaning solution is water with an added cleaning agent such as Struers Concentrated Soap Solution.
- Unplug the power cord before filling or emptying the tank. On units equipped with drains, open the drain plug and empty the liquid into a sink or a waste disposal container. Any liquid contains an amount of dissolved gas, which tends to absorb a part of the ultrasonic power. This dissolved gas should be stripped off whenever fresh liquid is filled in the tank. The stripping is done by simply running the ultrasonic cleaner for approx. 5 min. before starting the first cleaning cycle with fresh liquid.
- Excessive weight on the tank bottom dampens the ultrasonic energy and may cause damage to the transducer. Heavy samples must be placed with caution in the tank.
- It is faster and more efficient to run several small loads rather than a few big loads, i.e. use a large volume of liquid compared to the surface area of the parts being cleaned.
- Flush the cleaned items with Ethanol and dry with a hot air blower.
- The basket of Metason 120T fits into Drybox, so that a batch of specimens may be moved directly from cleaning to drying.

Temperature

The temperature of the aqueous solution in the tank will stabilise at approx. 50-60 °C when the ultrasonic cleaner is running continuously. Do not allow the temperature to exceed 70 °C.

Both heat and cavitation increase the chemical activity of the cleaning fluids. With Metason 200HT it is possible to stabilise the temperature of the cleaning solution. This improves the cleaning action but some materials may be damaged by this increased chemical activity. When in doubt, run test samples of items to be cleaned.

Cleaning Liquids

In order to obtain good cavitation which is essential for the cleaning mechanism, the cleaning liquid must have the following properties:

- The density of water, or slightly higher.
- Low vapour pressure at the required operating temperatures.
- Low viscosity at the required operating temperatures.

Reference Guide

Table of Contents	Page
1. Consumables and Accessories	8
2. Maintenance	8
3. Trouble-Shooting	9
4. Technical Data	10

1. Consumables and Accessories

Specification	Cat. Nr.
Concentrated Soap Solution for ultrasonic cleaning, 1l	49900000
<i>Basket for Specimens</i> for Metason 120T, 200HT (there is room for 2 baskets) and Drybox-2	04236904
<i>Drybox-2</i> Combined specimen drier and curing oven with adjustable temperature (approx. 40-90 °C) and low noise fan. With basket for specimens and slot for specimen holders 1 x 110-120 V / 50-60 Hz 1 x 220 V / 50 Hz 1 x 220 V / 60 Hz	03706118 03706133 03706134

2. Maintenance

The tank seldom requires hand cleaning because the ultrasonic action keeps it clean. If material is deposited on the tank bottom, it should be removed because it absorbs energy and reduces the cleaning effect. Use a non-abrasive cloth and water to clean the tank. Change the cleaning solution on a regular basis. Like most chemicals, solutions become spent over time.

Empty and clean the tank if the unit is out of operation over a long period.

The cabinet should be cleaned with a soft, damp cloth and a mild soap solution.

3. Trouble-Shooting

Error	Cause	Action
Unit does not start	Unit not properly plugged	- Check the line cord and the plug - Check the electrical supply outlet
	Mechanical timer not on	Turn timer clockwise
	Faulty switch	Call Struers service technician
	Faulty power supply fuse	Replace fuse
Decreased ultrasonic activity	Cleaning liquid has not been degassed	Let the unit run for 5 min
	Solution low or dirty	Exchange the solution
	Liquid level too low	Refill with liquid
	Solution cold	Allow solution to heat until activity peaks
	Tank bottom covered with soil particles	Clean tank with warm water and wipe with a non-abrasive cloth

WARNING

High voltage inside - dangerous shock hazard
Unqualified personnel should not disassemble or repair the unit.

4. Technical Data

Mains connection

220 V, 50/60 Hz.

Metason 60T Power Supply

80 W, HF-output power 35 W

Dimensions

Tank

Width: 140 mm
Depth: 150 mm
Height: 100 mm

Overall dimensions

Width: 250 mm
Depth: 290 mm
Height: 295 mm
Tank capacity 1.9 l
Weight: 5 kg

Metason 120T Power Supply

125 W, HF-output power 70 W

Dimensions

Tank

Width: 240 mm
Depth: 140 mm
Height: 100 mm

Overall dimensions:

Width: 300 mm
Depth: 335 mm
Height: 295 mm
Tank capacity 2.8 l
Weight: 6.4 kg

Metason 200HT Power Supply

469 W, HF-output power 140 W, Heating power 284 W

Dimensions

Tank

Width: 290 mm
Depth: 240 mm
Height: 150 mm

Overall dimensions

Width: 400 mm
Depth: 395 mm
Height: 370 mm
Tank capacity 9.5 l
Weight: 9.6 kg

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
Metason 200HT

Drybox-2

User-friendly accessories specially designed for materialographic specimens

Marking

Marking of specimens for subsequent identification

If a specimen is not properly marked, the preparation efforts may prove a waste of time if the sample can not be properly identified later.

- Universal engraving point, suited for the marking of practically all metals, plastics, ceramics and glass
- Very simple operation: You hold the Engraver like you would hold a pencil. The vibration amplitude is regulated by means of a control button

Cleaning

Effective ultrasonic cleaning of metallographic, ceramic and mineralogical specimens

- Fast and efficient cleaning
- Cleans pores and cavities
- Removes grinding and polishing debris and impurities
- Sturdy and reliable

The Metason ultrasonic cleaners produce a high-frequency voltage that is transformed into ultrasonic waves which spread through the cleaning liquid. These waves cause rapid formation and implosion of microscopic bubbles (cavitation) which provide the cleaning effect.

The Metason series consists of three versions with different capacities. A 60 min mechanical timer is standard on each model. Metason 200HT is also equipped with a heating which increases the effect of the cleaning solution. The temperature is adjustable to 60°C and automatically switches off if the tank runs out of cleaning solution. Specimen holders can be suspended into Metason 200HT, and individual specimens are placed in a wire basket for cleaning.

Drying

Drying of specimen holders and single specimens, and curing of cold mounted specimens

- Continuously variable adjustment of the air temperature between 40°C and 90°C
- Considerable reduction of curing time for cold mounted specimens
- Dries and cures up to 24 specimens per batch

Drybox-2

Combined drying and curing apparatus

Drybox-2 comprises a blower mounted in a steel case, producing a strong air flow that can be heated to about 90°C with a 2000/1500 W adjustable heating element. A maximum of 24 materialographic specimens may be dried in one batch.

Drybox-2 may considerably reduce the curing time for cold mounted specimens. Curing starts due to a chemical reaction when the two components are mixed. Heat advances the start of the reaction, but once initiated, the reaction requires cooling. Drybox-2 keeps the specimen temperature at the right level throughout the entire curing process.

Storing

Practical, convenient and space saving storage for materialographic laboratories

Struers has designed convenient and space saving storing cabinets for MD-consumables.

- Designed for materialographic laboratories

Storage cabinet for MD-consumables

The storage cabinets have shelves for storing each individual MD preparation surface. The storage cabinet allows the disc to dry effectively, and protects the preparation surface from contamination. Identification labels are delivered together with the storage cabinet. The labels can be mounted at the front of each shelf to allow individual disc information.

The cabinet version for MD preparation surfaces 350 mm dia. has shelves that pull out to the front. The shelves of the versions for MD preparation surfaces 300, 250 and 200 mm dia. can swing out either to the left or the right, so as to provide access, no matter where the storage cabinet is placed.

Struers A/S
 Pederstrupvej 84
 DK-2750 Ballerup, Denmark
 Phone +45 44 600 800
 Fax +45 44 600 801
 struers@struers.dk
 www.struers.com

Technical Data

Subject	Specification		
Metason	Metason 200HT	Metason 120T	Metason 60T
HF-output power	140 W	70 W	35 W
Tank dimensions (W x D x H)	290 mm x 240 mm x 150 mm	240 mm x 140 mm x 100 mm	140 mm x 150 mm x 100 mm
Tank capacity	9.5 l	2.8 l	1.9 l
Overall dimensions (W x D x H)	400 mm x 395 mm x 370 mm	300 mm x 335 mm x 295 mm	250 mm x 290 mm x 295 mm
Weight	9.6 kg	6.4 kg	5.0 kg
Drybox-2			
Output	2000 W at 220 V power supply 1500 W at 110-120 V power supply		
Rate of airflow	5 m/s		
Dimensions (W x D x H)	290 mm x 292 mm x 172 mm		
Weight	4.5 kg		
Struers Engraver			
Strokes per second	100 at 50 Hz power supply 120 at 60 Hz power supply		

Specifications

Specifications	Cat. no.
Metason 200HT Apparatus for ultrasonic cleaning of specimen holders, 140, 160 and 200 mm dia. and for single specimens. With heating to 60°C and timer for up to 60 min. Tank capacity 9.5 l. Including cover equipped with integrated grip for specimen holders. Baskets for single specimens (04236904) are ordered separately. (There is room for 2 pcs.) 1 x 115 V / 50-60 Hz or 1 x 220 V / 50-60 Hz	05146333
Metason 120T Apparatus for ultrasonic cleaning of specimens. With timer for up to 60 min. Tank capacity 2.8 l. Including wire basket (04236904) for specimens and cover 1 x 115 V / 50-60 Hz or 1 x 220 V / 50-60 Hz	051462xx
Metason 60T Apparatus for ultrasonic cleaning of specimens. With timer for up to 60 min. Tank capacity 1.9 l. Including wire basket for specimens and cover 1 x 115 V / 50-60 Hz or 1 x 220 V / 50-60 Hz	051461xx
Wire Basket for Specimens For Metason 120T, Metason 200HT and Drybox-2	04236904
Drybox-2 Combined specimen drier and curing apparatus with adjustable temperature (approx. 40-90°C) and low noise fan. With wire basket for specimens (04236904) and slot for specimen holders 1 x 110-120 V / 50-60 Hz, 1 x 220 V / 50 Hz or 1 x 220 V / 60 Hz	037061xx
Struers Engraver For marking of specimens 1 x 110 V / 50-60 Hz or 1 x 220 V / 50-60 Hz	001461xx
Engraver Tip Tungsten carbide tipped point for Engraver	40900043
Storage Cabinet For 10 MD preparation surfaces in 350 mm dia.	05666001
MD-Concert Storage cabinet for 10 MD preparation surfaces in 300, 250 or 200 mm dia.	05306101
MD-Concertino Storage cabinet for 8 MD preparation surfaces in 200 mm dia.	05306102

MD-Concert and MD-Concertino

Struers' products are subject to constant product development. Therefore, we reserve ourselves the right to introduce changes in our products without notice

AUSTRALIA & NEW ZEALAND
Struers Australia
 27 Mayneview Street
 Milton QLD 4064
 Australia
 Phone: +61 7 3512 9600
 Fax: +61 7 3369 8200
 info.au@struers.dk

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CANADA
Struers Ltd.
 7275 West Credit Avenue
 Mississauga, Ontario L5N 5M9
 Phone +1 905-814-8855
 Fax +1 905-814-1440
 info@struers.com

CHINA
Struers Ltd.
 No. 1696 Zhang Heng Road
 Zhang Jiang Hi-Tech Park
 Shanghai 201203, P.R. China
 Phone +86 (21) 6035 3900
 Fax +86 (21) 6035 3999
 struers@struers.cn

CZECH REPUBLIC
Struers GmbH
 Organizační složka
 Havlíčkova 361
 CZ-252 63 Roztoky u Prahy
 Phone +420 233 312 625
 Fax: +420 233 312 640
 czechrepublic@struers.de

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Struers Inc.