

# SELECTION GUIDE FOR COLD MOUNTING ACCESSORIES

Specifications			Contents	Cat. no.
Mounting cups	<b>FixiForm</b> Two part mounting polypropylene (PP) cup to be used with all Struers cold mounting resins. Can be re-used several times.	25 mm / 1" dia. 30 mm dia. 40 mm dia. 50 mm / 2" dia. 1¼" dia. 1½" dia.	10 pcs. 10 pcs. 10 pcs. 10 pcs. 10 pcs. 10 pcs.	40300085 40300086 40300089 40300090 40300087 40300088
	<b>FlexiForm</b> Flexible silicone rubber mounting cups for acrylic cold mounting resins	25 mm dia. 30 mm dia. 40 mm dia. 1¼" dia. 1½" dia.	5 pcs. 5 pcs. 5 pcs. 5 pcs. 5 pcs.	40300018 40300019 40300022 40300020 40300021
	<b>FlexiForm</b> Translucent silicone rubber mounting cup for mounting of rectangular specimens. To be used with acrylic cold mounting resins where they can be reused many times. Epoxies reduces the lifetime as the resin over time affects the silicone rubber.	68 x 37 x 35 mm, fits to specimen holder MAXDI, 02606920 90 x 50 x 35 mm, fits to specimen holder MAXOT, 02606922 120 x 60 x 45 mm	3 pcs. 3 pcs. 2 pcs.	40300082 40300083 40300084
	<b>SeriForm</b> Polypropylene mounting cups for acrylic and polyester cold mounting resins	25 mm dia. 30 mm dia. 40 mm dia.	10 pcs. 10 pcs. 5 pcs.	40300007 40300008 40300009
	<b>FlangeForm</b> Mounting cups with flange for AccuStop	30 mm dia. 40 mm dia.	Set of 3 Set of 3	40300012 40300014
	Fixation clips	<b>MultiClips</b> Multiple plastic clip for holding up to 5 small and thin specimens when mounting		50 pcs.
<b>Metal spring clips</b> For holding small and thin specimens when mounting		6 mm dia. 9 mm dia.	100 pcs. 100 pcs.	40300026 40300025
<b>Taper Section Angle</b>		Aluminium Copper Steel	50 pcs. 50 pcs. 50 pcs.	40300070 40300071 40300069
<b>Mould Insert</b>		Mould insert for use with 40 mm dia. mounting cups Max specimen size 35 x 20 mm	50 pcs.	05756912
Measuring and Mixing		<b>Cold mounting Consumables Kit</b>	400 mixing cups, 400 stirrers	
	<b>Consumables Kit for CitoVac</b>	100 dispensing tubes and one chamber protector		40300080
	<b>Epovac Consumables Kit</b>	80 mixing cups, 80 stirrers, 80 tubes, 2 rubber plugs		40300030
	<b>Disposable syringes</b>	5 ml 10 ml 20 ml	100 pcs. 100 pcs. 100 pcs.	40300047 40300048 40300049
	<b>Electric mixer for mixing of epoxy</b>	1 mixer, 1 disposable propeller		05696101
	<b>Disposable propeller</b>		20 pcs.	40300072
	<b>Silicon oil, release agent</b>		100 ml	40300076
	Dyes	<b>EpoDye</b>		20 g
<b>AcryDye</b>		Blue, red and yellow dye	3 x 20 ml	40300081
Protection	<b>Protecting caps</b> Plastic cap for cylindrical specimen/mount	For 25 mm dia. specimens/mounts For 30 mm dia. specimens/mounts For 40 mm dia. specimens/mounts For 50 mm dia. specimens/mounts	100 pcs. 100 pcs. 100 pcs. 100 pcs.	49900033 49900034 49900035 49900036
	<b>Protecting Lacquer</b>	In spray can	400 ml	49900012
CitoVac	CitoVac with built-in vacuum ejector. Vacuum impregnation unit. Compressed air supply. 4.5-6 bar required.	1 x 100-240 V / 50-60Hz		05926119
	CitoVac for external pump. Vacuum impregnation unit. External vacuum pump required.	1 x 100-240 V / 50-60Hz		05926219
Wupty	Tool for pressing mounts out of FixiForm			05696901
Drybox-2	Combined specimen drier and curing oven with adjustable temperature (approx. 40-90°C) and low noise fan. With basket for specimens (04236904) and slot for specimen holders	1 x 110-120 V / 50-60 Hz		03706118
		1 x 220 V / 50 Hz		03706133
		1 x 220 V / 60 Hz		03706134